1.1. Tabla de contenidos

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1.2. General information

- Wattage: 3000 W
- Voltage:
D-4 LARGE 230 V
- Intensity:
D-4 LARGE 13 A
- Output: 4 lts/h
- Water consumption: 55 l / min - 1 l / min
- Temperature distilled water: 25 °C
- Temperature draining water: 50 °C
- Conductivity (at 20 °C): 2.3 µS/cm in optimal conditions

ATTENTION

Since distilled water is in contact with the atmosphere before the conductivity measurement can be carried out, its value is actually much less (up to 10 times less) than what it is actually measured. The measurement of any ultrapure water in contact with atmosphere has a conductivity of 1-2 µS/cm and a pH of 5, due to a small amount of CO₂ (0.05 ppm) absorbed from the atmosphere.

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1.3. Safety devices

- Magnetic switch that avoids water still worked in case insufficient water level.
- Safety thermostat protects still in the event of water supply failure.

1.4. Packing list

- Water still.
- Water inlet tube.
- One filter.
- One brace.
- Silicone tubes set.
- Two nuts Ø 32 with liner
- Two nuts Ø 18
- Two bent olives

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1.5. Recommendations

Electric:

- Electrical wiring must be, at least, of 4 mm² section.
- It must have ground.
- Electrical system must let intensity of 13 A pass.

Outlet:


- The tap water should have screw mouth of 3/4"

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1.6. Assembly

1. For your security, the water still is supplied separately: the metallic base and boiler and condenser, water inlet level and security level.
2. When you assemble the condenser (1) over the boiler (2), make it carefully fitting on the o-rings. The distilled water outlet (3) should be as you see in the bottom figure.
3. Assemble the water inlet level threading the nut (4) onto the boiler and the silicone tube onto the condenser.
4. Assemble the security level threading the nut (4) onto the boiler on its back side and connect the wire (5) to the base.
5. Close the draining tap (7) if it was opened - it is located on the back side of the water still -.
6. Connect the apparatus to the water mains with the supplied tube (6).
7. Connect the apparatus to the mains with the supplied electrical wire (8).

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
1.7. Controls

1. Switch ON / OFF (1).
2. Lighting pilot of NOT ENOUGH WATER (2)
3. Magnetic switch (3).

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1.8. Operation

1. Connect the equipment to the suitable power supply.
2. Switch the apparatus on.
3. Turn on the tap water. The filling of the apparatus starts.
4. When the water reaches the heaters, the equipment starts distilling when the water reaches the security switch.
5. The apparatus can work continuous during all the day.

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1.9. Maintenance

- You should clean the still regularly. So, in this way, the results will be better and you can get a still more durable.

- When you finished using the water still, it is recommended to empty the used water and fill again with clean water. You will save many cleanings.

- The periodicity of the cleanings depends on salts and hardness of the supplied water.


- It is strongly recommended to clean the filters placed in the water tube.

Steps to follow:

1. Unscrew the thread (1) and the magnetic switch from the base (2).

2. Add a softener (as LITOPOL) through the cleanout and clean the embedded impurities. Finally, turn the draining tap on (3).

3. Screw the thread (1) and connect the switch (2) to the base.

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

1.10 Troubleshooting

1. The water stills does not work:
 - Check the plug.
2. The red pilot "Not enough water" is light:
 - Check that there is water flow in the tubes.
 - Check that the dirt is not blocked the water flow.
3. The apparatus is filled with water and it does not starts working:
 - Check that the water flow is enough.
 - Check that the security switch is working properly.
4. For whatever queries, call manufacturer Technical Service.

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