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Warning

Warning alert you to a possibility of personal injury.

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Caution

Caution alert you to a possibility of damage to the equipment.

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Note

Notes alert you to pertinent facts and conditions.

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Hot surface

Hot surface alert you possibility of burning injury by hot instrument surface.

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Electrics: 230 V, Int. 26 A

Production: 4 L/hour

Conductivity: It is possible to reach in right conditions 1 J.S/cm.

Water consumption: 120 L/h.

You have to have wire section in your laboratory of, at least, 6 mm.

Safety presostat protects still in the event of water supply failure.

Automatic device to cut water supply in the event of electrical failure

Device to fill a container without worrying when it finished. It cuts water and electrical supply to avoid useless expenses of water and electricity.

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
1.4. Unpacking

1. Take out the heaters and unpack them.
2. Unpack the condensers.
3. Remove the salt from the condensers with water.
4. Take the water still with another person's help.

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
- 1- Feed and cooling water tubes boiler 1
- 2- Cooling water tube boiler 2
- 3- Monodistilled water outlet and feed of the 2nd bidistiller boiler
- 4- Bidistilled water outlet

1. Lift the front blue cover and take it out.
2. Open the lateral doors with a coin and push the door from inside.
3. The water still is supplied disassembled in order to lower the risk of breakage.
4. Please, place the labelled tubes as you can see on photo n. 1
5. Place the condensers inclined to the left.
6. Both condensers carry PTFE sleeves on the joints. You should not remove them.

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7. Place silicone tube to the outlet hose of the condenser on the left side.

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8. Place element heaters in the holes indicated.


Fig. 3.

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9. Place the two heaters in their holes and fix them with the screw thread (piece 1). When you have finished this operation, you must connect each plug.



Fig. 4.

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10. Close the lateral doors and the front blue cover.

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11. For the tap of general electrical outlet (5), make installation more convenient to the location. Be sure that electrical installation is suitable with electrical consumption.

12. Connect water hose with screw thread to the indicated place (see below 4) and the other end and filter to the main water.

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1. Connect the switches (6 and 7). They and the pilot (8) will illuminate. Turn the tap on until the pilot (8) turns off.

2. Turn slowly the tap on until the pilot (8) turns off. The distillation of water will start automatically.


Fig. 7.

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3. The water distilled in the first boiler moves to the second boiler. When the water reaches the level set by the level switch (10) the bidistillation will start.

4. This process will keep if there is enough flow of water and electricity.

5. The tubes (11) drive the draining water to the drain.

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6. Insert the device (see below piece 13) into selected vessel to rake the distilled water. When this vessel is full, the equipment will be disconnected.



Fig.8.

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1.7. Maintenance

- You should clean the still regularly. So, in this way, the results will be better and you can get a still more durable.

- The cleaning periodicity depends on the salts which water supply contains.

Steps to follow:

1. Loosen the nuts (14 and 15) and put acetic acid, hydrochloric acid or whichever another product into the boiler. The water should be hot.
2. You have to wait for the attack to the residues and empty and rinse the boiler.
3. Switch off the cleaning switch (7) and rinse with clean water a few times.
4. You should clean the filters of the tubes whenever you clean the still.



Fig. 9.

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ATTENTION

Since distilled water is in contact with the atmosphere before the conductivity measurement can be carried out, its value is actually much less (up to 10 times less) than what it is actually measured. The measurement of any ultrapure water in contact with atmosphere has a conductivity of 1-2 µS/cm and a pH of 5, due to a small amount of CO₂ (0.05 ppm) absorbed from the atmosphere.

At a general level, features on distilled water are:

	Conductivity (uS/cm)	Resistivity (MOhm · cm)	Silicates (mg/l)	Heavy metals (mg/l)	Reduced Permanganate (min)	Sodium (mg/l)	Calcium hardness (mg/l)	Ammonium (mg/l)	Bacteria (u.f.c./ml)	pH (at 25°C)
Once distilled	10 – 2	0,1 – 0,5	1 – 0,5	1 – 0,5	30	5 – 2	3 – 1	0,01	<10	5 – 6,5
Bidistilled	2 – 1	0,5 – 1	0,7 – 0,1	0,8 – 0,1	30	1 – 0,5	0,3 – 0,1	0,01	<10	5 – 6,5