VACIOTEM-T
VACIOTEM-TV

4001489
4001490

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GENERAL INFORMATION

- 1) *Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.*
- 2) *Do not install or use the equipment without reading this handbook before.*
- 3) *This handbook must always be attached to the equipment and it must be available for all users.*
- 4) *If you have any doubt or enquiries, please contact your supplier or J.P. Selecta's technical service.*
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY EQUIPMENT TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) *If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 89/655/CEE, the manufacturer is not responsible for the damage that can occur.*
- 7) **DO NOT USE THE EQUIPMENT WITH LIQUIDS WHICH CAN GIVE OFF VAPOURS CAPABLE OF MAKING EXPLOSIVE MIXTURES.**

PACKING LIST

The standard equipment consist of the following components:

Description	Code	Quantity
Vaciotem (T/TV) oven	4001489 / 4001490	1
Shelves		2
Instruction Manual	80152	1
Pump connection hose	46147	1

ACCESSORIES

Denomination	Code
VACUUM PUMP «VACUUM»	5900621

Rotary blade pump with built-in oil return.

Motor with thermic protection device and starting up switch.

EQUIPMENT DESCRIPTION

Vacuum drying ovens VACIOTEM-T and VACIOTEM-TV are designed for drying at adjustable temperatures from 35°C up to 200°C with a temperature fluctuation of $\pm 1^\circ\text{C}$, and $\pm 2^\circ\text{C}$ of homogeneity. Temperature regulation is made by an electronic controller which allows delayed start up and temporized run. The heating elements surround the working chamber. The electronic controller of the VACIOTEM-TV have furthermore the vacuum level control and a RS-232 input/output.

The cylindrical inner chamber is made in AISI 304 stainless steel with removable anodized aluminium shelves. Inner door in tempered glass is assembled on an elastic silicone gasket to absorb all the contractions and expansions which may occur.

TECHNICAL FEATURES

Voltage supply 115-230V 50/60 Hz according to the characteristics plate indications.

Code		4001489/4001490
<i>Maximum Vacuum</i>		10^{-2}mmHg
<i>Capacity (litres)</i>		47
<i>Working volume</i>		70%
<i>Inner dimensions (cm)</i>	\varnothing	34
	<i>Depth</i>	52
<i>Overall dimensions (cm)</i>	<i>Height</i>	54
	<i>Width</i>	76
	<i>Depth</i>	70
<i>Shelves</i>		2
<i>Consumption (W)</i>		2200
<i>Weight (kg)</i>		73

INSTALLATION

Put the oven on a flat, horizontal and level surface leaving a space of 10 cm at the back and on each side of the equipment.

Connect the vacuum pump or vacuum net inlet to the outlet placed at the rear part of the oven (30) and turn the manual valve on (29) (the handle have to be positioned perpendicularly to the equipment). Connect the plug at the socket which is on the rear part of the oven (28). Turn the vacuum switch on.

If the work product can produce vapours, you should only use with the oven, water jet pumps or hydrotative pneumatic pumps. If an air ballast pump is used, a water refrigerated condenser to deposit gases and do not damage the pump has to be assembled between the pump and the oven. Besides, a condensation collector has to be added on each side of the pump. (suction and expulsion)

For maximum efficiency, do not place the oven in direct sunlight or other sources of heat.

Connect the oven at the mains. Be sure that the voltage supply is the same as the one indicated on the characteristics plate. Do not use the equipment if it is NOT earthed.

If you have to change the plug cap, keep in mind the following:


Blue cable: Neutral.

Brown cable: Phase.

Yellow/green cable: Earth.

OPERATION

CONTROL PANEL (VACIOTEM-TV):

1. RS-232 socket.
2. Air connection inlet.
3. Ventilation valve.
5. Vacuum set indicator led.
6. Breaking vacuum at the end of the cycle led.
7. Vacuum run status indicator led.
8. Pump indicator led.
9. Vacuum level reading display.
10. Vacuum selector push button.
11. End of cycle electrovalve push button.
12. Increasing value push button.
13. Reducing value push button.
14. Vacuum START/STOP push button.
15. Temperature Mode indicator led.
16. Time Mode indicator led.
17. Run status indicator led.
18. Standby time indicator led.
19. Time/temperature reading display.
20. Temperature selecting push button.
21. Time selecting push button.
22. Increasing value push button.
23. Reducing value push button.
24. Temperature START/STOP push button.
25. Main switch (POWER).
26. Safety thermostat indicator lamp (red).

CONTROL PANEL (VACIOTEM-T):

- 2. Air connection inlet.
- 3. Ventilation valve.
- 4. Vacuum gauge.
- 15. Temperature Mode indicator led.
- 16. Time Mode indicator led.
- 17. Run status indicator led.
- 18. Standby time indicator led.
- 19. Time/temperature reading display.
- 20. Temperature selecting push button.
- 21. Time selecting push button.
- 22. Increasing value push button.
- 23. Reducing value push button.
- 24. Temperature START/STOP push button.
- 25. Main switch (POWER)
- 26. Safety thermostat indicator lamp (red).
- 27. Vacuum pump switch.

REAR PART (T/TV MODELS):

- 27. Air connection valve (only TV model)
- 28. Pump power socket.
- 29. Vacuum handle valve.
- 30. Vacuum connector.
- 31. Safety thermostat.
- 32. Safety thermostat reset push button.
- 33. Fuses.

STARTING UP (TV MODEL):

1. Check that the installation protocol has been correctly followed.
2. Check that the ventilation valve (3) is closed.
3. Load the material that you have to work. Never load more than 70% of the space in the chamber. Be sure that hole in the rear part of chamber is free.
4. Close the oven door.
5. Switch on the main switch (25). The displays (9) and (19) will light.
6. Set vacuum parameters (see "Vacuum configuration parameters")
7. Set temperature parameters (see "temperature configuration parameters")

Vacuum configuration parameters:

The reading display (9) shows absolute pressure expressed in mbar (HPa). So when the oven door is opened, it shows the current atmospheric pressure (no 0).

1. Check that the RUN led indicator (7) does not light. Otherwise, push the START/STOP push button (14).
2. Push "mbar" (10) and at same time push (12) or (13) to set the work value. While you are doing this operation "mbar" led indicator (5) will light.
3. Push START/STOP (14) and the vacuum process will start. The led indicator (8) will light when the vacuum pump would be working.
4. The vacuum process can be stopped by pushing START/STOP (14). Otherwise, the process will continue until the end of the time set in the temperature controller (see "Temperature configuration parameters").

NOTE: At the end of the process, there are two possibilities:

- Maintaining the vacuum (led "6" off)
- Breaking the vacuum (led "6" on)

This can be chosen by pressing the key (11) which enables or disables this possibility in any moment of the process.

Temperature and time configuration parameters:

The heating process can run with or without vacuum. So, it is important the following:

During the heating, the manual valve (3) placed on the control panel, have to stay opened to avoid overpressure inside

The oven can work in the following modes:

MANUAL MODE:

1. Push (10) and, at same time, push (22) or (23) to set the work value. While you are doing this, the operation led indicator (15) will light.
2. Push START/STOP (24) and the heating will start. The oven start up and the RUN indicator led will light.
3. To disconnect the heating push again START/STOP (24) push button, the RUN indicator led (17) will turn off.

TEMPORIZED RUN MODE:

1. Follow 1 item of MANUAL MODE.
2. Push (21) and, at the same time, (22) or (23) to prefix the run time needed. During this operation, the indicator led TIME MODE (16) will light.

NOTE: The time range (between 1 minute and 9 h 59 min) is expressed in sexagesimal units. From 10 hours, time is expressed in hours and decimal fractions of them (each decimal equal 6 minutes)

3. The run time will begin to discount when chamber temperature would reach the set temperature.
4. Push START/STOP (24). Oven will start up and RUN indicator led (17) will blink until chamber temperature reach the set one. Then RUN indicator led (17) will light permanently.
5. When prefixed time is thoroughly elapsed, the process will stop and the RUN indicator led (17) will turn off. If user would select "End of cycle electrovalve push button" (see "Vacuum parameters configuration") the work vacuum level will continue until user would push START/STOP (14). Otherwise the chamber will automatically pressurize.

NOTE:

* In any moment user can stop the process by pushing START/STOP (24).

* If during work time you need to see set time you must push (21).

* If you need to see passed time or remaining time, you must push (22) or (23) keys.

DELAYED START UP MODE:

1. Follow 1 item of MANUAL MODE.
2. Push (21) twice. Maintain it pushed and at the same time push (22) or (23) to set delay time needed (between 1 and 24 hours). During this operation, the indicator led TIME MODE (16) will light.
3. If furthermore users need to set a run time, they have to follow the instructions "Temporized time configuration".
4. Push START/STOP (24). The RUN indicator led (17) and TIMER (18) will light. The remaining time will be able to see on the display (19). Each decimal digit will express the tenth part of an hour (6 minutes).

STARTING UP (T MODEL):

1. Check that the installation protocol has been correctly followed.
2. Check that the ventilation valve (3) is closed.
3. Load the material that you have to work. Never load more than 70% of the space in the chamber. Be sure that the hole in the rear part of chamber is free.
4. Close the oven door.
5. Switch on the main switch (25). The display (19) will light.
6. Switch on the vacuum switch VACUUM PUMP (27).
7. When the vacuum level is the needed one (see vacuum gauge) (4), disconnect the vacuum pump or the vacuum net with the VACUUM PUMP switch (27).

NOTE: In this vacuum gauge you can read absolute pressure so with the door the reading is 0.

8. Set temperature parameters (see "temperature configuration parameters" of instruction STARTING UP TV MODEL).
9. When temperature process ends, it is recommended to pressurize the chamber with the control panel manual valve (3) in order to avoid that the pressure induction put oil from vacuum pump into the chamber.

SAFETY THERMOSTAT

Vacuum ovens Vaciotem T and Vaciotem TV are equipped with a safety thermostat (31) which prevents that heaters temperature never reach at dangerous levels.

Never change the thermostat adjustment. It could cause unexpected stops of drying processes.

Thermostat have a reset push button (32) that allows to let the equipment ready to work elapsed some time. Nevertheless, if thermostat disconnects many times the process, user should call up to technical service.

RS-232 COMMUNICATION

RS-232 COMMUNICATION OVERVIEW

- 9600 bauds, 7 dat bits, 1 bit stop, without parity no handshake.
- Up to 15m of cable.
- Cable conection DTE to DTE 9 pins male.
- Protocol ASCII with LRC against transmission errors
- MASTER-SLAVE Communication (PC computer as a master).
- The device never sends anything if computer does not ask for it.

Start frame	Address		Data	LRC		End frame
#	N	N	Message (n bytes)	K	K	CR LF
	H	L		H	L	

Start frame	#	=	ASCII (23)	1 byte
Address	NN	=	01...99	2 bytes
Data	Message	=	Command + arguments	n bytes
LRC	KK	=	00..FF	2 bytes
End frame	<CR><LF>	=	ASCII (0D)+ASCII (0A)	1 byte

- Messages start with character #.
- Device address must be included.
This field goes from 00 to 99, size always 2 bytes in decimal-ASCII format.
- Data field is the message.
- LRC is a calculated value to check good transmission.
- At the end of the frame always a CR (0D) followed by LF (0A).
- Slave device answers always by sending its address.

LRC ERROR CHECKING FIELD

Frame includes a transmission error checking field called LRC. The field size is 2 bytes. This field is calculated by the sender device and includes it on the frame. The receiver device also calculates LCR value. If calculated value is equal than the transmitted LCR value, transmission was good.

LRC calculation procedure:

1. Convert each character on its hexadecimal value.
2. Add.
3. Take the smallest byte (two first characters starting from right hand side).
4. Subtract from FF on previous result.
5. Add 1.

ASCII	HEX
#	23
0	30
1	31
S	53
V	56
T	54
<BL>	20
+	2B
0	30
9	39
8	38
,	2C
0	+ 30
	2C9

LRC CALCULATION EXAMPLE

#01SVT<BL>+098,037<CR><LF>

		FF
	-	C9
		36
	+	1
		LRC = 37

- «Get» commands always ends by «?»
- Characters always in capital letters.
- Use spaces only when indicated by <BL>.
- Decimal separator is «,».

CONNECTION CABLE

It could be found in any electronic components specialized shop:

Cable Null - 9 pins modem

PC	Device
1	7, 8
2	3
3	2
4	6
5	5
6	4
7	1
8	

Connectors in both cable sides are 9 pins female.

RS-232 INSTRUCTIONS FOR VACIOTEM T AND VACIOTEM TV OVENSVaciotem T / Vaciotem TV

PVT? Get probe temperature PV (Point value):
 Command: #01PVT?KK<CR><LF>
 Answer: #01+000KK<CR><LF>

Vaciotem TV

PVP? Get pressure value
 Command: #01PVP?KK<CR><LF>
 Answer: #01000KK<CR><LF>

Vaciotem T / Vaciotem TV

SVT? Get Set point temperature:
 Command: #01SVT?KK<CR><LF>
 Answer: #01+000KK<CR><LF>

SVT Set point temperature
 Command: #01SVT<BL>+000KK<CR><LF>
 Answer: #01OKKK<CR><LF>

Vaciotem TV

SVP? Get Set point temperature:
 Command: #01SVP?KK<CR><LF>
 Answer: #01000KK<CR><LF>

SVP Set point temperature
 Command: #01SVP<BL>+0000KK<CR><LF>
 Answer: #01OKKK<CR><LF>

Vaciotem T / Vaciotem TV

RUN? Get status variable
 Command: #01RUN?KK<CR><LF>
 Answer:
 Temperature stop: #01STPTKK<CR><LF>
 Pressure stop: #01STPPKK<CR><LF>
 Temperature run: #01RUNTKK<CR><LF>
 Pressure run: #01RUNPKK<CR><LF>
 Alarm status: #01ALARKK<CR><LF>
 Standby time: #01STBYKK<CR><LF>

Model Vaciotem TV gives compound answers. For exemple:

#01STPT<BL>STPPKK<CR><LF>

Vaciotem TV

RUN Starts run:
 Command: #01RUNKK<CR><LF>
 Answer: #01OKKK<CR><LF>

Vaciotem T / Vaciotem TV

RUNT Starts temperature control:
 Command: #01RUNTKK<CR><LF>
 Answer: #01OKKK<CR><LF>

Vaciotem TV

RUNP Starts pressure control:
 Command : #01RUNPKK<CR><LF>
 Answer : #01OKKK<CR><LF>

Vaciotem TV

STOP Stops device:
 Command : #01STOPKK<CR><LF>
 Answer : #01OKKK<CR><LF>

Vaciotem T / Vaciotem TV

STPT Stops temperature control:
 Command : #01STPTKK<CR><LF>
 Answer : #01OKKK<CR><LF>

Vaciotem TV

STPP Stops device:
 Command : #01STPPKK<CR><LF>
 Answer : #01OKKK<CR><LF>

Vaciotem T / Vaciotem TV

ID? Gets device identification:
 (Software version and device code)

Command : #01ID?KK<CR><LF>
 Answer : #014001489V01KK<CR><LF> (Vaciotem-T)
 #014001490V01KK<CR><LF> (Vaciotem-TV)

7 characters for code + V + 2 for the version

Vaciotem T / Vaciotem TV

SAL? Gets alarm status:
 Command : #01SAL?KK<CR><LF>
 Answer : #01ALARM1KK<CR><LF>

* See table of the alarm status.

Vaciotem T / Vaciotem TV

TBY? Gets standby time:
 Command : #01TBY?KK<CR><LF>
 Answer : #0100,0KK<CR><LF>

SBY Sets standby time:
 Command : #01TBY?KK<CR><LF>
 Answer : #0100,0KK<CR><LF>
 Time in hours and tenth parts of an hour
 Between 0 and 24 hours

TBR? Gets standby set time:
Command : #01TBR?KK<CR><LF>
Answer : #0100,0KK<CR><LF>

Vaciotem T / Vaciotem TV

TRU? Gets run time:
Command : #01TRU?KK<CR><LF>
Answer : #010,00KK<CR><LF>
Answer : #0100,0KK<CR><LF>
Time until 9h59mim (9,59)
Time from 10h until 9.99 (99,9)

TRU Sets run time:
Command : #01TRU<BL>0,00KK<CR><LF>
Answer : #01TRU<BL>00,0KK<CR><LF>
Answer : #01OKKK<CR><LF>
Time until 9h59mim (9,59)
Time from 10h until 9.99 (99,9)

Vaciotem TV

INT Inhibits/enable key pad (temperature/pressure):
Command : #01INT<BL>1KK<CR><LF>
Answer : #01OKKK<CR><LF>
1 = inhibit 0 = enable

Vaciotem T/Vaciotem TV

ITT Inhibits/enable key pad (temperature/pressure):
Command : #01INT<BL>1KK<CR><LF>
Answer : #01OKKK<CR><LF>
1 = inhibit 0 = enable

Vaciotem TV

ITP Inhibits/enable key pad (temperature/pressure):
Command : #01INT<BL>1KK<CR><LF>
Answer : #01OKKK<CR><LF>
1 = inhibit 0 = enable

Vaciotem TV

SRLK Opens pump valve :
Command : #01SRLKKK<CR><LF>
Answer : #01OKKK<CR><LF>

RRLK Closes pump valve :
Command : #01RRLKKK<CR><LF>
Answer : #01OKKK<CR><LF>

Set commands			
COMMAND	CONCEPT	MESSAGE	ANSWER
SVT	Set Temperature (°C)	#NNSVT<BL>+/-000,0KK<CR><LF>	#NNOKKK<CR><LF> *
SVP	Pressure (mbar)	#NNSVP<BL>0000KK<CR><LF>	#NNOKKK<CR><LF> *
RUN	Run	#NNRUNKK<CR><LF>	#NNOKKK<CR><LF> *
RUNT	Run temperature	#NNRUNTkk<CR><LF>	#NNOKKK<CR><LF> *
RUNP	Run pressure	#NNRUNPKK<CR><LF>	#NNOKKK<CR><LF> *
STOP	Stop	#NNSTOPKK<CR><LF>	#NNOKKK<CR><LF> *
STPT	Stop temperature	#NNSTPTkk<CR><LF>	#NNOKKK<CR><LF> *
STPP	Stop pressure	#NNSTPPKK<CR><LF>	#NNOKKK<CR><LF> *
TBY	Standby time	#NNTBY<BL>00KK<CR><LF>	#NNOKKK<CR><LF> *
TRU	Run time	#NNTRU<BL>00,0KK<CR><LF>	#NNOKKK<CR><LF> *
SRLK	On pump	#NNSRLKKK<CR><LF>	#NNOKKK<CR><LF> *
RRLK	Off pump	#NNTRUKKK<CR><LF>	#NNOKKK<CR><LF> *
INT	Inhibit key pad	#NNINT<BL>0KK<CR><LF>	#NNOKKK<CR><LF> *
INTT	Inhibit temperature key	#NNINTT<BL>0KK<CR><LF>	#NNOKKK<CR><LF> *
INTP	Inhibit pressure key	#NNINTP<BL>0KK<CR><LF>	#NNOKKK<CR><LF> *

Read commands			
COMMAND	CONCEPT	MESSAGE	ANSWER
PVT?	Temperature PV (°C)	#NNPVT?KK<CR><LF>	#NN+000KK<CR><LF>
PVP?	Pressure PV (mbar)	#NNPVP?KK<CR><LF>	#NN0000KK<CR><LF>
SVT?	Temperature set (°C)	#NNPVT?KK<CR><LF>	#NN+000KK<CR><LF>
SVP?	Pressure Set (mbar)	#NNSVP?KK<CR><LF>	#NN0000KK<CR><LF>
RUN?	Device status	#NNRUN?KK<CR><LF>	#NNRUNKK<CR><LF>
ID?	Device identity	#NNID?KK<CR><LF>	#NN>KKKKKKKVKK<CR><LF><CR><LF>
SAL?	Alarms status	#NNSAL?KK<CR><LF>	#NNALARM0KK<CR><LF>
TBY?	Standby time	#NNTBY?KK<CR><LF>	#NN00,0KK<CR><LF>
TBR?	Standby Run time	#NNTBR?KK<CR><LF>	#NN00,0KK<CR><LF>
TRU?	Run time	#NNTRU?KK<CR><LF>	#NN0,00KK<CR><LF>

* The write error commands are on the error table.

Alarm status		
ALARM	CONCEPT	ANWER
0	Any alarm	#NNALARM0KK<CR><LF>
1	Overtemperature	#NNALARM1KK<CR><LF>
2	Power fail	#NNALARM3KK<CR><LF>

SPARE PARTS

To guarantee the safety of the equipment, the spare parts must be purchased from J.P.SELECTA, s.a.u.

Code	Description	4001489	4001490
21066	Chamber gasket	x1	x1
16003	Vacuum gauge	x1	x0
16214	Air inlet electrovalve	x0	x1
16215	Vacuum electrovalve	x1	x1
39005	Heating element	x2	x2
43017	Thermostat	x1	x1
36027	Rolid Relay	x1	x1
29268	Temperature control board	x1	x0
29269	Temperature-vacuum control board	x0	x1
15551	Fuse	x2	x2

MAINTENANCE

Before removing the cover, disconnect the equipment from the mains.

The manipulation of the internal electronic circuits of the oven by unauthorized personnel can cause irreparable damage. Take it to one of the J.P.SELECTA, s.a.u. authorized technical services.

CLEANING:

For the cleaning of the different parts, use the following products:

Cleaning of stainless steel: Alcohol

Cleaning of plastic: Alcohol with cotton duster.

GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.

Any manipulation of the equipment by unauthorized personnel by J.P. SELECTA, s.a.u. cancels automatically the guarantee.