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***INCUBATOR FOR ANAEROBIC CELL  
AND TISSUE CULTURES***

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# Table of contents

<b>1 Packing list</b> .....	<b>14</b>
1.1 Standard package.....	14
<b>2 Accessories</b> .....	<b>14</b>
<b>3 Safety</b> .....	<b>14</b>
3.1 Safety icons .....	14
3.2 Risks that the operator can be subjected to:.....	15
3.3 Qualified users .....	15
3.4 Modifications .....	15
<b>4 Description</b> .....	<b>15</b>
<b>5 Installation</b> .....	<b>15</b>
5.1 Suitable location.....	15
5.2 Connection to power supply .....	16
5.3 Connection to CO2.....	16
<b>6 Technical features</b> .....	<b>16</b>
<b>7 Operation</b> .....	<b>17</b>
7.1 Control panel.....	17
7.2 Rear part.....	17
7.3 Starting up .....	18
7.4 Parameter configuration .....	18
<b>8 Alarms</b> .....	<b>20</b>
<b>9 RS-232 Communications</b> .....	<b>20</b>
<b>10 Maintenance</b> .....	<b>21</b>
10.1 Spare list.....	21
10.2 Microfilter repalcing.....	21
10.3 Cleaning CO2 regulation jet.....	22

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## 3 Safety

This equipment has safety devices.

This manual indicates areas where risks are possible.

### 3.1 Safety Icons

We have marked with a grey line the risky situations and the safety measures that should be adhered to.

#### Risk of danger

Danger risk.  
Beware, follow the instructions as indicated.

#### Electrical risk

Risk of electrical shock if accessing zones that are shown by this icon.  
Beware, follow the instructions as indicated.

#### Risks of burns in high temperature areas

The temperature in the zone indicated with this icon can exceed 60°C. Use thermal protective gloves to complete described functions.

Beware, follow the instructions as indicated.

#### Important information

- Important information in obtaining the best results or optimum performance of the equipment.
- Important information to extend the life of the equipment and maintaining its optimal performance.

### 3.2 Risks that the Operator can be Subjected to

- Possibility of touching surfaces that are over 60°C.
- Electrical risk

### 3.3 Qualified Users

This equipment must only be used by qualified personnel.

This equipment must only be used by personnel that has read and understood this manual.

### 3.4 Modifications

Modifications of the safety devices are not allowed by the manufacturer and may cause unexpected risks.

**DO NOT OPEN THE DOOR OR TOUCH THE GLASS DOOR DURING THE STERILIZATION PROCESS**


## 4 Description

The external enclosure is epoxy coated and it has a thermal isolation around the chamber. The chamber is made of stainless steel with shelves which can be easily removable, for an easy clean. It has a double door, an inner one made from tempered glass with a silicone join and another external door with a magnetic lock and a heating system to avoid condensations.

Users will be doctors or their assistants under their supervision.

The incubator has a CO<sub>2</sub> inlet with a microfilter on the rear part. And an outlet on the control panel to test the CO<sub>2</sub> concentration.

Electronic control of the temperature and of the CO<sub>2</sub> concentration. Humidity level constant at 98% produced by the evaporation of the water inside the optional tank.


## 5 Installation

Once the equipment is unpacked, make a visual inspection, check the packing list and make sure that the equipment does not got any damages during the transport.

Keep the packing for some days just in case it was necessary to return the equipment.

### 5.1 Suitable Placement

Place the equipment on a leveled and flat surface, enough strong to support its weight, leaving a free space of 20 cm around the incubator.

## 5.2 Connection to Power Supply

Place the equipment near the socket of 16 A.

Due to safety reasons, the socket must be connected to the ground.

Before connecting the plug to the socket, please check that the voltage marked on the unit plate matches the voltage of the socket.

## 5.3 CO<sub>2</sub> Connection

The CO<sub>2</sub> cylinder must be equipped with two pressure gauges. The first one, the nearest to the cylinder, has to indicate the pressure of the gas inside it, and it must be checked often in order to avoid exhausting the CO<sub>2</sub> during work. The second one must be adjusted by a pressure reducer in order to maintain the pressure between 0.5 and 1.5 bar.

Connect the CO<sub>2</sub> cylinder through a gas reducer to the rear part inlet (6) using the supplied hose.


## 6 Technical Features

Overall dimensions (cm):	
height	89
width	64
length	60
Chamber dimensions(cm):	
height	65
width	50
length	46
Capacity of shelves:	9
Capacity (litres):	150
Incubation power (W):	1100
Sterilization power (W):	2025
Weight (kg):	110


## 7 Operation

### 7.1 Control Panel

- 1 Single button control
- 2 CO<sub>2</sub> sample outlet
- 3 RS-232 connector
- 4 Screen
- 5 Alarm indicator led

### 7.2 Rear Part

- 6 CO<sub>2</sub> inlet
- 7 Microfilter
- 8 24V DC 100 mA alarm outlet
- 9 Safety thermostat reset button
- 10 Fuse carrier
- 11 Power supply cord
- 12 Main switch

### 7.3 Stating Up

- 1 Check that the equipment is connected to the power supply, connected to the CO<sub>2</sub> cylinder that outlet pressure of the CO<sub>2</sub> is adjusted at a pressure between 0.5 to 1.5 bar by the gas pressure reducer.
- 2 If working with humidity is needed, fill in the water tank and put it inside, at the bottom of the chamber.
- 3 Turn on the main switch (12). During the first 120 seconds, the microprocessor does not make CO<sub>2</sub> readings. During this time, the display shows "--" instead of CO<sub>2</sub> %. The temperature reading is visible immediately. Nevertheless, during this period, the working parameters can be modified and the incubator can be started on.


°C	CO <sub>2</sub>	STATUS
19.5	--	STOP


- 4 To change from STOP mode to RUN mode, press the single button (1) for a while. The incubator starts to regulate the temperature and the CO<sub>2</sub> concentration at the last used values.
- 5 When the incubator is in RUN mode, press the single button (1) to change the incubator to SHOW mode. It shows the set temperature and the CO<sub>2</sub> concentration for about 5 seconds without stopping the operation. Then it goes on with the normal operation.
- 6 To stop the incubator, press the single button (1) for at least 5 seconds. The display will show the stop menu. Select "Yes" by turning the single button (1) and pressing on it.

STOP?	<NO>
YES	

### 7.4 Configuration of Parameters

- 1 To access the configuration menu, check that the equipment is in STOP mode.
- 2 Then press the single button (1) for 5 seconds.
- 3 The menu options are as follows:

### Limit of the Parameters:

#### SETTINGS:

TEMPERATURE: Between 10.0 and 50.0 °C  
 CO<sub>2</sub>: Between 0,0 and 20,0 %

#### STERILIZATION:

TEMPERATURE: Between 140 and 180 °C  
 TIME: Between 1 and 30 minutes

#### CONFIGURATION:

ALARM:            MAXIMUM  
                     TEMPERATURE    Low limit 55 °C and < SET

                    MINIMUM  
                     TEMPERATURE    Low limit 10 °C and > SET

                    MAXIMUM  
                     % CO<sub>2</sub>                    Low limit > SET and 21%

                    MINIMUM  
                     % CO<sub>2</sub>                    Low limit 0% and < SET

DATE:            Date and time

LANG:            Spanish and English

PRINTER:        PRINT  
                     TIME            Between 1 h and 24 h

OFFSET TEMP    Between -5 and +5 °C

OFFSET CO<sub>2</sub>    Between -2 and +2 °C



## 8 Alarms

Some fails can happen during the operation, and will be shown as a beep and a message in the display.

To avoid unnecessary alarms, errors 01, 04, 05 and 06 must be checked after 90 minutes from the starting or from the opening.

Alarms are established as follows:

- 01 Open door for more than 20 seconds. While the door is open the CO<sub>2</sub> supply is cut. When alarm 01 is produced, the heating of the chamber is cut but not the heating of the door. The error disappears when closing the door.
- 02 The temperature of the chamber exceeds the alarm temperature set by the user for more than 7 minutes. The consequence is that all the main heating is cut. The beep stops when pushing the single button (1)
- 03 Low CO<sub>2</sub> supply pressure for more than 20 seconds. This is an informative alarm. It does not have any effect on the function. When pushing the single button (1), the alarm stops but it appears again if the problem persists.
- 04 The measured CO<sub>2</sub> concentration value exceeds the maximum fixed by the user. Close the CO<sub>2</sub> solenoid valve. The audible alarm stops after some minutes but the message only disappears when pushing the single button (1).
- 05 The chamber temperature is lower than the minimum programmed value. This is an informative alarm. It does not have any effect on the function. The error message will disappear when the machine returns to the normal levels by pushing the single button (1).
- 06 The measured CO<sub>2</sub> concentration value is lower than the minimum fixed by the user. This is an informative alarm. It does not have any effect on the function. The error message will disappear when the machine returns to the normal levels by pushing the single button (1).
- 07 CO<sub>2</sub> valve opened for more than 7 minutes uninterrupted. It is considered that the chamber might have a leakage loss. This is an informative alarm. It does not have any effect on the function. The audible alarm stops after some minutes but the message only disappears when pushing the single button (1).
- 08 The difference between the two probes is higher than 2°C for more than 5 minutes. This alarm leaves the incubator in STOP status. This alarm indicates that one of the probes is not working correctly or that it is not well connected.

The alarm status produces an intermittent acoustic signal for 15 seconds. Furthermore, it gives an electrical signal of 24V DC 100 mA in the terminal block (8) on the rear part where it is possible to connect a buzzer, a lamp or other device the user could have.

## 9 RS-232 Communications

Through the connector RS-232 (3), it is possible to:

- Record the temperature, CO<sub>2</sub> concentration level and status in a PC.
- Manage the control panel from a PC locking the panel of the incubator.


Ask JP SELECTA s.a.u. for more information.

## 10 Maintenance

### 10.1 Spare Parts List


Code	Description
07319	Screen cord
07323	Encoder cord
07184	Printer cord
13007	Contactor (relay)
15520	Microfilter 0.2 micromilimeters
16058	Solenoid valve NC
16062	Encoder
16220	Pressure switch
20017	Main switch
21211	Joint glass door
21212	Joint external door
24017	Fan motor
29232	RS-232 circuit board cord
29325	Control and power circuit board
29348	Screen
34167	CO <sub>2</sub> regulation jet
36028	Solid relay
39121	Chamber heater element
39292	Door heater element
43017	Safety thermostat
43065	Chamber Pt100 probe
43066	Door Pt100 probe
43136	CO <sub>2</sub> sensor
45011	Mains filter 6.2 A
45076	Tranformator 115-115V/20-14-9V
45078	Transformator peltier
47188	Glass door

### 10.2 Microfilter Replacing

Each time the CO<sub>2</sub> cylinder is replaced, it is advisable to replace the microfilter (7) too. To do it:

- Cut the clamps.
- Extract the both ends of the microfilter from the silicone hoses.
- Insert the new microfilter.
- Lock it again with clamps.

### 10.3 Cleaning CO<sub>2</sub> Regulation Jet

If the incubator is not able to achieve the desired CO<sub>2</sub> concentration despite any CO<sub>2</sub> in the bottle and the valve itself, you should proceed to clean the jet:

- Loosen the CO<sub>2</sub> inlet with a spanner.
- Extract the part with the regulation jet.
- Dismantle the regulation jet from the inlet part.
- Clean the CO<sub>2</sub> regulation using a thin wire (<0.2mm) or pressured air in counter way that the CO<sub>2</sub> in the normal use. As indicated on the picture.

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