**CINTA CALEFACTORA DE GRAN SUPERFICIE  
LARGE SURFACE HEATING MANTLE**

**4000200**

## INDICE

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## GENERAL INFORMATION

- 1) Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.
- 2) Do not install or use the equipment without reading this handbook before.
- 3) This handbook must always be attached to the equipment and it must be available for all users.
- 4) If you have any doubts or enquiries, please contact with your supplier or J.P. Selecta's technical service.
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY EQUIPMENT TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 2009/104/CE, the manufacturer is not responsible for the damage that can occur.
- 7) Do not use the equipment with liquids which can give off vapours capable of making explosive mixtures.

## PACKING LIST

The standard equipment consists of the following components:

Piece	Code	Qty.
Heating mantle	4000200	1
Instruction manual	80415	1

## TECHNICAL FEATURES

Voltage supply 115-230V 50/60 Hz according to the characteristics plate indications.

Power consumption per band	Weight	Overall dimensions cm		Heating power
		Length	Width	W/cm <sup>2</sup>
1000 W	1.5 Kg	171	15	0.39

## DESCRIPTION

The heating mantles are flexible silicone bands insulated in mineral wool that allow, by means of a spring fastening system, to adapt to 200-liter metal drums with a diameter of 60cm to heat them before use.

The bands are equipped with a Schuko plug to connect directly to the power supply. The maximum temperature that can be reached on the band inner surface is 210°C.

## INSTALLATION

It does not require any special installation. The bands are mounted by simply hugging the drum, making sure that the band is as flat as possible for a better thermal contact.

Once installed, plug them into the mains.

**CAUTION!!! IMPORTANT FOR YOUR SAFETY**

**Be sure that the voltage supply is the same as the one indicated on the characteristic plate of the equipment.**

**Do not use the equipment if it is NOT earthed.**

**If you have to change the plug, bear in mind the following:**

**Blue cable: Neutral.**

**Brown cable: Phase.**

**Yellow/green cable: Earth.**

**Note:**

The equipment does not have neither a control thermostat nor a safety thermostat.

It is recommended never to leave the installation unattended.

**Note:**

For a correct operation, when assembling the heating bands they cannot overlap each other, as it would generate a high temperature at that point and could damage the bands.

## GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.

Any manipulation of the equipment by unauthorized personnel by J.P. SELECTA, s.a.u. cancels the guarantee automatically.

### CAUTION

The ceramic tube envelope can reach temperatures of about 80°C

### Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.