**V-1100****4120025****Espectrofotómetro de rango visible  
Spectrophotometer visible range**

# Indice

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# Index

<b>1 INTRODUCTION .....</b>	<b>16</b>
1.1 Package contents .....	16
1.2 Appearance.....	16
1.3 Operating Panel .....	17
1.4 Descripción del Teclado .....	18
<b>2 SAFETY .....</b>	<b>18</b>
<b>3 OPERATING PRINCIPLE .....</b>	<b>19</b>
<b>4 INSTALLATION .....</b>	<b>19</b>
<b>5 OPERATION .....</b>	<b>20</b>
5.1 Local Software Structure .....	20
5.2 Basic Operation .....	20
<b>6 MEASUREMENT.....</b>	<b>21</b>
6.1 Measure the Absorbance .....	21
6.2 Measure the Transmittance.....	22
6.3 Two-point Method.....	23
6.4 Coefficient Method.....	23
<b>7 SPECIFICATIONS.....</b>	<b>25</b>
<b>8 WARRANTY .....</b>	<b>25</b>
<b>9 TROUBLESHOOTING.....</b>	<b>26</b>
<b>10 CE CONFORMITY DECLARATION / DECLARACIÓN DE CONFORMIDAD .....</b>	<b>27</b>

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## **ENGLISH VERSION**

Contact with  
dangerous  
substances

Electrical accident

## 1 INTRODUCTION

Spectrophotometers V-1100 described in this manual are designed to be used by properly trained personnel in a suitable equipped laboratory. For the correct and safety use of this equipment, it is essential that laboratory personnel follow generally accepted safe procedures in addition to the safety precautions called for in this manual.

The covers on this instrument may be removed for servicing. However, the inside of the power supply unit is a hazardous area and its cover should not be removed under any circumstances. There are no serviceable components inside this power supply unit.

### 1.1 Package contents

Carefully unpack the contents and check the materials against the following packing list to ensure that you have received everything in good condition.

#### Packing List

Description	Quantity
Spectrophotometer .....	1
Mains Lead .....	1
Glass cuvettes .....	4
Manual .....	1


### 1.2 Appearance

Please see Fig. 3-1:

Fig. 3-1

- 1 Lid of the compartment
- 2 Cell Holder
- 3 Pole
- 4 Operating Panel
- 5 Wavelength Knob
- 6 USB Port
- 7 Print Port
- 8 Power Socket 85-265VAC
- 9 Power Switch

### 1.3 Operating Panel

Please see the panel of V-1100 (Fig. 3-2):


Fig. 3-2

- 1 LCD (128 X 64)
- 2 Key Button

Switch the Photometric Mode

Confirm/Print

Decrease Number/Set Zero

Increase Number/Set 100%T

## 2 SAFETY

Spectrophotometer VR-1100 is an electrical test instrument widely used in laboratories.

This equipment has been designed to be used mainly in laboratories.

It is important that the laboratory staff know the safety general procedures, together with the safety precautions explained in this manual.

Some of the chemical products used in this spectrophotometer are corrosive and/or inflammable and the samples can be radioactive, toxic or potentially infectious.

Be careful to follow the safety standardized procedures in laboratories for samples and chemical products manipulation.

DANGER:  
contact with  
dangerous  
substances

DANGER:  
Electrical accident

#### Electrical connection:

Power requirements of V-1100 are from 85V to 265VAC. Make sure that the local power supply is within this range.

The power cord shall be inserted in a socket provided with a protective earth contact.

#### Warning:

Any interruption of the protective conductor inside or outside the equipment is likely to make it dangerous.


NEVER touch or handle the power supply on V-1100 due to the high voltage in the inside, even when it is disconnected from the mains.

### 3 OPERATING PRINCIPLE

V-1100 equipments are composed of five parts:

- An halogen lamp to supply the light.
- A monochromator to isolate the wavelength of interest and to eliminate unwanted second order radiations.
- A compartment for sample cuvettes.
- A photodetector to receive light transmittance and to convert it to electrical sign.
- A digital display where absorbance or transmittance.

Block diagram for the spectrophotometer.

### 4 INSTALLATION

After carefully unpacking the contents, check the materials with the packing list to ensure that you have received everything in good condition.

Keep the packaging for some days. If there is something missing, please contact your distributor as soon as possible.

Place the instrument in a stable level surface, suitable for the equipment weight and away from direct sunlight.

Set the unit up in an area that is free of dust, corrosive gases and strong vibrations.


Leave a free space of 30cm around the equipment for a good ventilation.

Turn on the spectrophotometer V-1100, it will begin to self test.

After a 20 minutes' pre-warm, you could take any readings.

## 5 OPERATION

### 5.1 Local Software Structure

### 5.2 Basic Operation

#### 1) Select Test Mode

Press **MODE** to select the measure mode.

#### 2) Set Wavelength

Turn the Wavelength Setting Knob to select the wavelength you want, and the wavelength value can be displayed on the screen in real-time.

#### 3) Input the Coefficient

When the system prompts you to input C, K or B, press

button to change the value till it displays the one you want, then press

to confirm.

Note: The system will always memorize your last input, until you input another value.

#### 4) Calibrate Zero

Close the lid of the compartment and press

to calibrate Zero T%.

#### 5) Calibrate 100%

Pull the Reference in the light path, press

to calibrate 100%TO.

#### 6) Print the result

Press the **ENTER** button to print the test result.


Fig. 4-1

## 6 MEASUREMENT

### 6.1 Measure the Absorbance

MODE

1) Press to choose «A» mode (Fig.4-1).

2) Turn the wavelength knob till it displays the wavelength value you want.

to set 100%T


Fig. 4-2

3) Pull the reference in the light path, press (Fig.4-2).

4) Pull the unknown concentration sample in the light path and then record the displayed value (Fig. 4-3).

Fig. 4-3

ENTER

5) Press to print the test result.

6) Repeat steps 4) and 5) to test other unknown concentration samples.

### 6.2 Measure the Transmittance

MODE

1) Press to choose the «T» mode (Fig. 4-4).

2) Turn the wavelength knob to set the wavelength at the point you want.

3) Pull the Reference in the light path to calibrate 100%T (Fig. 4-5).


Fig. 4-4

4) Pull the unknown concentration sample in the light path, then the value displayed on the screen is the one you need (Fig. 4-60).

ENTER

5) Press to print the test result.


Fig. 4-5

6) Repeat steps 4) and 5) to test other samples.


Fig. 4-6

WL.	546 nm
<input type="radio"/> A	Input Conc. : <b>1000</b>
<input type="radio"/> T	
<input checked="" type="radio"/> C	
<input type="radio"/> F	

Fig. 4-7

**6.3 Two-points Method**

This method allows measuring a concentration from two pattern concentrations.

1) Choose «A» or «T» mode and the wavelength.

2) Insert the «blank» cuvette, press

to get 100%T.

3) Insert the «pattern» cuvette and press

to choose «C» mode.

4) Press

or

to input the standard concentration, and press

to confirm. (Fig. 4-7).

5) Repeat to insert the second pattern.

6) Set the unknown concentration sample in the cell holder, then its concentration will display on the screen. (Fig.4-8).

7) Press

to print the test result.

8) Repeat step 6) and 7) to test other samples.

WL.	546 nm
<input type="radio"/> A	<b>985C</b>
<input type="radio"/> T	
<input checked="" type="radio"/> C	
<input type="radio"/> F	

Fig. 4-8

**6.4 Coefficient Method**

1) Set the wavelength.

2) Press

to choose «F» Mode.

3) Press

or

to set the coefficient value of «K» and «B», followed

with

ENTER
▲

pressed to confirm. (Fig. 4-9).  
From now on, you can test the samples.

4) Press

ENTER
▲

to set 100%T transmittance or zero absorbance (Fig. 4-10).

5) Set the samples in the cell holder. Its concentration will be displayed on the screen automatically (Fig.4-11).

6) Press

ENTER
▲

to print the test result.

WL.	546 nm
<input type="radio"/> A	INPUT K: <b>1000</b>
<input type="radio"/> T	
<input type="radio"/> C	
<input checked="" type="radio"/> F	

WL.	546 nm
<input type="radio"/> A	INPUT B: <b>1000</b>
<input type="radio"/> T	
<input type="radio"/> C	
<input checked="" type="radio"/> F	

WL.	546 nm
<input type="radio"/> A	BLANKING... PLEASE WAIT...
<input type="radio"/> T	
<input type="radio"/> C	
<input checked="" type="radio"/> F	

WL.	546 nm
<input type="radio"/> A	<b>1252</b>
<input type="radio"/> T	
<input type="radio"/> C	
<input checked="" type="radio"/> F	

## 7 SPECIFICATIONS

<u>Model</u>	<b>V-1100D</b>
<u>Wavelength Range</u>	325-1000nm
<u>Spectral Bandwidth</u>	4nm
<u>Optical System</u>	Single Beam, Grating 1200lines/mm
<u>Wavelength Accuracy</u>	±2nm
<u>Wavelength Repeatability</u>	1nm
<u>Wavelength Resolution</u>	0.5nm
<u>Photometric Accuracy</u>	±0.5%T
<u>Photometric Repeatability</u>	±0.3%T
<u>Photometric Range</u>	-0.3-3A, 0-200%T
<u>Stray Light</u>	0.3%T@360nm
<u>Stability</u>	±0.004A/h @500nm
<u>Display</u>	Graphic LCD (128X64 dpi)
<u>Keyboard</u>	4 Button keypad
<u>Photometric Mode</u>	T, A, C, F
<u>Detector</u>	Photodiode
<u>Sample Compartment</u>	Standard 10mm path length cuvette holder (with optional holder)
<u>Light Source</u>	Tungsten lamp
<u>Output</u>	USB Port Parallel Port (printer)
<u>Power Requirement</u>	85V ~ 265VAC
<u>Dimensions (W x D x H)</u>	480 x 360 x 160mm
<u>Weight</u>	10 kg

## 8 WARRANTY

This product is guaranteed for one year from the date of purchase and it covers defective components and assembly.

Please keep the purchase ticket as it is the only valid document for the purchase date.

This guarantee does not cover damage caused by an incorrect use of the equipment. This equipment must be properly validated before its use in the value sample analysis for the user.

Any return without the authorization of J.P. SELECTA, s.a.u. is not accepted.

## **9 TROUBLESHOOTING**

<b>PROBLEM</b>	<b>Possible Causes</b>	<b>Solution</b>
No functioning after power is on	<ul style="list-style-type: none"> <li>-Power cord not connected to outlet</li> <li>-Dead Power outlet</li> <li>-Internal fuse melted or defective electronic component</li> </ul>	<ul style="list-style-type: none"> <li>-Plug instrument in</li> <li>-Change to a different outlet</li> <li>-Change fuse or call an authorized service engineer.</li> </ul>
Instrument cannot set 100%T (0A)	<ul style="list-style-type: none"> <li>-Light beam blocked: Holder misaligned</li> <li>-Lamp is old or defective</li> <li>-Lamp is off alignment</li> <li>-Defective electronic component</li> </ul>	<ul style="list-style-type: none"> <li>-Check sample holder</li> <li>-Replace lamp</li> <li>-Readjust the lamp</li> <li>-Call an authorized service engineer</li> </ul>
Unstable display	<ul style="list-style-type: none"> <li>-Insufficient warm up time</li> <li>-Poor grounding</li> <li>-Excessive vibration, strong air current near light source, or strong external light.</li> <li>-Lamp old or defective</li> <li>-Unstable power supply</li> <li>-Defective or dirty detector or defective electronic component</li> </ul>	<ul style="list-style-type: none"> <li>-Warm up at least 20 minutes</li> <li>-Check ground connection</li> <li>-Improve working conditions</li> <li>-Replace with a new lamp</li> <li>-Equip a manostat to make the voltage stable</li> <li>-Call an authorized service engineer</li> </ul>
Incorrect readings obtained	<ul style="list-style-type: none"> <li>-Insufficient sample volume</li> <li>-Wrong wavelength setting Failed to blank (0A/100%T)</li> <li>-Stray sample preparation vapours</li> <li>-Bubbles or particles in solution</li> <li>-Instrument out of electronic calibration</li> </ul>	<ul style="list-style-type: none"> <li>-Fill cuvette with more samples</li> <li>-Check analytical procedure and wavelength setting.</li> <li>-Prepare sample away from instrument. Use proper ventilation</li> <li>-Check sample preparation and analytical procedure</li> <li>-Call an authorized service engineer</li> </ul>