pH-metro digital

Digital pH-meter

pH-2003

4120300

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
1001552

Probe holder

Upto 4 prober on the same holder.
Multiposition arm.

Code 1001552

7001638

Magnetic stirrer «Agimicro»

Regulable speed from 40 to 1700 r.p.m.
Enginereeng material chemical resistant.
Included, magnet teflon coated Ø8 x 32 mm.

Code 7001638

Buffer solutions for calibration

Buffer	pH=4.00 (20°C)	pH=4.01 (25°C)
Buffer	pH=7.00 (20°C)	pH=6.98 (25°C)
Buffer	pH=9.23 (20°C)	pH=9.18 (25°C)
Electrolite	KCl 3M	

Code	4120107
Code	4120108
Code	4120120
Code	4120109

Warning symbol

Important
information
symbol

2 Safety

Across this manual, some risk of danger situations are remarked. Proceed as described to avoid any hazard.

- **Warning symbol of a dangerous operation:** Proceed according indicated instructions.
- **Important information symbol:** Important information to obtain good results or for an optimal operation of the equipment or to extend the life of the equipment or to avoid the degradation of some of its components.
- When using a pH-meter, the operator could be exposed to chemical hazard.
- This equipment must only operated by personnel trained an suitable qualified to work in a chemical laboratory and its common hazards.
- This equipment must only operated by personnel who have read and understood this manual or have been qualified on its manipulation and operation.
- This equipment its not designed and not manufactured to work in an explosive ambient area.

3 Overview.

pH-2003 is an portable pH-meter suitable for measures on the field or in the laboratory.


It does not measure the temperature, but this can be entered for compensation.

The transport case, the electrode of resistant engineered material and its small size do it indicated for field measures.

It is delivered with a case to facilitate the transport of pH-meter and its accessories.

4 Technical specification

Power supply:	2 alkaline batteries 1.5 V size AA
pH range:	0..14
pH precision:	± 0.03
pH resolution:	0.01
mV range:	±1400
mV precision:	±2
mV resolution:	1
Temperature correction range:	0 ... 60°C (keypad entry)
Input impedance:	5 x 10 ¹¹ Ohm
Calibration points:	2
Ambient working conditions:	5 ... 35°C humity (Hr) lower than 80%
Dimensions:	170 x 75 x 30 mm
Weight:	0.5Kg

5 Description

All unit functions are selected from the front keypad. The display shows the pH reading, the temperature set and the pH-meter status.

1. BNC Electrode connector.
2. Display.
3. On/Off.
4. Mode push button
5. Enter push button
6. Up push button.
7. Down push button.
8. Back light display push button
9. Electrode holder

6 Operation

6.1 Placing batteries

- Place the batteries on the rear side holder.
- Check correct polarity.
- If the pH-meter will not be used for a long time, remove the batteries.
- When message «LoBat» appears on the display, change the batteries.

6.2 Electrode connection

The electrode must be stored with its immersion cap filled with electrolyte.

Transport the pH-meter in the briefcase with the electrode unplugged, and the pH-meter protected with the short-circuit connector.

To connect the electrode, remove the short-circuit connector from the pH-meter and to connect the electrode.

Note: The BNC connectors are extracted and placed by turning 1/4 turn.

Keep the electrode with its immersion cap until use.

CAUTION:

The immersion cap contents KCl 3M as electrolyte which is highly corrosive.

Take suitable cautions to avoid electrolyte to be in contact with skin or wear. Avoid any contact with eyes.

6.3 Keypad

Press Start/Stop to switch on the pH-meter.

Press «Mode» to select the different operations:

- «°C»: To enter the measuring temperature.
- «STD1»: Input standard 1.
- «STD2»: Input standard 2.
- «MEAS»: Selection between «pH» or «mV».

6.4 Changing mV to pH

pH-2003 allows measure on mV or pH.

On starting (On/Off), measure is in mV.

Press «Mode» and select «MEAS» to measure pH.

Measure in mV could be used to measure with an ion selective electrode.

6.5 Measure temperature

pH-2003 does not measures temperature. But it calculates a temperature compensation according to a temperature value entered from keypad.

In other words you can measure the sample temperature by a thermometer and then enter it on the pH-2003.

To enter temperature:

- Select «°C» by pressing on «Mode» (°C units blink on the display).
- Use keys «Up» or «Down» to select temperature.
- Validate with «ENTER»

Note: If the temperature set it quite different from the real temperature on sample, the pH measured could have an error.

6.6 Calibrating pH-2003

El pH-2003 must be calibrated on 2 points.

Choose this two points by considering:

- One point is always pH=6.98
- If you are going to measure an acid solution, choose pH=4.00
- If you are going to measure an alkali solution, choose pH=9.18

Calibrating pH-2003:

- Select «STD1» by pressing «Mode»
- Immerse the electrode cap on the pH=6.98 buffer.
- Wait until the mV reading is stable.
- Press «ENTER»

- Clean and rinse the electrode cap with clean water.

- Select «STD2» by pressing «Mode»
- Immerse the electrode cap on the pH=4.00 or 9.18 buffer.
- Wait until the mV reading is stabilized.
- Press «ENTER»

- Clean and rinse the electrode cap with clean water.

pH-meter is calibrated.

Note: Usually a pH-meter must be calibrated:

- First time use
- Daily.
- After changing the electrode
- After measuring strong acids (pH<2)
- After measuring strong alkalis (pH>12)
- After measuring very dense solutions that can obstruct the membrane.

15 Warranty

The warranty covers the defective components or faulty assembling of the unit for 1 year since the purchasing date.

Keep purchasing bill as the only valid document for the purchasing date.

The warranty does not covers damages caused for a unit malfunction. The unit must be suitably qualified before been used for valuable samples processing.

Return good are not accepted without previous J.P. SELECTA, s.a.u. approval.

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.

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