Conductímetro portatil

Portable conductivity meter

CD-2004

4120410

Indice

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David Pecanins
Rble. Calidad

1 Contents list and accessories

First, check that the contents is complete and is in a good state.
Keep the packaging for some days.

1.1 Standard contents

- Conductivity-meter (1)
- Battery 1.5 AA (4 units)
- RS-232 communication cable (6)
- Transport briefcase
- Conductivity sensor K=1
- Inox temperature sensor Rt2252
- Instruction manual

1.2 Accessories

Conductivity probes (glass)

Fig.	Code	K	Temp. probe	Size	Uses
1	4120220	1	No	Ø12 x 120	General use on the lab up to 200uS
2	4120221	5	Yes	Ø12 x 120	Up to 20 ms
3	4120222	0.1	No	Ø12 x 120	Pure water up to 20uS
4	4120223	10	No	Ø12 x 120	High conductivity up to 200mS

Universal arm support

Code 1001552

With a device for 4 electrodes or temperature probes. Metallic arm support of high stability. Especial for working at different positions and heights.

1001552

5504300

Magnetic stirrer «agimicro»

Code 5504300

Features:

Regulable speed from 0 to 1400 r.p.m.
Excellent resistant to chemical agents.
Rear hole for support bar (accessory)
Magnet PTFE coated Ø8 x 32 mm included.

Certified conductivity standards (250ml)

Features:

±0.05	1278uS/cm	(20°C)	1413uS/cm	(25°C)	Code	4120160
±0.025	4.915mS/cm	(20°C)	5.446mS/cm	(25°C)	Code	4120161
±0.06	11.67uS/cm	(20°C)	12.88uS/cm	(25°C)	Code	4120162

4120113

Serial printer

Code 4120113

Serial printer with ink cartridge (no thermal).
Direct connection to CD-2004
Paper width: 52mm

2 Safety

Risky situations in the conductivity meter are indicated with the following icons:

- **Warning symbol of a dangerous operation:** Proceed according the indicated instructions.
- **Important information symbol:** Important information to obtain good results and an optimal operation of the equipment.
- The use of the conductivity meter expose the operator to dangerous chemical hazard.
- This equipment must only be operated by staff qualified to work in a safety way in a chemical laboratory.
- This equipment is not designed and not manufactured to work in an explosive ambient area.

Warning symbol

Important information symbol

• Standard solutions

Attention

- Conductivity standard solutions are irritating chemical products. Avoid any contact with eyes and use a lab coat to protect your wear.
- In case of skin contact, please quickly clean with abundant water.
- Once open, please keep the solutions in a closed flask. Do not use the solutions after 1 month the flask is open.



3 Overview

The CD-2004 pH-meter is a portable conductivity meter protected from water splashing. It's suitable for outdoor measurements (lakes, rivers,..) thanks to its capability of data storing. You will be able to measure simultaneously conductivity and temperature.

The CD-2004 also measures temperature but if there's no temperature sensor, sample temperature can be entered by using the keypad and also there's an automatic temperature compensation applied.

4 Technical specifications

Power supply:	4x 1.5V AA size. alkaline batteries
Conductivity range:	0 ... 1.999x10⁵ μS/cm in 5 ranges:
	0.00 ... 19.99 μS/cm
	20.00 ... 199.9 μS/cm
	200.0 ... 1999 μS/cm
	2.00 ... 19.99 mS/cm
	(Constant sensor = 5 or 10) 20.00 ... 199.9 mS/cm
TDS range:	0 ... 19990mg/L in 5 ranges::
	0.00 ... 10.00 mg/L
	10.00 ... 100.0 mg/L
	100.0 ... 1000 mg/L
	1.00 ... 10.00 g/L
	(Constant sensor = 5 or 10) 10.00 ... 19.99 g/L
Salinity range:	0.00 ... 8.00 %
Temperature range:	0.0 ... 40.0°C
Temperature measure accuracy:	±0.6°C
Conductivity measure accuracy:	±1.6% (FE)
Salinity reading accuracy:	± 0.35% (FE)
Working conditions:	5 ... 40°C humidity (Hr) lower than 85%
Dimensions:	201 x 100 x 45 mm
Weight:	0.5Kg
Protection degree:	IP65
Temperature sensor:	Rt 2252

4.1 Factory set parameters

Parameter:	Valid range	Factory set
Sensor constant:	0.01; 0.1; 1; 5; 10	1.00
ATC factor:	0.000 ... 0.040	0.020
TDS factor:	0.40 ... 1.00	0.50
ASD automatic stop:	ON/OFF	ON

5 Description

5.1 Display

- Conductivity reading/TDS/Salinity
- Units
- Actual measurement temperature
- Status indicator

5.2 Rear panel

- Conductivity sensor connector.
- Mini BNC connector for temperature sensor connection
- Serial RS-232 connector

6 Operation

Once the equipment is unpacked, please visually check that all goods are OK. Keep the packaging for some days.

6.1 Placing batteries

- Place the batteries on the case, in the rear side of the pH-meter.
- Pay attention polarity is correct.
- If you are not going to use the unit for some time, please remove the batteries.
- Replace the batteries when «LoBat» appears on the display.

6.2 Sensors connection

Connect the conductivity and the temperature sensor on the rear panel.

Do not move the unit with the sensors connected.

Note: The BNC connectors must be removed and fixed by turning 1/4 turn.

6.3 Panel

Switch on the pH-meter by pressing the «ON/OFF» key.

Wait until the equipment replies before pressing any other key.

Keys:

- «ENTER/PRINT»: Validate input. Prints the reading.
- «MODE /MEAS»: Changes from MEASuring status to MODE status and vice versa.
- «▲ /C/T/S»: Up. Change from conductivity, TDS or salinity.
- «▼ / SAVE»: Down. Stores the reading.

6.4 Measurement mode and configuration mode

The equipment can be configured in two different modes:

«Measurement mode»: when the equipment measures conductivity, TDS or salinity.

«Configuration mode»: when the equipment is used for configuration or calibration.

There are some different «states» inside each «mode»:

6.4.1 «Measurement mode» states:

DISPLAY indication	Meaning
CON	When measuring conductivity
TDS	When measuring TDS
SALT	When measuring salinity
COM	When transferring data RS-232
PRN	When sending data to the printer
SAVE	When saving data to memory
OVER	When it's out of range

6.4.2 «Configuration mode» states:

Indicators during **conductivity** configuration.

DISPLAY indication	Meaning
CONT/ADJS	Sensor constant selection and adjustment
COEF	Temperature compensation coefficient
CAL/CA-M	Sensor constant calibration

Indicators during **TDS** configuration.

DISPLAY indication	Meaning
CONT/ADJS	Sensor constant selection and adjustment
COEF	Temperature compensation coefficient
CAL/CA-M	Sensor constant calibration
TCAL/CA-M	Transfer coefficient selection and adjustment
TCOE	Transfer coefficient calibration

Indicators during **Salinity** configuration.

DISPLAY indication	Meaning
CONT/ADJS	Sensor constant selection and adjustment
CAL/CA-M	Sensor constant calibration

Common parameters to all measurement modes.

DISPLAY indication	Meaning
PRN/ SURE/ NULL/ PRN-	Data printing dialog
VIEW/ N001/ NULL/	Data visualization dialog
DEL/ SURE/ DEL-	Data erase dialog
ASD/ ON/ OFF	Equipment automatic stop configuration

6.5 Conductivity measurement initial display


When switching on the CD-2004 («ON/OFF» key) the conductivity and temperature measurement initial display appears.

Push «▲/C/T/S» key to change measurements in conductivity, TDS and salinity.


Note: In order the measurements are precise, it's necessary to configure and calibrate the conductivity meter.

Conductivity measurement display

TDS measurement display

Salinity measurement display

6.5.1 Conductivity measurement: (CON)

In order this measurement is correct, you must use the configuration mode to set up the following parameters:

- Sensor constant.
- Temperature compensation coefficient.
- Conductivity sensor constant calibration.

6.5.2 Total Diluted Solids measurement: (TDS)

In order this measurement is correct, you must use the configuration mode to set up the following parameters:

- Sensor constant.
- Temperature compensation coefficient.
- «TDS transferring coefficient»
- Conductivity sensor constant calibration
- Coefficient calibration: «TDS transferring coefficient»

6.5.3 Salinity measurement: (SALT)

In order this measurement is correct, you must use the configuration mode to set up the following parameters:

- Sensor constant.
- Conductivity sensor constant calibration

It's not necessary to insert a temperature compensation coefficient, as the equipment performs this compensation automatically.

6.6 Suitable conductivity sensor selection

Sensor constant to measure in a conductivity range:

Conductivity	Resistivity (Ohm)	Constant (cm ⁻¹)
0.05 ... 2 µS/cm	20M ... 500K	0.01 ; 0.1
2 ... 200 µS/cm	500K ... 5K	0.1 ; 1
0.2 ... 2 mS/cm	5K ... 500	1
2 ... 20 mS/cm	500 ... 50	1 ; 5 ; 10
20 ... 200 mS/cm	50 ... 5	5 ; 10

Sensor constant to measure in a TDS range:

TDS	Conductivity	Constant (cm ⁻¹)
0 ... 1000 mg/L	0 ... 2 mS/cm	1
1000 ... 10000 mg/L	2 ... 20 mS/cm	1 ; 10
1000 ... 19990 mg/L	20 ... 40 mS/cm	10

Sensor constant to measure in a salinity range:

A 10 (cm⁻¹) constant sensor is normally used to measure salinity; if the salinity is lower than 1%, we can also use a 10 (cm⁻¹) constant sensor.

6.7 Parameter configuration

In order the measurements are precise, it's necessary to configure and calibrate the conductivity meter.

6.7.1 Sensor constant configuration.

Every conductivity sensor is marked by its «constant» value. (For instance 0.998).

The sensor «constant» must be entered in the conductivity meter as for conductivity, TDS or salinity measurement.

The sensor constant configuration has two phases:

- Enter the «constant» nominal value: (0.01; 0.1; 1; 5; 10)
(By following the example, the nominal value is 1)
- Adjust this constant to the sensor marked value.
(By following the example, the marked value is 0.998)

The procedure is as follows:

- Push «MODE». The equipment changes to configuration mode.
- Push «▲ /C/T/S» or «▼ /Save» till «CONT» appears (on the lower left corner).
- Push «ENTER». The constant current value will appear.
- Push «▲ /C/T/S» or «▼ /Save» to select the nominal value.
- Push «ENTER» to accept the start. (ADJS will appear)
- Push «▲ /C/T/S» or «▼ /Save» to select the exact value of the sensor marked constant.
- Push «ENTER» to accept the entry. (CONT will appear)

6.7.2 TDS transfer coefficient configuration.

This coefficient is only needed to measure TDS .

It's not necessary to measure conductivity nor salinity.

The procedure is as follows:

- With the equipment set up to measure TDS. (TDS on the lower left corner)
- Push «MODE». The equipment changes to configuration mode.
- Push «▲ /C/T/S» or «▼ /Save» till «TCOE» appears (on the lower left corner).
- Push «ENTER». The coefficient current value will appear.
- Push «▲ /C/T/S» or «▼ /Save» to select the value.
- Push «ENTER» to accept the entry. (TCOE will appear)

6.7.3 Temperature coefficient configuration.

This is the coefficient which has to be entered in the conductivity meter to measure conductivity and TDS.

The procedure is as follows:


- With the equipment set up to measure conductivity or TDS.
- Push «MODE». The equipment changes to configuration mode.
- Push «▲ /C/T/S» or «▼ /Save» till «COEF» appears (on the lower left corner).
- Push «ENTER». The coefficient current value will appear.
- Push «▲ /C/T/S» or «▼ /Save» to select the value.
- Push «ENTER» to accept the entry. (COEF will appear)


Nominal value inlet
of sensor constant

Exact value inlet
of sensor constant

Exact value inlet
of sensor TDS
transfer coefficient

Temperature
coefficient inlet

6.8 Calibration.

In order the measurements are precise, it's necessary to calibrate the conductivity meter.

6.8.1 Sensor constant calibration.

In the course of time and use, the sensor «constant» as well as the TDS transfer factor change their value. So they need a regular calibration.


6.8.1.1 Conductivity calibration standards selection.

According to the sensor constant value, select a calibration standards in accordance with the following table:

Sensor constant:	Advised standards conductivity
0.01 ... 0.1	100 µS/cm to 25°C
1	1413 µS/cm to 25°C
5 ... 10	12.88 mS/cm to 25°C

The procedure is as follows:

- Connect the conductivity meter sensor. Temperature sensor **switch off**.
- Standards must be at 25°C.
- Clean the sensor with distilled water.
- Immerse the sensor in the standards.
- If the sensor constant has not entered, see 6.7.1.
- Push «▲ /C/T/S» or «▼ /Save» till «CAL» appears (on the lower left corner).
- Push «ENTER». The absolute conductivity reading will appear (without temperature compensation).
- Wait till the reading will be stable.
- Push «ENTER» and «CA-M» will appear.
- Push «▲ /C/T/S» or «▼ /Save» to adjust the reading to the exact value of the conductivity standards.
- Push «ENTER» to accept the entry. The equipment will re-calculate the constant value to adjust it to the standards.
- The new constant value is shown in the display for 5 seconds. It comes out automatically from the calibration menu.


Automatic calibration of sensor constant

6.8.2 TDS transfer coefficient calibration.

In the course of time and use, the TDS transfer factor change its value. So it needs a regular calibration.

The procedure is as follows:

- Connect the sensor to the conductivity meter. Temperature sensor **switch off**. Set up for TDS measurements.
- Standards must be at 25°C.
- Clean the sensor with distilled water.
- Immerse the sensor in the standards.
- If the sensor constant is not entered, see 6.7.2.
- Push «▲ /C/T/S» or «▼ /Save» till «TCAL» appear (on the lower left corner).
- Push «ENTER». The absolute TDS reading will appear (without temperature compensation).
- Wait till the reading will be stable.
- Push «ENTER». «CA-M» will appear.
- Push «▲ /C/T/S» or «▼ /Save» to adjust the reading to the exact value of the TDS standards.
- Push «ENTER» to accept the entry. The equipment will re-calculate the TDS transfer coefficient value to adjust it to the standards.
- The new constant value is shown in the display for 5 seconds. It comes out automatically from the calibration menu.

Automatic calibration of sensor transfer coefficient

6.9 Save readings.

Push «▼ / SAVE» directly from the measurement display to store the reading. (Conductivity, TDS or salinity)

Conductivity, TDS and salinity measures are saved, printed and deleted separately.

Up to 250 conductivity readings, 250 TDS readings and 250 salinity readings can be stored (total 750). If you store more than this quantity, the new ones overwrite the first one stored.

6.10 Delete readings.

Deleting the stored readings:

Note: All the readings of the selected measure will be deleted (conductivity, TDS or salinity). It's not possible to delete just one reading.

- First select conductivity, TDS or salinity. They will be deleted in a separate way.
- Push «MODE/MEAS».
- Push «▼ / SAVE» till «DEL» appears in the display.
- Push «ENTER/PRINT» to confirm.

6.11 Print readings.

The readings could be printed if the CD-2004 has a printer connected (RS-232 serie connection, Baud rate 9600).

For an instant print of the current reading: Push «ENTER/ PRINT».

For a print of all the stored data: (see 6.9).

- First select conductivity, TDS or salinity. They are separately printed.
- Push «MODE/MEAS».
- Push «▼ / SAVE» till the «PRN» display.
- Push «ENTER/PRINT» to confirm.

6.12 Visualize readings.


Visualize the stored readings:

- First select conductivity, TDS or salinity. As they are separately visualized.
- Push «MODE/MEAS».
- Push «▼ / SAVE» till the «VIEW» display.
- Push «ENTER/PRINT» to confirm.

6.13 Automatic switch off.

Configure the CD-2004 to automatically switch off after 2 hours of no use.


- Push «MODE/MEAS».
- By pushing «ENTER/PRINT».
- Select with «▼»: and «▲»
ON : Automatic switch off.
OFF: Without automatic switch off.


Data deletion

Print stored readings

View stored readings

Configuration display of automatic switch off

6.14 RS-232 Communication

Use the communications software to register measurements on a PC. Connect the communications cable to a free COM port. Speed 9600.

6.15 Conductivity sensor manipulation.

- Follow the own sensor instructions, if there is some.
- Keep the electrode with its protection.
- Any hit to the sensor could modify the sensor geometry and so its constant.
- After measuring concentrated or aggressive solutions, please clean the sensor with distilled water.
- After any measurement, clean thoroughly to avoid interferences.

7 Warranty

The warranty covers the defective components and the assembly faults for 1 year from the purchase date.

Please keep the purchase invoice as the unique valid document of the purchase date.

The warranty does not cover the damage due to a defective operation of the equipment. The equipment must be validated in a suitable way before its use in the valuable samples analysis for the users.

No goods return are accepted without the previous authorization of J.P. SELECTA, s.a.

8 CE Declaration of conformity

The manufacturer: J.P. SELECTA, s.a.
Ctra. NII Km 585.1
08760 ABRERA (BARCELONA)
SPAIN

Declares that the equipment: CD 2004 Code: **4120410**

Complies with the following european directives :

73/23/CEE	Electrical safety directive.
89/336/CEE	Electromagnetic compatibility directive.


Complies with the following standards:

EN 61326:2002
Electrical material for measurement, control and laboratory use. Electromagnetic compatibility requirements (CEM).

EN 61010-1:2002
Safety requirements of electrical equipments of measurement, control and laboratory use. Part 1: General regulations.

Abrera, May 2007

 Ramon Ramon
 Technical Dtor.

 David Pecchini
 Quality Resp.