pH-metro digital

Digital pH-meter

pH-2005


4120500

Indice

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Index

1 Packing list and accessories	12
1.1 Packing list.....	12
1.2 Accessories	12
2 Safety	13
3 Overview	13
4 Technical specification	13
5 Description	14
5.1 General.....	14
5.2 Rear part	14
5.3 Keyboard	14
6 Operation	15
6.1 Location and connection to th power supply	15
6.2 Probe suppor holder.....	15
6.3 Plugging electrode and temperature probe	15
6.4 Measures	16
6.5 Calibration.....	16
6.5.1 Calibration at one point	17
6.5.2 Calibrations at two points	17
6.6 Potential measurement	17
6.7 Preparing buffer solutions.....	17
6.8 Electrode handling.....	18
7 Warranty	18

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
ENGLISH VERSION

1 Packing list and accessories

Once unpacked the equipment, check visually that everything is in good condition and the contents are correct.

1.1 Packing list

- 1 pH-meter
 - 2 Probe holder
 - 3 Short-circuited connector
 - 4 Dust keeper
 - 5 USB communication wire
 - 6 Power cord
 - 7 Buffer solution powder
- Instruction manual

1.2 Accessories

pH combined electrodes:

Fig.	Code	Range pH	Range (°C)	Temp. Probe	Electrolite	Size	Material Body	Application
1	4120102	0 - 14	0-80°	No	Liquid	Ø12 x 120	Glass	Laboratory general use.
2	4120105	0 - 14	0-60°	Si	Gel	Ø12 x 120	Epoxi	Anti-chock body. Laboratory general use and outdoors.
3	4120101	0 - 14	0-60°	No	Gel	Ø12 x 120	Epoxi	Anti-chock body. Laboratory general use and outdoors.
4	4120114	0 - 10	0-60°	No	Gel	Ø12 x 120	Epoxi	Suface sensor. Skin, meat, wet semisolids, culture.
5	4120104	0 - 14	0-60°	No	Gel	Ø12 x 120	Epoxi	Sharper sensor. Semisolid food, fruits, meat.
6	4120103	0 - 14	0-80°	No	Gel	Ø8x 60	Glass	For in tube measures. Aqueous an solutions.
7	4120106	0 - 14	0-80°	Si	Liquid	Ø12 x 120	Glass	Fast response on aqueous means.

Probe holder

Up to 4 probe on the same holder.
Multiposition arm.

Code 1001552

Certified buffer solutions for calibration

Buffer	pH=4.00 (20°C)	pH=4.01 (25°C)	Code	4120107
Buffer	pH=7.00 (20°C)	pH=6.98 (25°C)	Code	4120108
Buffer	pH=9.23 (20°C)	pH=9.18 (25°C)	Code	4120120
Electrolite	KCl 3M		Code	4120109

Warning symbol

Important
Information
symbol

2 Safety

Across this manual, some risk of danger situations are remarked. Proceed as described to avoid any hazard.

- **Warning symbol of a dangerous operation:** Proceed according indicated instructions.
- **Important Information symbol:** Information important to obtain good results or for an optimal operation of the equipment or to extend the life of the equipment or to avoid the degradation of some of its components.
- Using a pH-meter, the operator, could be exposed to chemical hazard.
- This equipment must only be operated by personnel trained and suitable qualified to work in a chemical laboratory and its common hazards and who have read and understood this manual or have been qualified on its manipulation and operation.
- This equipment its not designed and not manufactured to work in an explosive ambient area.

3 Overview

The pH-2005 is a table-top meter with a big character display for easy reading. Suitable for indoor measures on the laboratory. It allows to read pH and temperature simultaneously.

The pH-2005 can measure temperature but if there is any temperature probe connected, sample temperature can be entered through the keyboard in order it makes the automatic temperature compensation. Function of automatic buffer solution recognition for 4.00, 6.86 and 91.18. Also can be done measures in mV with electrodes of selective ion.

4 Technical specification

Power supply:	See feature label
pH range:	0.14 (Display -2...18.00)
pH precision:	± 0.02 pH
pH resolution:	0.01 pH
mV range:	±1999 mV
mV precision:	±1 mV
mV resolution:	1 mV
Temperature correction range:	0 ... 99.9 °C
Temperature precision:	+0.5°C
Temperature resolution:	0.1°C
Temperature compensation precision:	> 0.01 pH ± 1 bit
Input impedance:	> 1 x 10 ¹² Ohm
Calibration points:	2
Room working conditions:	5 ... 40°C humidity (Hr) lower than 85%
Dimensions:	290 x 210 x 95 mm
Weight:	1.5Kg

5 Description


5.1 General


1. Keyboard.
2. Display.
3. Mooring rod for probe holder

5.2 Rear part

4. pH or mV electrode socket
5. Reference electrode socket
6. Temperature probe socket
7. USB port
8. Fuse holder
9. Main switch
10. Power supply socket

5.3 Keyboard

11. pH or mV selection key
Leave from the calibration mode
12. Manual temperature set key
13. Manual second calibration point correction (SLOPE)
14. Manual first calibration point (STD)
15. «ENTER» key

6 Operation

6.1 Location and connection to the power supply

Place the pH-meter on a flat surface, clean of dust, not to direct sun light exposition, and near to a mains plug. For an ergonomic use, let 30cm free space round the pH-meter.

Use the power cord to connect pH-meter to mains plug provided with an earth conductor.

6.2 Probe support holder

- Pull the mooring rod platform (3) until the needed size
- Insert the probe holder on the rod mooring
- Fix holder with the nut
- Adjust the position tighten its several joints

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To connect the electrode, the short-circuiting plug must be extracted from the socket (4) and then, connect the electrode in the same socket.

It is advised to keep the short-circuited plug in the socket (4) while it is not in use.

The electrode must be stored with its end preserved by its protective cap previously filled with electrolytic solution.

Note: The BNC connectors have to be extracted and fixed doing 1/4 of turn.

Keep the protective cap on the electrode until just before to use it.

ATTENTION:

The protective solution of the electrode contain potasic chloride 3M.
Be careful: potasic chloride 3M can cause skin irritations and or cloth damages.
Avoid contact with eyes.

6.4 Measures

To do measures follow the steps below:

- Make sure to have connected the device to a suitable power supply
- Make sure to have installed the required electrode (pH o mV) in the socket (4)
- If it is available, make sure to have the electrode connected to the socket (6)
- Switch the main switch (9) on
- Select pH or mV with the key (11)
- Remove the protective caps from the electrode
- Rinse the electrode in distilled water
- Introduce the electrode in the liquid to measure

If you work with temperature probe it must keep introduced inside the liquid

In order to improve the accuracy it is advised to do a preheating 30 minutes.

Normally this device does not need the manual introduction of temperature but, if it is needed, temperature can be measured with an external thermometer and then enter it through the keyboard.

It is advisable to have a probe connected to the equipment in order to read the temperature of the sample.

6.5 Calibration

It is advisable to calibrate the device before to do the measurements or, at least, once a day in case of continuous working.

The device recognizes automatically three different kinds of buffer solutions: pH 4.00, pH 6.86 and pH 9.18. In these cases is not necessary to do any adjustment after press the key «STD» or «Slope» simply press the key «ENTER» to end the calibration.

Also can be done calibrations with other kind of buffer solutions or doing calibrations at two points.

Before to introduce the electrode in the buffer solution it must be rinsed with distilled water.

After doing the calibration in one point the old calibration date will be erased automatically.

It is advisable to know the value of the buffer solution at several temperatures.

6.5.1 Calibration at one point

The sequence is the following:

- Rinse the electrode with distilled water
- Introduce the electrode in the buffer solution
- Measure the liquid solution with a thermometer. Press the key «TEMP» and correct it if necessary
- Press the key «STD», the display will show «Std, Yes», then press the key «ENTER» in order the device recognizes automatically solutions pH 4.00, pH 6.86 or pH 9.18, in case of other kind of buffer correct manually
- Press the key «ENTER» to end the calibration at one point.


6.5.2 Calibration at two points

The sequence is the following:

- Follow the instructions from the section «calibration at one point»
- Rinse the electrode once again and introduce it in the second buffer solution.
- Measure the temperature of the liquid and press the key «TEMP» until the display show the current temperature.
- When the reading be stable pres the key «SLOPE», the display will show «Std YES» then press «ENTER» so that the display the current pH and temperature.
- Press the key «ENTER» to end the calibration at two points

If the device has a temperature probe connected, it must keep introduced in the liquid in order the equipment adquire automatically its current value. In this case is not necessary to introduce manually the temperature.

6.6 Potential measurement (mV)

- Place the selective ion electrode and the reference electrode on the support holder.
- Rinse the electrode first with distilled water and then with sample solution.
- Connect the selective ion electrode to the plug (4)
- Connect the reference electrode to the plug (6)
- Introduce both electrodes in the sample solution and stirr it. Then the display will show the potential (mV) of this selective ion electrode with its belonging polarity (+/-).
- If the potential of the sample solution is out of range, the display will show «Err»

6.7 Preparing buffer solutions

Together to the equipment several powder satchets are supplied in order to prepare buffer solutions:

- pH4.00. 10.12 gr of powder to solve in 1000 ml of bidistilled.
- pH6.86. 3.387 gr of powder to solve in 1000 ml of bidistilled
- pH9.18. 3.80 gr of powder to solve in 1000 ml of bidistilled.

6.8 Electrode handling

- Follow the instructions which comes with the electrode if there are.
- Keep the electrode with its end cap filled with electrolytic liquid KCl 3M when it is not used.
- In order to prevent splashes on skin or clothes preventions must be taken. KCl 3M can produce skin irritations and attacks clothes. Avoid eyes contact.
- Avoid than the electrode left more than 1h its end in the air. If it becomes dry it can lose its features. It is advised to keep it immersed in distilled water.
- Rinse with plenty of distilled water after each measurement.

7 Warranty

The warranty covers the defective components or faulty assembling of the unit for 1 year since the purchasing date.

Keep purchasing invoice as the only valid document for the purchasing date.

The warranty does not covers the damages caused by an inappropriate use. The device must be suitable validated before to be used to analyse valuable samples.

J.P. SELECTA, s.a. will not accept any return of goods unless previous approval.

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