Conductímetro digital de sobremesa

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
English version

1 Packing list and accessories

Once unpacked the equipment, visually check that everything is in good condition and its contents is correct. Keep the packaging for some days.

1.1 Standard contents of the packaging

- Conduct-meter (1)
- Probe holder (2)
- Dust keeper (3)
- Instruction manual (4)
- USB communication wire (5)
- Power supply cord (6)
- Fuses (7)

1.2 Accessories

Conductivity probe (glass):

Fig.	Code	K	Temp. probe	Size	Uses
1	4120220	1	Yes	Ø12 x 120	General use on the lab up to 200µS
2	4120221	5	Yes	Ø12 x 120	up 20 ms
3	4120222	0.1	Yes	Ø12 x 120	Pure water 20µS
4	4120223	10	Yes	Ø12 x 120	High conductivity up to 200mS

Probe holder

Cat. number 1001552

It is able to hold up to 4 electrodes or temperature probes. High stability because of its metallic base. Suitable to work at several heights and positions.

Magnetic stirrer «Agimicro»

Cat. number 5504300

Features:

Adjustable speed from 0 to 1400 rpm
 Excellent resistance to chemical reagents
 Back hole to place a holder bar (optional accessory)
 It is supplied with teflon coated. Dimensions Ø8 x 32 mm.

1001552

5504300

Certified conductivity patterns (250ml)


Features:

±0.05 1278 µS/cm (20°C) 1413 µS/cm (25°C) Cat. number **4120160**
 ±0.025 4.915 mS/cm (20°C) 5.446 mS/cm (25°C) Cat. number **4120161**
 ±0.06 11.67 µS/cm (20°C) 12.88 µS/cm (25°C) Cat. number **4120162**

2 Safety

Warning symbol

Important information symbol



- **Warning symbol of a danger operation:** Proceed according the indicated instructions.
- **Important information symbol:** Important information to obtain good results or an optimal operation of the equipment or to extend its life or to avoid the degradation of some of its parts.
- Tests and determinations reported in this manual can expose the user to a chemical reagents contact.
- This device have to be used by personnel properly qualified to work in a chemical laboratory safely.
- Not suitable to work in flammable or explosives atmospheres.

3 Overview

The CD-2005 is a table top conductivity meter with a big display in order to make easier the reading. Suitable to measures in laboratories. It allows to read conductivity and temperature simultaneously. Furthermore it makes an automatic temperature compensation. The temperature reading can be visualised into η S/cm or into salinity units TDS (mg/L).

4 Technical specifications

Power	230V 50/60 Hz
Conductivity range	0...20.00 μ S/cm 0...200.0 μ S/cm 0...2000 μ S/cm 0...10.00 mS/cm
TDS range	0...1000 mg/L
Reading precision	± 1.5 end of scale + 1 digit
Temperature correction range	0...60 $^{\circ}$ C
Temperature precision	± 0.5 $^{\circ}$ C
Temperature compensation factor	2%
Calibration points	1
Room working conditions	5...40 $^{\circ}$ C humidity (Hr) lower than 85%
Dimensions	300x200x75 cm
Weight	0.5 kg

5 Description

5.1 General

1. Keyboard
2. Display
3. Probe holder

5.2 Display



4. Conductivity or TDS reading
5. Units (μ S/cm or mg/L)
6. Temperature measure
7. Status indicator.

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6 Operation

6.1 Locate and connection to the power supply

Place the conductivity meter on a surface even, plane, free of dust, out of the direct sunlight and near to a power supply plug. To an ergonomic use, left about 30 cm free from objects around.

Before to connect the device to the mains check that the voltage is the indicated on the feature label attached on the device. Use the supplied connection cord plugged to a socket properly hearted.

Keep in mind:

- The **probe constant** is a probe parameter that can be:

0.01, 0.1, 1 or 10. If the probe constant is specified with an exact number, enter the nearest value. (Example: 0.97 select constant = 1)

- The **Correction factor** is calculated every time that the unit must be calibrated.

Note:

If the temperature entered manually is far from the real temperature, the conductivity reading could have a big error.

6.2 Connection of the electrode and temperature probe

Temperature probe: Connector BNC mini: It has to be extracted and locked by turning 1/4 of turn.

Connector of conductivity probe: Hold the el connector by its rear part, press and turn.

Disconnect: Pull back from the front of the connector.

6.3 Control panel

The device starts up by the main switch on the rear part (18).

All functions of the equipment can be selected by the keyboard:

- Key "COND/TDS" (8): Pressing the key «COND/TDS» can be selected the reading mode "ηS/cm" or salinity "mg/L".
- Keys "CONS" (9) (probe constant): With the keys "▼" and "▲" this value varies: 0.01, 0.1, 1 and 10.
- Keys "ADJ" (10) (adjust): With the keys "▼" and "▲" this value varies in order to adjust the value and makes it meet with a pattern.
- Keys "TEMP" (11) (temperature): With the keys "▼" and "▲" the working temperature can be introduced manually when there are not any temperature probe connected.
- Key "ENTER" (12): Key to confirm values.
 - °C: To manually enter the temperature if there's no probe connected. Use the keys «▼» & «▲» & ENTER.
 - **Probe constant and correction factor**
Use the keys: ENTER, «▼» & «▲» .

6.4 Temperature constant

The temperature constant is a feature from the conductivity probe given by the manufacturer. It can be 0.01, 0.1, 1 or 10. If the probe specifies a different value must be selected the nearest one, for example, if the values is 0.97 se must be selected 1.

With the keys "▼" and "▲" (8) can be modified the values from: 0.01, 0.1, 1 and 10 and validate by pressing the key "ENTER" (12).

6.5 Calibration

The CD-2005 can be calibrated in one point. It is necessary to use a calibration pattern near the reading point.

Calibration of the CD-2005:

- Having conductivity pattern with a determined value (P)
- Check that the conductivity constant of the probe is correct: 0.01, 0.1, 1 or 10
- Introduce the conductivity cell inside the pattern liquid
- Wait that the reading will be stable
- Make note of the reading (L)
- Calculate the correction factor: $C = P / L$
- Press one of the keys «ADJ» (10), the display will show the current value
- With the keys «▼» and «▲» enter the calculated value "C"
- Press «ENTER» (12) to validate the new correction value

The conductivity meter is calibrated.

Note: Usually the conductivity meter must be calibrated:

- On starting up
- When the probe is changed

6.6 Conductivity cell handling

- Follow the own instructions from the conductivity cell if it has.
- Avoid to hit the cell.
- Clean with water after every measurement.

7 Warranty

The warranty covers the defective components or faulty assembling of the unit for 1 year since the purchasing date.

Keep the purchasing invoice as the only valid document for the purchasing date.

The warranty does not covers the damages caused for a unit malfunction. The unit must be suitably qualified before to be used for valuable samples processing.

Return good are not accepted without previous J.P. SELECTA, s.a.u. approval.

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