# MEDIDOR MULTIPARAMÉTRICO MULTI-PARAMETER ANALIZER

**MP-2006** 

4120700

 $\epsilon$ 

# Contenido

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# **Content**

- 1. General
- 2. Main Specification
- 3. Construction
- 4. Operation
- **5.** Maintenance
- 6. Complete set of the Meter
- 7. Appendix

#### 1. General

MP-2006 multi-parameter analyzer (hereafter refers to meter) is a new laboratory analyzer which contains ion measuring part, conductive measuring part, dissolved oxygen measuring part and temperature measuring part. It allows inspecting the relevant parameters of above parts. It also allows users to select individual parameter for measuring according to the practical needs. In which ion measuring part allows to measure: potential value, pH value (or pX value) and ion concentration. Conductive measuring part can measure: conductivity, resistivity, TDS and salinity. Dissolved oxygen measuring part can measure: electrode current, concentration of dissolved oxygen and saturation of dissolved oxygen. Temperature measuring part can measure temperature value of present solution.

#### The meter has following character:

#### A. Ion Measuring Part

- a. To measure Potential value, pH value (or pX value) and ion Concentration.
- b. To measure several general ions. Meter provides several general ion modes: H<sup>+</sup>, Ag<sup>+</sup>, Na<sup>+</sup>, K<sup>+</sup>, NH<sub>4</sub><sup>+</sup>, Cl<sup>-</sup>, F<sup>-</sup>, NO<sub>3</sub><sup>-</sup>, BF<sub>4</sub><sup>-</sup>, CN<sup>-</sup>, Cu<sup>2+</sup>, Pb<sup>2+</sup>, Ca<sup>2+</sup>, etc...for easy use. User can match relevant ion selective electrode and reference electrode to measure corresponding ion concentration directly. When measuring finishes, use can easily read the concentration between different concentration units.
- c. Apart from the ion part provided by the meter, in case you need to measure other ion, you can also set your own ion mode to measure other ion.
- d. Meter has the ability of automatic recognition standard buffer solution. You can select several pH standard buffer solutions to calibrate electrode. Total 10 different standard buffers are provided.
- e. Set own standard solution groups.
- f. With one-point, two-points and multiple-points calibration (mostly 5 points) function.
- g. With several concentration measuring modes: Direct Reading concentration measuring mode, Known Addition measuring mode, Sample Addition measuring mode and GRAN Method measuring mode.
- h. On the basis of high speed micro-processor, with high precision A/D, matched with precise measuring electrode, meter can ensure the high precision. The valid resolution of potential can be reached to 0.01mV.

#### **B.** Conductive Measuring Part

- a. Meter can measure Conductivity, Resistivity, Total Dissolved Solid (TDS) and Salinity.
- b. Within full measuring range, meter has automatic temperature compensation, automatic calibration, automatic measuring range and automatic frequency switching functions.
- c. With calibration function. User could calibrate Cell Constant or TDS Factor.

#### C. Dissolved Oxygen Measuring Part

- a. Meter can measure Dissolved Oxygen's Concentration, Saturation and Electrode current.
- b. With automatic temperature compensation function.
- c. With calibration function which can make Zero Oxygen calibration, Full Scale calibration, Air Pressure calibration and Salinity Compensation calibration.

#### D. Other Characteristic


- 1. Big 5.7" touch screen for easy use.
- 2. Apart from ion concentration, meter allows to measure all parameters simultaneously. Of cause, it also allows user to select necessary parameters for measuring.
- 3. Meter supports GLP:
- a. Meter requests user to set operator number and then can record down the procedure of all user.
- b. Can record down and print calibrating data.
- c. Can store 200 sets measuring data in each of pH measuring data, Conductivity measuring data, TDS measuring data, Salinity measuring data and Dissolved Oxygen measuring data if all are meet GLP. Meter can also store 6 kinds of ion measuring value including corresponding pX value and ion concentration value 100 sets of each.
- 4. Can view, print and delete measuring data
- 5. Can view the using parameter under current measuring mode as well as previous calibration data.
- 6. There are three measuring modes: Continuous measuring mode, Timed Reading measuring mode and Balance Reading measuring mode to meet different user's measuring requirements.
- 7. With standard RS232 interface.
- 8. With electricity cut off protection function which can protect the stored measuring data, calibrating data and set parameter to be saved when meter is turned off or abnormal electricity cut off occurred.

# 2. Main Specifications

- **1. Measuring Range:** a. mV:  $(-1999.99 \sim 0)$ mV,  $(0 \sim 1999.99)$ mV;
  - b. pH/pX:  $(-2.000 \sim 19.999)pH/pX$ ;
  - c. Concentration: various concentration values in correspondence to potential measuring range and indicating electrode.
  - d. Conductivity:  $0.000\mu\text{S/cm} \sim 199.9\text{mS/cm}$
  - e. Resistivity:  $5.00\Omega$ .cm  $\sim 20M\Omega$ .cm
  - f. TDS:  $0.000mg/L \sim 19.99g/L$
  - g. Salinity:  $(0.0 \sim 8.00)\%$
  - h. Dissolved Oxygen Concentration: (0.00 ~ 19.99)mg/L
  - i. Dissolved Oxygen Saturation:  $(0.0 \sim 199.9)\%$
  - j. Temperature:  $(-5.0 \sim 135.0)$  °C
- **2. Resolution :** a. pH/pX: 0.01/0.001pH
  - b. mV: 0.1/0.01mV
  - c. Dissolved Oxygen: 0.01mg/L
  - d. Saturation: 0.1% e. Temperature: 0.1°C
- **3. Electronic Unit Accuracy:** a. pH/pX:  $\pm 0.002$ pX  $\pm 1$  bit
  - b.  $mV: \pm 0.03\%$  (FS)
  - c. Ion concentration:  $\pm 0.5\% \pm 1$  bit
  - d. Conductivity:  $\pm 0.5\%$  (FS)  $\pm 1$  bit
  - e. TDS:  $\pm 0.5\%$  (FS)  $\pm 1$  bit
  - f. Salinity:  $\pm 0.1\% \pm 1$  bit
  - g. Dissolved oxygen: ±0.10mg/L
  - d. Dissolved oxygen saturation: ±1.0%
  - e. Temperature:  $\pm 0.3 \,^{\circ}\text{C} \pm 1 \, \text{bit}$
- **4. Input Impedance:**  $>3 \times 10^{12} \Omega$
- **5. Output Method:** touch 5.7" LCD screen with standard RS232 interface.
- 6. Normal working condition:
  - a. Ambient temperature:  $(5.0 \sim 35.0)$
  - b. Relative humidity:  $\leq 85\%$
  - c. Power supply: DC (9VDC, 800mA, inside positive, outside negative)
  - d. No performance-affected vibrations nearby.
  - e. No corrosive gases in surrounding
  - f. No strong electromagnetic interference except the geomagnetic fie
- 7. Outside Dimension: 1 x b x h: mm: 290 x 200 x 70
- 8. Weight: about 1kg


## 3. Construction

#### A. Front panel of meter


- 1. Main body
- 2. Touch display screen

#### B. Rare panel of meter

- 3. Switch
- 4. Socket
- 5. RS-232 interface
- 6. DO electrode socket
- 7. Conductive electrode socket
- 8. Reference electrode socket
- 9. pH electrode socket
- 10. Grounding socket
- 11. Temp electrode socket


#### C. Accessories


- 12. pH combined electrode
- 13. Electrode cover
- 14. Q9 short circuit plug


15. DJS-1C conductive electrode


16. DO-958-S dissolved oxygen electrode


17. T-818-B-6 temperature sensor

18. Connecting lines of printer

19. RS-232 communication line


20. 9V, 800mA power adaptor (inside is positive, outside is negative)

21. Multi-function supporter

# 4. Operation

#### 4.1. Installation

- a). Put Main body (1) and multi-function supporter (21) on table.
- b). If you need to measure solution's pH value, conductivity, dissolved oxygen etc..., please separately select relevant electrode attached with meter (or you can select other electrode). Please use pH combined electrode (12) to measure pH value, use DJS-1C conductive electrode (15) to

measure conductivity, TDS and salinity, use DO-958-S dissolved oxygen electrode (16) to measure dissolved oxygen, dissolved oxygen saturation, use T-818-B-6 temperature sensor (17) to measure temperature. Put relevant electrode on multi-function electrode supporter (21) and pull off the cover (13).

c) Pull off Q9 short circuit plug, insert pH combined electrode (12) into pH electrode socket (9), insert DJS-1C conductive electrode (15) into conductive electrode socket (7), insert DO-958-S dissolved oxygen electrode (16) into dissolved oxygen socket (6), insert T-818-B-6 temperature sensor (17) into temperature socket (11), insert 9V,800mA power adaptor (inside positive, outside negative) (20) into power socket (4)

If user has serial TP-16 printer (for detailed installation, setting and operating, please refer to the instruction of printer or relevant chapter of the meter), please insert printer's connection line (18) into RS-232 interface (5)

If user needs the communication connection, Please insert RS-232 communication line (19) into RS-232 interface (5).

d). When all connections are properly connected, turn on the meter and you can start to measure according to following procedure.

#### 4.2 Switch on/off the meter

Insert the power apparatus (9VDC 800mA power adapter, inside is positive and outside is negative) matched with meter into power socket then turn on the meter, meter displays 'MP-2006 multi-parameter analyzer'. Wait for a while, meter enters into start mode, user can use the meter now.

When you finish using, please insert short circuit plug into pH electrode socket, and then turn off the meter.

Note: Please follow below regulations and use the meter carefully to avoid any damage to the meter!


- a. Meter should be pre-heated for 0.5h before measuring
- b. For the measuring in ion part, in order to guarantee the precision of the measurement, suggest make zero point potential calibration after pre-heated for 0.5h. Usually when pH electrode socket is connected with short circuit plug, user could find the displayed potential value excurse from

zero point potential, at this time you need to calibrate zero point potential (zero point potential should be 0.00mV, you can select pH measuring parameter and view displayed zero point potential value in measuring mode). For detailed zero point potential calibration, please refer to the relative chapter about pH or pX measuring mode, usually when you enter into measuring mode, press 'SETUP' key, select 'Calibrate zero point potential' and then press 'Enter'.

- c. To protect and to use the meter better, before turn on the meter, please check the pH electrode socket on backside of meter and make sure it is connected with measuring electrode or short circuit plug, otherwise the high resistance apparatus of the meter will be damaged.
- d. If meter is not used, please connect the short circuit plug to avoid damage of the meter.

#### 4.3 Start Mode


The starting mode is shown as following chart, the upper left side of the screen shows present system time. On the bottom is present measuring mode including measuring parameter. On the right side is 'SETUP', 'VIEW' and 'MEAS' key. Press 'SETUP' on upper right side, you can set 'Measuring Mode', 'Ion Mode', 'System Time', 'Manual Temperature' and 'System Setup' etc... Press 'VIEW' in the middle, you can view 'pH STD Data', 'pX STD Data', 'Conductive Measuring Parameter', 'TDS Factor', 'Dissolved Oxygen measuring parameter', 'View Stored Data'. On lower right side is 'MEAS' key, you can start measuring by pressing this key.

Display chart of starting mode

#### **4.4 Set Function in Start Mode**

In this function, you can set 'Measuring Mode', 'Ion Mode', 'Time and Date', 'Manual Temperature' and 'System Setup' etc... Under starting mode, pressing 'Setup', meter displays setup menu which shown as follows:

User could press relevant menu item and then press 'ENTER' to select relevant function part. Press 'CANCEL' button or press 'X' (quit) to exit from Setup Menu.

- **Set Measuring Mode**: to set current measuring mode (Continuous measuring mode, Timed Reading measuring mode, Balance Reading measuring mode or Ion Concentration measuring mode) and measuring parameters.
- Set 'Ion Mode': to exchange ion mode, for example to change from measuring C1 to F.
- **Set 'System Time'**: to set time of the meter.
- **Set 'Manual Temperature'**: If meter is not connected to temperature sensor, meter would use manual temperature value.
- Set 'system Setup': to set necessary GLP selection.
- **Set 'Measuring Para.'**: from this part, you can set Cell Constant, TDS Factor, Temperature Coefficient, Air Pressure value of DO, Salinity compensation etc...


#### 4.4.1 Set Measuring Mode

Apart from ion concentration, meter can run ion measuring part, conductive measuring part and dissolved oxygen part simultaneously. Each part has one individual measuring display window to display three measuring parameters at same time. In which the upper one is the main measuring parameter, the other two are the auxiliary measuring parameter. You can only select main measuring parameter to start measure. The auxiliary measuring parameter is contained automatically by meter.

- The ion measuring part takes pH/ pX measuring parameter and Ion Concentration as main parameter while potential and temperature value as auxiliary parameter.
- Conductive part takes Conductivity, TDS and Salinity as main parameter while resistivity and temperature value as auxiliary parameter.
- Dissolved oxygen part takes Dissolved oxygen and Saturation as main parameter while dissolved oxygen current and temperature value as auxiliary.

Meter also allows user to select individual parameter to measure. Meanwhile, meter supports three measuring modes: Continuous measuring mode, Timed Reading measuring mode and Balance Reading measuring mode. Ion concentration supports four measuring modes: Direct Reading, Known Addition, Sample Addition and GRAN Method measuring mode. For easy use, meter has the function of setting measuring mode to support above operation.


Press 'SETUP' button and then press 'ENTER' button to set measuring mode which is shown as follows:

On left side, it is present selected measuring parameters, i.e. main measuring parameter chart. In the middle is measuring mode including Continuous Mode, Timed Reading measuring mode and Balance Reading measuring mode. When selected parameter is ion concentration, measuring mode would be changed to concentration measuring mode respectively which includes Direct Reading Mode, Known Addition mode, Sample Addition mode and GRAN Method mode, etc...

On the right side of the screen is press key, press 'MEAS' to start measurement. In case you want to start the same measurement next time, you can press 'ENTER' key, meter will save all present setting automatically and then return to Start Mode. Or you can press 'CANCEL' to exit and return to Start Mode.

On the bottom of screen, meter displays present measuring result which is shown as follows:


Measuring display chart when select conductivity, pH, DO measuring parameter.

The actual display will vary from the value of selected main measuring parameter. Each measuring part has one individual measuring window which displays the name of main measuring parameter, measuring value of main parameter and measuring value of auxiliary parameter. There is a '>' button on the left side of the window, you can select directly the other main measuring parameter under present part. For example, press '>' button beside conductivity window, meter will switch main measuring parameter to TDS measuring, press again, switch to salinity measuring, press once more, switch back to conductivity measuring (note: this method do not suitable for ion measuring part).

There is a "▼"or "▲" key on the right side of each measuring window, press this key to switch the displaying position of main measuring parameter (note: this method do not suitable for individual main measuring parameter).

#### 4.4.1.1 Select Measuring Parameter

Press 'Parameters' displayed on the left side of screen, meter shows 'Select Parameters' window which shown as follows:

Display chart of Select Parameters'


You can select ion measuring part, conductive measuring part and DO measuring part simultaneously or you can select individual parameter to measure. The left side of the screen is the selectable measuring parameters. The right side is selected measuring parameters. Press selected parameter display area once you can delete one selected parameter. Press relevant measuring parameter in 'Select parameter' area on left side, you can select measuring parameter. When selecting finishes, press 'X' (quit) to exit 'set menu'.

#### 4.4.1.2 Select Ion Concentration Measuring Mode

When you select measuring ion concentration, meter will shield other measuring parameter automatically, i.e. measure ion concentration should be done solely and can not be measure with other parameter at same time.

When you select measuring ion concentration, you need to select 'ION Conc' measuring mode including Direct Reading mode, Known Addition mode, Sample Addition mode and GRAN Method measuring mode etc... which is shown as following chart, user could select suitable concentration measuring mode according to measuring requirement.

For detailed ion concentration measuring, please refer to chapter 'ion concentration measuring'


Display chart of set ion concentration measuring mode

#### 4.4.1.3 Select Measuring Mode

Meter supports three measuring modes: Continuous measuring mode, Timed Reading measuring mode, Balance Reading measuring mode which can be selected directly at current measuring mode display area.

#### 4.4.1.3.1 Continuous Measuring Mode

It is the most normal measuring mode, when you start measurement, meter will measure, calculate and display measuring result continuously. Please refer to following display chart. During measurement, user could view measuring parameter, calibrate electrode, store or print measuring result. When measurement finishes, you can press 'X' (quit) exit the measuring mode and return to start mode.


Display chart of continuous measuring mode

#### 4.4.1.3.2 Timed Reading Measuring Mode

It is particularly set for easy use. For example, if you need to measure conductive data continuously in 30 minutes, you can select 'Timed Reading' measuring mode. When measurement starts meter will automatically measure, calculate and display measuring result till the set time up and then to save the measuring data automatically (if it is connected with serial printer, meter will print automatically) and start next measurement.

Press 'X'(quit) and 'Enter' key to exit from 'Timed Reading' measuring mode.

If user selects 'Timed Reading' measuring mode and needs to set interval, the default interval is 10 minutes which displayed as follows:


Display chart of 'Timed Reading' measuring mode

#### **4.4.1.3.3** Balance Reading Measuring Mode

It is another measuring mode. Firstly set balance condition (**please refer to system setup chapter for detailed operation**). When measurement starts, meter will automatically measure, calculate and display measuring result. As soon as the measuring meets the balance condition set, measurement finishes.

During measurement, user could view measuring parameter, calibrate electrode etc... When measurement finishes, user could store, print measuring result and press 'X' (quit) on upper right corner to exit the measuring mode or to start next measuring by pressing 'CONT'.


Display chart of 'Balance Reading' measuring mode

#### 4.4.2 Set Ion Mode

It is designed mainly for easy use. Meter provides about 10 general kinds of ion mode in correspondence with different ion measurements and user could select relevant ion mode to make concentration measurement. When concentration measurement finishes, user could view current ion concentration value according to different ion concentration unit.

Press 'SETUP' button, select 'Ion mode' and then press 'ENTER' to enter into 'set ion mode' part which is shown as follows:

Display chart of 'Set Ion Mode'

Meter divides ion into four valences: one-valence ion mode, two-valence ion mode, three-valence ion mode and custom ion mode. Meter provides some general ion modes. One-valence ion (including one-valence cation, one-valence anion) has 10 kinds. They are: H<sup>+</sup>, Ag<sup>+</sup>, Na<sup>+</sup>, K<sup>+</sup>, NH<sub>4</sub><sup>+</sup>, Cl<sup>-</sup>, F<sup>-</sup>, NO<sub>3</sub><sup>-</sup>, BF<sub>4</sub><sup>-</sup>, CN<sup>-</sup>, etc...(because H<sup>+</sup> is always allowed, it is not displayed in 'Set Ion Mode'). Two-valence ion is three kinds of cation: Cu<sup>2+</sup>, Pb<sup>2+</sup>, Ca<sup>2+</sup>. On the right side of screen, there is a

four direction key, you can select relevant ion by pressing these four keys (display highlight). In the middle of the screen is the ion name and molecular weight currently selected, press 'ENTER', meter will select this ion mode as actual measuring ion mode, press 'X' key to exit from 'Set Ion Mode' function and return to start mode.

User could set his own ion mode when he has relevant ion select electrode. User could also measure ion concentration according to normal operating procedure of ion mode.


Three keys on the bottom of screen are 'NEW', 'EDIT', 'DELETE' key. They are designed to match custom ion mode. Press 'NEW' button, user could define ion mode. Meter will request user to input ion name and relevant molecular weight. When setting finishes, meter will add this ion mode in ion mode storage automatically. Next time you can just select this ion mode to measure. User can also edit custom ion mode or delete custom ion mode.

You could only start concentration measuring when you select correct ion mode. If you select different ion mode, it will cause incorrect measuring result. For example, if you want to measure concentration of sodium ion, then you should select 'Na' ion mode in 'Set Ion Mode' and start to measure concentration of sodium ion.

Meter supports maximum 20 custom ion modes.


#### 4.4.3 Set System Time

Press 'SETUP' key, select 'Set Time and Date' and press 'ENTER'. Meter enters into set system time part. User can set present system time. On the upper side of the window is present system time while at the bottom is the corresponding time key including 'year', 'month', 'day', 'hour', 'minute', 'second' key. It is displayed as follows:

Display chart of 'Set System Time'

Press corresponding time key to set relevant time. For example, if user wants to set present month, you can operate as follows: press 'month' key once, meter will add one on present figure. At the same time meter displays another two press keys, i.e. 'Enter' and '...' key. It is displayed as follows:

Display chart of 'Set System Time'


Press 'month' key repeatedly till your needed month and press 'Enter' to finish setting. (Or you can input the figure: press 'Input' key, there is an input window jump out for you to input necessary month. When input finishes, press 'Enter' to exit from input window. Press 'Enter' again, you can fulfill the setting of present month).

Your can set other time according to the same procedure.

When setting finishes, press 'X' (quit) key to exit 'Set System Time' mode and return to start mode.


#### **4.4.4 Set Manual Temperature**

When temperature electrode socket is connected with temperature sensor, meter will automatically adopt the temperature of temperature sensor. Otherwise, meter will adopt manual temperature set by user as current temperature. Press 'SETUP' key, select 'Set Manual Temperature' and press 'Enter', meter enters into 'Set Manual Temperature' part. It is displayed as follows:

Display chart of 'Set Manual Temperature'

Press once at 'Manual Temperature' displaying area, meter shows an input window which is displayed as following chart. You can input relevant manual temperature value in input window. When input finishes, press 'Enter' to finish the input. In case there is any mistake occurred, press 'Delete' button to input again. Press 'Cancel' to quit the input. When setting finishes, meter returns back to start mode.

Display chart of 'Set Manual Temperature'


Note: All parameters needed to amend will use the similar way to achieve.

All chapters about setting manual temperature here after will be operated according to this.

#### 4.4.5 System Setup

This function is GLP regular setting includes setting electrode calibration interval, Balance Reading measuring condition, operator's number, whether to delete stored data automatically and whether open GLP function.

Press 'SETUP' key, select 'System Setup' item, press 'Enter' to enter into system setting mode. It is displayed as follows:

Display chart of 'System Setup'

Press relevant item to amend corresponding parameter. When amendment finishes, press 'X' (quit) to exit, then return to start mode.

Electrode calibration interval means meter will remind user the time interval of calibrating electrode. Meter will automatically calculate the time from previous calibration till now. In case previous calibration time exceeds the time interval set by user, meter will jump a reminding window to remind user take care and re-calibrate electrode. The unit of electrode calibration interval is hour. The displayed pH electrode calibration interval is 4 hours.

The operator's number is a three digit number, the range is between  $000 \sim 200$ , all operation records contain operator's number.

Automatic delete stored data function means when stored data reaches the amount that the meter

allows, whether you agree to cover or repeat the storage. For example, the permitted stored pH measuring data is 200 sets, if you want to store number 201 data, in case automatic delete stored data function is opened, meter will store the number 201 data in first position, that means meter will store data from the beginning, otherwise meter will give up storing current data. Please pay special attention on it.


GLP function means meter has two modes, i.e. either select supporting GLP or do not select supporting GLP. There will be different pattern in search and print stored data under these two modes.

'Balance Reading' measuring condition is in correspondence to the Balance Reading condition of each measuring parameter, displayed pH balance condition is 0.01pH. When change amount of pH measuring value is less than 0.010pH, the measuring is valid. This function is mainly for 'Balance Reading' measuring mode. When user selects 'Balance Reading' measuring mode to measure and if all measurement within set balance time are accord with balance condition, the current measurement finishes.

Balance time is only valid in 'Balance Reading' measuring mode and the unit is second (s), range is  $1 \sim 200$  seconds.

#### 4.4.6 Set Measuring Parameter

Select 'Set Measuring Para.', meter jumps out a selectable menu which displayed as follows:


Displayed chart of 'Set Measuring Parameter'

User can enter into relevant parameter set function part according to practical need and press 'CANCEL' to quit.

#### 4.5 View Function

For easy use, meter allows user to view present measuring parameter at random including previous calibrating data and present using parameter, allows user to edit measuring parameter directly. Meter also allows view stored data.


In start mode, press 'View' key, select relevant item to view, edit measuring parameter and view stored data which displayed as follows:

Display chart of 'View' in start mode

#### 4.5.1 View pH STD Data

Under this function, user could view current pH calibration data. In start mode, press 'View' key, then press 'Enter' to view 'pH STD Data'. Displayed as follows:

Display chart of 'View pH STD Data'

On top of the screen is last calibration data including calibration time, calibration point, calibration data and current slope in percentage in correspondence to each STD. At the bottom is the condition of STD buffer group: on right side is "<=>"(move) key, 'PRINT' key, 'STD' key and 'SETUP' key.

If last calibration is multi-point calibration, user can press "<=>"(move) key to show detailed data of each standard solution point.

Here, user could calibrate electrode slope (please refer to calibration chapter for detailed calibration procedure) directly by pressing 'STD'.

Press "SETUP' key, user could set new STD buffer group.

Press 'X'(quit) key to exit 'View pH STD Data'.

If user wants to print present calibration data, you could connect the connecting line (attached with meter) with printer TP-16 (should be bought separately) and press 'PRINT' to print out calibration data.


We suggest using TP-16 printer. If you want to use other serial printer, the printing content will be stagger. Please note.

- 1. The communication agreement of meter's RS232 is: 9600, N, 8, 1, i.e. baud ratio 9600bps, no odd/even, 8 bit, one stop bit.
- 2. When to connect serial printer, please set baud ration according to instructions and set shake hand method (busy signal) to standard control method.
- 3. Only connect the printer when the power is cut off and printer is turned off.

#### 4.5.1.1 Setup pH STD Buffer Group

Meter has automatic distinguishing function, can distinguish 10 standard buffer solutions. You can make multi-points calibration but mostly not exceeds 5 points. Because pH range of 10 standard buffer solutions maybe duplicated, in order to ensure the precision of measuring result, user needs to establish relevant STD buffer group. For example, if you want to adopt 4.003pH, 9.182pH standard buffer solution to calibrate electrode slope, you can set STD buffer group as 4.003pH and 9.182pH. In normal practice, meter can automatically distinguish these two standard buffer solutions.

Press 'SETUP' key to enter into current STD buffer group, displayed as following chart. In which there are 10 figures in correspondence to 10 standard buffer solutions in 'Setup pH STD Buffer Solution' window. These 10 standard buffer solutions are 1.679pH, 3.557pH, 3.776pH, 4.003pH, 6.864pH, 7.000pH, 7.413pH, 9.182pH, 10.012pH and 12.454pH. Under the figure is standard pH value of STD, for example the figure above 3.56 is in correspondence to 3.557pH STD. The figure with tick means this standard buffer solution has been selected while no tick on figure means the solution has not been selected. Opposite color means present STD is operational. Press the figure with opposite color once more, you can select or delete current STD. The relevant STD range will be shown under the window.

'Setup pH STD Buffer Group' in 'pH STD Data' mode

For example if user wants to select 3.557pH standard buffer solution, press the figure in correspondence to 3.557pH STD buffer, the figure which indicates 3.56pH will display the opposite color, press this figure with opposite color again, the figure indicates 3.56pH will be ticked which means it has been selected.


To avoid duplication of pH value between standard buffer solutions, you should select the STD buffer which you actually needed and delete all the other STD buffer which you actually don't needed

When selection finishes, press 'X' (quit) key to exit 'Set pH STD Buffer Group' and return to 'View pH STD Data' mode.

(Please refer to appendix for distinguishable 10 standard buffer solutions.)

#### 4.5.2 View 'pX STD Data'

In start mode, press 'VIEW' key, select 'pX STD Data', then press 'ENTER' to view 'pX STD Data'. Displayed as follows:

Display chart of view 'pX STD Data'

On top of the screen is 'Previous STD' including STD time, STD point and present slope in percentage in correspondence to STD buffer. At the bottom is the 'Current Ion Mode', i.e. all above pX STD Data are in correspondence to this ion mode: on right side is '<=>'(move) key, 'PRING' key, 'STD' key and 'SETUP' key.

If previous STD is multi-points calibration, user can press '<=>'(move) key to show detailed data of each STD point.

From here, user could calibrate electrode slope (please refer to calibration chapter for detailed calibration procedure) directly by pressing 'STD'.


Please press 'PRINT' to print out STD data when the meter is connected with model TP-16 printer.

Press 'SETUP' button, you could set ion mode (please refer to chapter 'set ion mode' in previous page). Press 'X'(quit) to exit view 'pX STD Data' mode.

#### **4.5.3 View Conductive Measuring Parameter**

Under start mode, press 'VIEW', select 'Cond Meas Para', then press 'ENTER' to view Conductive measuring parameter which displayed as follows.

On top of the screen is Previous STD, at the bottom is present conductive parameter value including 'Cell Constant' and 'Temperature Coefficient'. On right side is 'PRINT', 'STD' and 'SETUP' key for setting cell constant and temperature coefficient.

Display of View 'conductive measuring parameter'

Please press 'PRINT' to print out calibration data when the meter is connected with model TP-16 printer.

# 4.5.3.1 Setup 'Cell Constant'


Generally, each conductive electrode from our company has relevant conductive constant, you only need to set constant value on electrode once before using. After using a period of time, if you suspect the electrode constant is not precise, you can calibrate with conductive standard solution again. Press 'STD' key, you can calibrate electrode slope (please refer to chapter Calibration for detailed procedure).

So there are two ways to set present cell constant value: one is to re-calibrate with conductive standard solution, as soon as calibration finishes, meter will calculate cell constant automatically. Another ways is to set cell constant value directly.

You can only select one method. Presume you get cell constant value through calibration the previous time and now you want to input cell constant directly, meter will delete previous calibration data, so please take care.

Press 'SETUP' key beside 'Cell Constant', there will be an input window, please input new cell constant then press 'ENTER'. Displayed as follows:

Press 'X'(quit) to exit 'Set Cell Constant' mode.

Display chart of 'Set Cell Constant'

#### 4.5.3.2 Setup Temperature Compensation Coefficient

If meter needs precise measuring, the influence of temperature will cause inaccurate measuring of conductivity. At this time we need to set temperature compensation coefficient, the acquiescent value is 0.020.


In 'Cond Meas Para' mode, press 'SETUP' key beside temperature coefficient, meter will show an input window, you can input new temperature coefficient and press 'ENTER'.

Generally it is not necessary to set temperature compensation coefficient, the acquiescent temperature compensation coefficient is 0.020.

# 4.5.4 View Dissolved Oxygen Measuring Parameter

Under start status, press 'VIEW' key, select 'DO Meas Para' and then press 'ENTER' key to view measuring parameter of dissolved oxygen. The display is as follows.

On top of the screen is 'Previous STD', at the bottom is present 'Air Pressure' and DO 'Salinity Factor'. On right side is 'PRINT', 'STD' and 'SETUP' key for setting air pressure and DO salinity factor.

Display chart of 'Dissolved Oxygen Measuring Parameter'

#### 4.5.4.1 Setup 'DO Air Pressure'

Press 'SETUP' key beside 'Air Pressure' to set present air pressure, an input window jumps out and you can input new air pressure.

Generally it is not necessary to set air pressure. The acquiescent air pressure is 101.3kPa.

# 4.5.4.2 Setup 'DO Salinity'

Press 'SETUP' key beside 'Salinity Factor' to set salinity, an input window jumps out and you can input new salinity.

Generally it is not necessary to set salinity. The acquiescent salinity is 0.0g/L.

#### 4.5.6 View Stored Data

When measuring finishes, you can store measuring result for easy viewing and printing in future. Meter stores data according to measuring parameter and supports GLP. You can store 200 sets in each of pH value, conductivity, TDS, salinity and saturation etc...which accord with GLP. Apart

from hydrogen ion, all other ion can make pX measuring and concentration measuring. So the storage is according to pH measuring result and concentration measuring result. Also all the measuring result are stored according to present ion mode. For example, user select sodium ion mode to measure concentration of sodium ion, meter will automatically assign a storage space for user to store the measuring result. 100 sets pNa measuring results (pX measuring result) and 100 sets concentration measuring results can be stored. When user selects another ion mode to measure other ion concentration, meter will similarly assign the same area to this ion mode to store measuring result.

Note: Because of the limitation of storage, meter only allows to store 6 measuring results of ion mode.

Under start mode, press 'VIEW' key, select 'View Stored Data' and then press 'ENTER' key to view stored data which is displayed as follows.

On top of the screen is mode name and stored No. in correspondence to present view mode. Each page displays mostly 4 stored data and displayed condition differs from view mode and also related to ON/OFF of GLP mode.

• pH mode includes storing time, operator's number, slope, potential, pH value and temperature etc...


On right side of the screen is press button including 'Print', 'Delete', 'Front Page' and 'Next Page'. On left bottom is search mode selection button and right bottom is the amount of actual stored ion mode.

You can press directly the displayed data and select relevant stored data, press 'PgUp', 'PgDn' to view the stored data in other pages.

If user wants to print, delete stored data, press 'PRINT', 'DEL' key to do so.

Press 'Current Mode' key, you can view in sequence pH stored data, pX stored data, ion concentration stored data, conductivity stored data, TDS stored data, salinity stored data, DO stored data, saturation stored data, etc...

Press 'Ion Mode No' key, you can select other stored ion mode. The meter can store up to 6 ion modes.

Display chart of 'View Stored Data'

#### 4.6 Ion Measurement


The principle of ion measurement is: take various ion selection electrodes as indicate electrode, supplement with suitable reference electrode and insert all electrodes into measured solution to structure electrochemistry system for measuring.

#### 4.6.1 pH Measurement

It is the most normal ion measuring mode. We particularly pick up pH measuring mode for easy use and add automatic distinguishing pH STD buffer function in calibrating electrode.

In start mode, if there is 'pH measuring parameter' mode, then press 'MEAS' directly to start measurement. Otherwise press 'SETUP' button, select 'Setup Measuring Mode' and then select 'pH STD Data (please refer to chapter 'Set Measuring Mode' for details).

When meter enters into pH measuring mode, it is displayed as following chart. On top of the screen is current measuring mode, system time. In the middle is pH STD data, relevant potential and current temperature. On right side of the screen are press keys including 'View', 'STD', 'SORE', 'PRINT' and 'SETUP' keys. Timed reading measuring mode and Balance Reading measuring mode will shield 'STORE', 'PRINT' key automatically.

Display chart for pH measuring status

Press 'View' key to search pH STD data. Please refer to chapter 'View pH STD Data' for details. Press 'STD' key to re-calibrate pH electrode. Please refer to chapter 'pH Electrode Calibration' for details.

During measurement, press 'SETUP' key to set manual temperature, mV displayed resolution, pH displayed resolution and calibrated potential zero point.

If user selects **continuous measuring mode** (**please refer to chapter 'Setup Measuring Mode' for details**), when measuring result stable, user could store current measuring result by pressing 'STORE' key. If it is connected with TP-16 printer before turning on, meter will print current measuring data by pressing 'PRINT' key.

If user selects **Timed Reading measuring mode (please refer to chapter 'Setup Measuring Mode' for details),** meter will automatically measure, calculate and display measuring result. When set time interval reaches, meter will store measuring result automatically and start next measuring after printing finishes.

If user selects Balance Reading measuring mode (please refer to chapter 'Setup Measuring


**Mode' for details),** when measuring starts, meter will automatically measure, calculate and display measuring result. In case measuring in accordance with set balance condition, the measuring finishes! At this time, you can press 'STORE' key to store present measuring result. If it is connected with TP-16 printer before turning on, meter will print present measuring data by pressing 'PRINT' key.

Press 'X' key on top right corner and then press 'ENTER' to finish the measurement.

# 4.6.1.1 pH Electrode Calibration

- 1. Before each measurement, suggest re-calibrate the electrode. When calibration starts, previous calibration data will be recovered.
- 2. Electrode should be re-calibrated after a period of using
- 3. Meter will remind automatically when the electrode calibration time interval is set to suitable value.

Under pH measuring mode, press 'STD' key, select 'pH Electrode Calibration' and then press 'ENTER' key or in 'View pH STD Data' mode, press 'STD' key, you can calibrate electrode slope also. Meter reminds 'Calibrate Electrode?' press 'ENTER' key, meter enters into calibration mode which is displayed as follows:

Display chart for calibration function under pH ion mode

On upper half of the screen is current measuring data, meter displays current pH value (slope is set as 100.00%), potential value and temperature value. In the middle of the screen is 'Current STD Buffer Group' set by user. At the bottom is current 'STD Result'. On right side is 'AUTO' (or 'MANU') key, "SETUP' key and END' key.


Automatic distinguish (or manual distinguish) means whether present distinguish mode is automatic distinguish or manual distinguish, press 'AUTO' (or 'MANU') key, you can switch between automatic distinguish and manual distinguish quickly. Press 'SETUP' key to set current STD buffer group. In standard solution, when electrode displaying is stable, press 'ENTER' key to calibrate present solution, press 'CANCEL' key to quit.

Meter has automatic distinguishing function, can distinguish 10 standard buffer solutions. Because pH range of 10 standard buffer solutions maybe duplicated, in order to ensure the precision of measuring result, before measuring, you'd better check the set STD buffer group. For example, you want to use 4.003pH, 9.182pH STD buffer to calibrate electrode slope, then you can set STD buffer

group as 4.003pH and 9.182pH, otherwise meter could not distinguish these two standard buffer solution automatically which will cause inaccurate of calibration result or even cause the mistake of calibration. In calibration mode, press 'SETUP' key beside 'Current STD Buffer Group' to set current STD buffer group (please refer to chapter 'Setup pH STD Buffer Group' for details). Calibration procedure is as follows:

- a. Before calibration, please prepare 1 to 5 standard buffer solution (could be normal STD buffer or STD buffer prepared by user himself), place the solution in constant temperature for certain period.
- b. Set STD buffer group according to the procedure stated before. If it is abnormal standard buffer solution, please select manual distinguish method.
- c. Clean pH measuring electrode, reference electrode and temperature sensor, and then put them in measured standard buffer solution.
- d. When display stable, press 'ENTER' key, meter displays 'Store STD Data...' and store the calibration data.
- e. Later, meter reminds user 'End Standardizing?' which is displayed as follows. If user has some other standard buffer solution to be calibrated, please select 'CONTINUE' key and repeat previous procedure to calibrate other solution till calibration finishes. During calibration, user could press 'END' key anytime to finish the calibration.

Meter supports mostly 5 point calibration. When the fifth solution is calibrated and confirmed, meter will finish calibration automatically.

Display chart when a certain standard solution calibration finishes.

For the standard buffer solution which pH ranges are duplicated, such as 6.864pH and 7.000pH STD, we suggest calibrate as follows:

When user calibrate 6.864pH standard buffer solution, set STD buffer group to 6.864pH only and then calibrate. When calibration of 6.864pH STD, STD buffer group again to 7.000pH only and calibrate.

Use manual distinguish method to calibrate, i.e. each time when you calibrate standard buffer solution, manually input pH value of current STD which is in correspondence to current

temperature. You can calibrate in this way but it is a little complicated.

For normal standard buffer solution, you can use automatic distinguish function. To match with previous STD buffer group set before, meter can automatically distinguish these standard buffer solutions. You don't need to change distinguish method to calibrate (if it can not distinguish, meter will remind user: wrong calibration and request user to exchange electrode or to re-set STD group, or to exchange automatic distinguish method to manual distinguish method, you can select either according to the normal practice.)

If you use your own standard buffer solution (not normal STD) to calibrate electrode, you must select manual distinguish method.

For example, user has one standard buffer solution and the known pH value at 25°C is 2.704pH, at 25.1°C is2.710pH, 25.2°C is 2.720pH, then we suggest stabilize the calibration temperature at 25°C. When calibration starts, firstly set distinguish method to 'MANU', wait till display stable, press 'ENTER' key, meter will request user to input pH value at present temperature, if present temperature is 25.2°C, please input 2.720pH and press 'ENTER' key. When input finishes, meter will store present calibration data.

If user either has normal standard buffer solution or has his own standard buffer solution then he can operate separately according to automatic distinguish method and manual distinguish method.

# 4.6.1.2 Setup mV Resolution

You can select mV displayed resolution, i.e. 0.1mV or 0.01mV. Press 'SETUP' key, select 'Setup mV Resolution' and then press 'ENTER' key to set mV displayed resolution. When setting finishes, meter returns to measuring mode.

All parts hereafter about 'setup mV Resolution' are operated according to this instruction.

#### 4.6.1.3 Setup pX Resolution

You can select pX displayed resolution, i.e. 0.01pX or 0.001pX. Press 'SETUP' key, select 'Set pX resolution' and then press 'ENTER' key to set pX displayed resolution. When setting finishes, meter returns to measuring mode.

All parts hereafter about 'pH, pX resolution' are operated according to this instruction.


#### 4.6.1.4 Calibrate Zero Point Potential

To ensure the high precision of meter, please calibrate potential zero point 0.5h after turning on the meter. Connect the meter with short circuit plug, when potential displayed value excurse from zero point potential (0.00mV) far, you need to calibrate zero point potential. Press 'SETUP' key, select 'Calibrate zero point potential' and press 'ENTER' key, meter reminds you 'Calibrate current zero point potential?', press 'ENTER' key to calibrate zero point potential. Otherwise, press 'CANCEL' to give up the operation and return to measuring mode.

All parts hereafter about 'Calibrate zero point potential' are operated according to this instruction.

#### 4.6.2 pX Measuring

You can make pX measuring in correspondence to each ion mode. When you select relevant ion mode, you can start pX measuring. For example, if you have sodium ion electrode and want to measure sodium ion Na in solution, you need to select present ion mode as Na ion mode in 'Setup Ion Mode' then to begin normal measuring. Please see 'Setup Ion Mode' for detailed operation Under start mode, if you have pX measuring parameter, then press 'MEAS' key directly to start measuring. Otherwise, press 'SETUP' key, select measuring mode and the select 'pX STD Data'. (Please see 'Setup Measuring Mode' for details)

pX measuring mode

The displayed interface is the same as pH measuring mode. On top of the screen is current measuring mode, system time, pX measuring result and corresponding potential and temperature value. On right side of the screen is press buttons including 'VIEW', 'STD', 'STORE', 'PRINT' and 'SETUP' key.

Press 'VIEW' key to view 'pX STD Data'. For details, please refer to 'View pX STD Data' Press 'STD' KEY, you can re-calibrate pX electrode, please refer to 'pX electrode calibration' for details.

During measuring, you can press 'SETUP' key to set manual temperature, mV displayed resolution, pX displayed resolution and calibrate zero point potential.

If user selects continuous measuring mode (please refer to chapter 'Setup Measuring Mode' for details), when measuring result stable, user could store current measuring result by pressing 'STORE' key. If it is connected with TP-16 printer before turning on, meter will print current measuring data by pressing 'PRING' key.

If user selects 'Timed Reading' measuring mode (please refer to chapter 'Setup Measuring Mode' for details), meter will automatically measure, calculate and display measuring result. When set time interval reaches, meter will store measuring results automatically and start next measuring after printing finishes.

If user selects 'Balance Reading' measuring mode (please refer to chapter 'Setup Measuring Mode'


for details), when measuring starts, meter will automatically measure, calculate and display measuring results. In case measuring is in accordance with set balance condition, the measuring finishes! At this time, you can press 'STORE' key to store current measuring result. If it is connected with TP-16 printer before turning on, meter will print current measuring data by pressing 'PRINT' key.

Press 'X' on top right corner and then press 'ENTER' to finish the measurement.

# 4.6.2.1 pX Electrode Calibration

Before each measurement, suggest re-calibrate the electrode. When calibration starts, previous calibration data will be recovered. Electrode should be re-calibrated after a period of using.

Under pX measuring mode, press 'SETUP' key, select electrode calibration and then press 'ENTER' key. Or in 'View Measuring Parameter' unit, you can press 'STD' to calibrate electrode slope also. Meter reminds 'Calibrate Electrode Slope?' press 'ENTER', meter enters into calibration mode which is displayed as follows. On upper half of the screen is present measuring data, meter displays present pX value (slope is set as 100.00%), potential value, temperature value and STD value (pX value) in correspondence to present solution. At the bottom of the screen is present STD result. On right side is 'ENTER', 'SETUP', and END' key.


Display chart for a certain ion mode is calibrated

Before calibration, please prepare 1 to 5 standard buffer solution (could be normal standard buffer solution or standard buffer solution prepared by user himself), place the solution in constant temperature for certain period.

- a. Clean relevant ion select electrode, reference electrode and temperature sensor, and then put them in measured standard buffer solution.
- b. Press 'SETUP' key and input present STD value (pX value) in correspondence to present standard solution.
- c. When display stable, press 'ENTER', meter displays 'Store STD data...' and then store the calibration data.
- e. Later, meter reminds user 'End standardizing?' which is displayed as follows. If user has some other standard buffer solution to be calibrated, please select 'CONTINUE' key and repeat

previous procedure to calibrate other solution till calibration finishes. During calibration, user could press 'END' key anytime to finish the calibration.

During calibration, each time when you change standard solution, please made sure to input relevant STD value (pX value) in corresponding to current standard solution.

Display when a certain ion mode calibrating a certain standard solution finishes

Note: When calibration of the fifth standard solution finishes, meter will stop calibration automatically.

# 4.6.2.2 Setup mV Displayed Resolution

The same as 'Setup mV Display Resolution' in pH measuring

#### 4.6.2.3 Setup pX Display Resolution

The same as 'Setup X Display Resolution' in pH measuring

#### 4.6.2.4 Calibrate Zero Point Potential

The same as 'calibrate zero point potential' in 'pH measuring'

#### **4.6.3** Ion Concentration Measurement

It is used to measure ion concentration of solution. Meter supports four concentration measuring modes, including direct reading concentration mode, known addition mode, sample addition mode and GRAN measuring mode.

# 4.6.3.1 Directly Reading Concentration Mode

The mode has following formula:

$$E_x = E_0 + S \times \log(C_x + C_b)$$

In which:  $E_x$  balance potential of measured sample.

 $E_0$ ~ zero potential value


 $S \sim$  electrode slope

 $C_x$  concentration value of measured sample

 $C_b \sim$  blank concentration value.

So, user only needs to calibrate relevant slope, gets slope and zero potential value to calibrate concentration of measured sample. If you need to calibrate the concentration value (blank concentration value) of blank standard solution, you can select 'Cblank'.

User selects relevant ion mode according to ion electrode used (for example, if you want to measure Ag ion concentration, you can select ion mode as Ag ion, please refer to chapter 'Setup Ion Mode' for details). Set measuring parameter as ion concentration and select direct reading concentration mode (please refer to chapter 'Setup Measuring Mode), when selection finishes, press 'MEAS' key to start the measuring or press 'ENTER' to return to start mode. Under start mode, press 'MEAS' key, you can also make direct reading concentration measuring which is displayed as follows:


Display chart for 'Direct Reading Concentration Mode' measuring

On top left side of the screen are current system time, present measuring potential, temperature value and relevant pX value. In the middle is current measuring result and relevant blank concentration value. At the bottom of the screen is previous STD result. On right side are pressing keys including 'ENTER', 'END', 'BLANK' and 'SETUP' key.

Press 'BLANK' key to calibrate blank concentration value of ion.

Press 'SETUP' key, you can set manual temperature, view calibration data, calibrate electrode, select concentration unit, set mV display resolution and set pX display resolution, etc...

Please clean the relevant ion selection electrode and put in measured solution, meter displays current measuring value. When reading figure is stable, press 'ENTER' key, meter will calculate measuring result which is displayed as follows:


Display chart when 'Direct Reading Concentration Mode' measuring finishes At this time, press 'STORE' key to store current measuring result, press 'PRINT' key to print out the result, press 'CONT' to continue next ion concentration measurement, press 'EEN' to quit direct reading concentration measuring mode and return to start mode.

If you want to select concentration unit of sample's concentration, you can operate according to following procedure: press 'Sample Cs' once, the sample's concentration unit will switch one time automatically.

#### 4.6.3.1.1 Blank Concentration Calibration in Direct Reading Mode Measuring


If you want to make blank concentration calibration, press 'BLANK' key and meter displayed as follows:

Display of Blank Concentration Selection

Press 'ENTER' key, Meter displays 'Cali Blank Conc?' press 'ENTER' to start blank concentration calibration. Select 'Clear Blank Conc', you can make present blank concentration zero.

When meter enters into blank concentration calibration, the display is the same as direct reading concentration measuring mode which is displayed as follows:

Direct read concentration calibrating blank concentration

When user has prepared blank standard solution, put relevant ion select electrode and temperature electrode into solution, wait till display stable, press 'ENTER' key, meter will calculate blank concentration value and store it automatically.

Press 'ENTER' key, meter returns to direct reading concentration measuring mode.

During Calibration, you could press 'END' key any time to quit the calibration and return to direct reading concentration measuring mode.

- 1. During blank calibration, all blank solution should be similar with the blank solution in chemical analyze.
- 2. There is blank calibration in direct reading concentration mode and known addition mode. Detailed operation is the same as this measuring mode.

#### 4.6.3.1.2 Clear Blank Concentration

If you want to clear previous blank concentration value, you can press 'BLANK' key under direct reading concentration measuring mode and select 'Cali Blank Conc? Clear Blank Conc' then press 'ENTER' to make blank concentration value to zero.

#### 4.6.3.2 Known Addition Mode

Known addition is also called standard addition. Firstly, measure balance potential value of system and then add standard solution with known concentration in measured system. Measure balance potential value of system again, and then calculates the concentration value of measured sample according to the changing value before and after adding. The formula is as follows:

$$Cx = \frac{\rho \times Cs}{(1+\rho) \times 10^{(E2-E1)/S} - 1} + \frac{\rho \times Cb}{(1+\rho) \times 10^{(Eb2-Eb1)/S} - 1}$$

In which:  $Cx^{\sim}$  concentration value of measured sample

Cs concentration value of standard solution (added solution)

 $S \sim$  electrode slope

Cb~blank standard concentration value

E 1∼Measured potential value before adding standard solution.

E 2~Measured potential value after adding standard solution.

 $\rho \sim \text{standard solution adding volume(Vs)/measured sample volume (Vx)}$ 


E b1 ~ Measured potential value before adding standard solution in blank calibration

Eb2~Measured potential value after adding standard solution in blank calibration

Before measuring, input standard solution's concentration value and adding volume, then input sample's volume. Later measure electrode potential value E1 before adding and then measure electrode potential value E2 after adding, meter can calculate sample's concentration value Cx according to above formula. If you want to make blank calibration, you can separately measure blank standard solution's electrode potential changing value before and after adding standard solution according to the similar method, i.e. Eb1, Eb2 and then calculate the blank concentration value of blank standard solution.

Similarly, you can select relevant ion mode according to ion electrode used (for example, if you want to measure Na ion concentration, you can select ion mode as Na ion, please refer to chapter 'Setup Ion Mode' for details). Set measuring parameter as ion concentration and select 'Known Addition' mode (please refer to chapter 'Setup Measuring Mode), When selection finishes, press 'MEAS' key to start the measuring. Or you can press 'ENTER' to return to start


mode. Under start mode, press 'MEAS' key, you can also make known addition measurement. Firstly you need to set some parameter which is displayed as follows:

Parameter setting in 'Known Addition Measuring Mode'

On upper side of the screen is the adding volume. The volume before adding is the sample volume. 'STD Conc' means the concentration value of added standard solution. The relevant concentration unit is listed on right side. The displayed concentration of standard solution is 100mol/L. In the middle of the screen is blank concentration value (Cblank). At the bottom is previous STD. On right side is 'ENTER', 'CANCEL', 'BLANK' and 'STD' key.

According to practical condition, set the parameter of present known addition measuring mode, if next adding volume is 1.5mL, then you can operate as follows: press once in displayed 'adding volume' area, meter will open the input window automatically which is displayed as follows:

Display chart for set adding volume


You can set other parameter in sequence according to above method, press once on 'ConcUnit' to select the concentration unit of standard solution. When you select it, the color of 'ConcUnit' is opposite (displayed concentration unit is 100.0mmol/L). When setting finishes, you can press 'ENTER' to start known addition measuring. Press 'CANCEL' to quit setting and return to start mode.

From here, if you want to re-calibrate electrode, you can press 'STD' key, meter displays 'Calibrate Electrode?' press 'STD' again to calibrate electrode. (Please refer to chapter 'pX Electrode Calibration' for detailed calibration method).

When all setting finish, you can start measuring of known addition concentration mode which is displayed as follows. In which on upper side of the screen is the present measuring data, at the bottom is the records of measuring result including potential, temperature before and after adding. On right side is 'ENTER', 'END' and 'SETUP' key.


Press 'SETUP' key, you can set manual temperature, select concentration unit, set mV display resolution, pX display resolution which is already been calibrated potential zero point.

Clean the relevant ion select electrode and put in measured solution, meter displays current potential, temperature value and relevant pX value.

Known addition concentration measuring mode

When display stable, press 'ENTER' key, meter stores present potential and displays 'Add Standard Solution', you can add standard solution according to set volume, when display stable again, press 'ENTER', meter displays 'Measuring Finish!' and calculate the concentration value of present sample which is displayed as follows:

Display chart when known adding concentration measuring mode finishes

At this time you can store present measuring result or print the result.

Similarly, if you want to select the concentration unit of sample's concentration, you can also operate according to following procedure: press once at display area of 'sample Cs', the concentration unit will switch for one time.

Press 'END' to finish known addition concentration measuring and return to start mode.

#### 4.6.3.2.1 Blank Concentration Calibration in Known Addition Mode Measuring

Before measurement, if you want to make blank concentration calibration, press 'BLANK' key, select 'Cali Blank Conc' and then press 'ENTER' key, meter displays 'Cali Blank Conc?', press 'ENTER' to start blank concentration calibration.

When meter enters into blank concentration calibration, the display and operation is absolutely the same as 'known addition Measuring Mode'. When measurement finishes, you can press 'END' to return to current parameter setting mode.

#### 4.6.3.2.2 Clear Blank Concentration

If you want to clear previous blank concentration value, you can press 'BLANK' key, select 'Cali Blank Conc?' and then press 'ENTER' to clear blank concentration value to zero.

#### 4.6.3.3 Sample Addition Mode

This mode is similar with 'Known Addition Mode'. The only difference is in standard adding method, to add standard solution into sample and measure the electrode potential changes caused by concentration changing of measured groups and enhance measure the concentration of sample. Similarly, if you add sample into standard solution, you can measure the concentration of sample through measuring the potential changes before and after adding. The formula is as follows:

$$Cx = Cs \times [(1 + \rho) \times 10^{(E2-E1)/S} - \rho]$$

In which:  $Cx^{\sim}$  concentration value of measured sample (adding solution)

C s ∼ concentration value of standard solution


 $\rho$  standard solution volume(Vs)/measured sample volume(Vx)

E 1∼system potential value when measured sample is not added.

E2 system potential when measured sample is added

 $S\sim$  electrode slope

You can select relevant ion mode according to ion electrode used (For example, if you want to measure Na ion concentration, you can select ion mode is Na ion, please refer to chapter 'Setup Ion Mode' for details). Set measuring parameter as ion concentration and select 'Sample Addition Mode' (please refer to chapter 'Setup Measuring Mode), When selection finishes, press 'MEAS' key to start the measuring. Or you can press 'ENTER' to return to start mode. Under start mode, press 'MEAS' key, you can also make concentration measuring which displayed as follows:

Parameter setting in Sample Addition

On upper side of the screen is the adding volume which is the sample's volume going to be added. Volume before adding is the volume of standard solution. 'STD Conc' means the concentration value of standard solution. The relevant concentration unit is listed on right side. The displayed concentration of standard solution is 100mmol/L, At the bottom is 'previous STD'. On right side is 'ENTER', 'CANCEL', and 'STD' key.

This mode needn't measure the blank concentration value of blank standard solution. For detailed operation, please refer to 'Known Adding Measuring Mode'.

#### 4.6.4 GRAN Method


Apart from normal measuring method, you can also use GRAN method to measure the sample with low content. According to mathematical principle of GRAN, you can get sample's concentration with following formula:

$$(Vs + Vx) \times 10^{E/S} = 10^{E0/S} \times (CxVx) + 10^{E0/S} \times (CsVs)$$

When measure, input standard solution' concentration (Cs), Volume (Vs) and measured sample's volume (Vx), then measure sample's electrode potential value when first standard solution is added. Repeat the measuring for three to eight times in sequence, meter can calculate the concentration value of measured solution.

You can select relevant ion mode according to ion electrode used (for example, if you want to measure Ag ion concentration, you can select ion mode is Ag ion, please refer to chapter 'Setup Ion Mode' for details). Set measuring parameter as ion concentration and select GRAN method (please refer to chapter 'Setup Measuring Mode), When selection finishes, press 'MEAS' to start the measuring. Or you can press 'ENTER' to return to start mode. Under start mode, press 'MEAS' key, you can also make concentration measuring

Similarly, GRAN method is the same as 'Known Addition Mode' which needs to set parameter such as addition volume, volume before adding and standard solution concentration, etc... (Please refer to chapter parameter setting of 'Known Addition Mode'). When setting finishes, you can start measuring. The displayed is as follows:


Display chart for GRAN method

On upper side of the screen is present measuring data and lower side is measuring result including adding times and measuring result in sequence. On right side is 'ENTER', 'END' and 'SETUP' key.

Clean the electrode and put it in measured solution, add standard solution according to set volume, meter displays present potential, temperature value and relevant pX value. When display stable, press 'Enter' key, meter will record down and display current measuring result and remind 'Add Standard solution...' you can continue to add standard solution.

During adding, if you press 'SETUP', you can set manual temperature, select concentration unit, set mV, pX display resolution.

Repeatedly add standard solution and measure potential after adding. When you finish adding three times the standard solution, meter automatically displays 'Continue to Add?' which is displayed as follows:

Display chart when GRAN method exceeds 3 times

If you still want add measuring, you can press 'CONTINUE' to continue the adding, otherwise you can press 'END' finish measuring and meter will calculate the concentration value of measured sample. Measuring finishes.

You can press 'END' any time to finish the measuring. If adding time exceeds 3 times, meter will also admit previous adding is valid and calculate the concentration value of measured sample, measuring finishes.

When measuring finishes, you can press 'STORE' to save present result, 'PRINT' to print result, Press 'END' to quit GRAN method and meter returns to start mode.

#### 4.6.4 Conductive Measuring

Meter allows measure conductivity, resistivity, TDS and salinity value. Within all measuring range, it has following switching functions: automatic temperature compensation, automatic calibration, automatic range and automatic frequency etc...

When measure conductivity and TDS, connect temperature electrode, according to set temperature coefficient, meter will automatically compensate conductivity to the value under 25°C. If you don't connect temperature electrode, meter displays the original conductivity of measured solution which is not been compensated.

When measure salinity, connect temperature electrode, meter will automatically compensate salinity to the value under 18°C. If you don't connect temperature electrode, meter displays the original salinity of measured solution which is not been compensated.


# 4.6.4.1 Conductivity Measuring

Please notice before measure conductivity:

- 1. Selection of conductive electrode
  - The conductive electrode with constant 1.0, 10 has two styles: 'bright' and 'platinum black', electrode cast with platinum is usually called platinum black electrode. Better measuring range of bright electrode is  $0 \sim 3000 \mu \text{S/cm}$ . If the range exceeds  $3000 \mu \text{S/cm}$ , measuring error will be big.
- 2. When measure high conductivity, usually adopts conductive electrode with big constant. When conductivity  $\geq 20.00$ mS/cm, you should adopt the electrode with constant 5 or 10. When conductivity  $\geq 100.00$ mS/cm, you should adopt the electrode with constant 10.
- 2. Recommended Conductivity Range and relevant electrode Constant

Conductivity range (µS/cm)	resistivity (Ω·cm)	Recommended electrode constant (cm <sup>-1</sup> )
0.05~2	20M~500K	0.01, 0.1
2~200	500K∼5K	0.1, 1.0
200~2000	5K~500	1.0
2000~20000	500~50	1.0, 10
2000~2×10 <sup>5</sup>	500~5	10

Under start mode, if you have conductivity measuring parameter, you can start measuring directly by pressing 'MEAS' key, otherwise you have to press 'SETUP', select measuring mode and then select conductivity measuring parameter. (**Please refer to 'Setup Measuring Mode' for details**). The display is as follows:

Conductivity measuring status

On upper side of the screen is present measuring mode and system time. Measuring main window displays present conductivity and relevant resistivity, temperature value. On right side is press button 'VIEW', 'STD', 'STORE', 'PRINT' and 'SETUP'. 'Timed Reading Measuring Mode' and 'Balance Reading Measuring Mode' will shield 'STORE', 'PRINT' key automatically.

Press 'VIEW' key to view conductivity measuring parameter, (Please refer to 'Conductivity Measuring Parameter' for details).

Press 'STD' key to re-calibrate conductive constant. (Please refer to 'Conductive Constant Calibration' for details).

Press 'SETUP' key to set manual temperature.

On the title of measuring main window displays present measuring parameter type. The displayed is conductivity. On left side, there is an arrow towards right direction, press this arrow to switch between conductivity, TDS and salinity.

If you select **continue measuring mode** (**Please refer to 'Setup Measuring Mode' for details**), When measuring result stable, you can press 'STORE' to save current measuring result. If it is connected with TP-16 printer before turning on, you can press 'PRINT' key, meter will print present measuring data.

If you select 'Timed Reading Measuring Mode' (Please refer to 'Setup Measuring Mode' for details), meter will automatically measure, calculate and display measuring result. When the time interval meets the time set by user, meter will automatically store measuring result. If it is connected with TP-16 printer before turning on, meter will automatically print measuring data and start next measuring when printing finishes.

If you select 'Balance Reading Measuring Mode' (Please refer to 'Setup Measuring Mode' for details), When measuring starts, meter will automatically measure, calculate and display measuring result. In case the measuring suits the set balance condition, the measuring finishes! At this time, you can press 'STORE' to save current measuring result. If it is connected with TP-16 printer before turning on, you can press 'PRINT' key, meter will print present measuring data. Press 'X' on top right corner and the press 'ENTER' to finish the measuring.


#### 4.6.4.1.1 Conductive Constant Calibration

Each electrode is marked electrode constant when ex-factory. If you suspect the correctness of electrode constant, you can re-calibrate the electrode according to following procedure:

Select suitable standard solution (see chart 1) according to electrode constant. For the compound, please see (chart 2), the relationship between standard solution land conductivity value, (see chart 3).

i. Connect the conductive electrode with meter, turn off temperature electrode (meter does not connect with temperature sensor), meter regards manual temperature as present temperature value, set manual temperature as 25°C, the conductivity vale displayed now is the absolute conductivity which is not be temperature compensated.

- ii. Wash the conductive electrode with distilled water.
- iii. Insert conductive electrode into standard solution.
- iv. Control solution's temperature as  $(25.0 \pm 0.1)\Box$  or  $(20.0 \pm 0.1)\Box$  or  $(18.0 \pm 0.1)\Box$  or  $(15.0 \pm 0.1)\Box$
- v. Press 'STD' key, select 'Cell Constant' and confirm to enter into cell constant standardizing status which is displayed as follows:

Cell constant standardization

- vi. Press 'SETUP' button, input the relevant figure in chart 3, i.e. present standard solution's conductivity.
- vii. When reading figure stable, you can press 'ENTER' key, meter will automatically calculate new electrode constant value, calibration finishes. Press 'END' key, meter stops the electrode constant calibration.

Chart 1 KCL standard solution for measuring cell constant

Electrode constant (l/cm)	0.01	0.1	1	10
KCL solution approximate	0.001	0.01	0.01or 0.1	0.1or 1
concentration (mol/L)				

**Chart 2 Compound of standard solution** 

Approximate	Volume concentration KCL(g/L) solution
concentration (mol/L)	(20□ in air)
1	74.2650
0.1	7.4365
0.01	0.7440
0.001	dilute100mL 0.01mol/Lsolution to 1L

Chart 3 Relationship between KCL solution approximate concentration and conductivity

appro. Temp.	15.0□	18.0□	20.0□	25.0□	30.0□
conc. (mol/L)					
1	12120	97800	101700	111310	131100
0.1	10455	11163	11644	12852	15353
0.01	1141.4	1220.0	1273.7	1408.3	1687.6
0.001	118.5	126.7	132.2	146.6	176.5

#### 4.6.4.2 TDS Measuring

When measure high TDS, conductive electrode with big cell constant is usually be adopted, When TDS  $\geq 10.00$ g/L, you must adopt the electrode with electrode constant 5 or 10.

Before measure TDS, you must refer to chapter 'Selection of Conductive Electrode' to confirm the conductive electrode to be used in measuring, set suitable electrode constant, temperature compensation coefficient and TDS switching coefficient. If you need to calibrate the electrode, you should also calibrate electrode constant under conductivity measuring parameter. You should note that under TDS measuring parameter, you can only calibrate TDS switching coefficient.

Under start mode, if you have already selected conductive measuring mode parameter, such as conductivity parameter, you can start measuring directly. Otherwise, you should enter into 'Setup Measuring Mode', select conductive measuring unit to measure. (Please refer to 'Setup Measuring Mode' for details).

For the display, measuring procedure and operating, please refer to chapter 'conductivity measuring'.

#### 4.6.4.1.1 Calibrate TDS Factor

Select suitable standard solution according to the character and measuring range of measured solution. For the relationship between TDS standard solution and conductivity, please see chart 4.

- a. Press '>'beside main measuring window, switch measuring status to conductivity measuring. Under conductivity measuring status, set conductive electrode constant properly through view conductive using parameter or re-calibrate conductive electrode constant through calibrate conductive electrode. When it is finished, press '>' to enter into TDS measuring mode.
- b. Wash the conductive electrode with distilled water.
- c. Insert conductive electrode into standard solution. Control solution's temperature as  $(25.0 \pm 0.1)\Box$
- d. Press 'SETUP' key, input relevant figure in chart 4, i.e. present standard solution's TDS value.
- e. When reading figure stable, you can press 'ENTER' key, meter will automatically calculate new TDS factor, calibration finishes. Press 'END', meter stops the calibration of TDS factor.

Chart 4 Relationship between conductivity and TDS standard solution

conductivity	TDS Factor					
μS/cm	KCl	(mg/L)	NaCl	(mg/L)	442	(mg/L)
23	11.6		10.7		14.74	
84	40.38		38.04		50.5	
447	225.6		215.5		300	
1413	744.7		702.1		1000	
1500	757.1		737.1		1050	
2070	1045		1041		1500	
2764	1382		1414.8		2062.	7
8974	5101		4487		7608	
12880	7447		7230		11367	
15000	8759		8532		13455	
80000	52168		48384		79688	

- 1. 442 stands for 40% Na2SO4 \ 40% NaHCO3 \ 20% NaCl \.
- 2. The value listed in above chart is the value under  $25 \square$

#### 4.6.4.3 Salinity Measuring

Before measure salinity, you must refer to chapter 'Selection of Conductive Electrode' to confirm the conductive electrode to be used in measuring, set suitable electrode constant.

When measure salinity, we usually adopt conductive electrode with constant 10. When salinity  $\geq 1.00\%$ , you must adopt the electrode with electrode constant 10.

Under start mode, if you have already selected conductive measuring mode parameter, such as conductivity parameter, you can start measuring directly. Otherwise, you should enter into 'Setup Measuring Mode', select conductive measuring unit to measure. (Please refer to 'Set Measuring Mode' for details).

For the display, measuring procedure and operating, please refer to chapter 'conductivity measuring'.


#### 4.6.4 Dissolved Oxygen Measuring

Meter could measure dissolved oxygen concentration, dissolved oxygen saturation and electrode current. Meter has automatic temperature compensation function and could make zero oxygen calibration, full scale calibration, atmosphere calibration and salinity calibration.

When meter is used for first time, clean the dissolved oxygen electrode with distilled water an insert into measured solution. When meter is turned on, you can make the measurement. Meter calculates, displays dissolved oxygen concentration (or saturation) and electrode current value.

If this is your first time using or long time not using the meter, please calibrate dissolved oxygen electrode first, you can refer to chapter Dissolved Oxygen Calibration for details.

In 'Setup Measuring Mode', select dissolved oxygen parameter such as dissolved oxygen, saturation etc... to start measuring of DO which is displayed as follows:

Display chart for dissolved oxygen measuring status

On upper side of the screen is present measuring data, system time, DO measuring result and relevant current DO value, current temperature value. On right side is press key including 'VIEW', 'STD', 'STORE', 'PRINT' and 'SETUP'.

Press 'VIEW' to view conductivity measuring parameter, set current air pressure and salinity etc... (Please refer to 'View solution's Measuring Parameter' for details).

Press 'STD' to re-calibrate DO's zero oxygen, full scale etc... (Please refer to 'Dissolved Electrode Calibration' for details).

During measurement you can press 'SETUP' to set manual temperature.

If you select **Continue Mode (Please refer to 'Setup Measuring Mode' for details),** When measuring result stable, you can press 'STORE' to save current measuring result. If meter is connected with TP-16 printer before turned on, you can press 'PRINT' key, meter will print present measuring data.

If you select Timed Reading Measuring Mode (Please refer to 'Setup Measuring Mode' for details), meter will automatically measure, calculate and display measuring result. When the time interval meets the time set by user, meter will automatically store measuring result. If meter is connected with TP-16 printer before turning on, meter will automatically print measuring data and start next measuring when printing finishes.


If you select 'Balance Reading measuring mode' (Please refer to 'Setup Measuring Mode' for details), When measuring starts, meter will automatically measure, calculate and display measuring result. In case the measuring suits the set balance condition, the measuring finishes! At this time, you can press 'STORE' to save current measuring result. If meter is connected with TP-16 printer before turning on, you can press 'PRINT' key, meter will print present measuring data.

Press 'X' on top right corner and the press 'ENTER' to finish the measuring.

#### 4.6.4.1 Dissolved Oxygen Calibration

In order to get precise measuring result, dissolved oxygen electrode should be calibrated before measuring. Meter has several calibration functions such as zero oxygen calibration and full scale calibration.

Under DO measuring mode, press 'STD' key, select 'DO electrode' and confirm. You can start to calibrate dissolved oxygen electrode which is displayed as follows:

Display chart for dissolved oxygen calibration

On upper side of the screen are present DO electrode current, STD value and current temperature value. In the middle are present Air Pressure and Salinity Factor. At the bottom is STD result. On right side are press keys including 'Zero STD', 'Full STD' and 'SETUP' for setting air pressure and salinity.

#### 4.6.4.1.1 Zero Oxygen Calibration

Clean DO oxygen electrode with distilled water and insert into 5% fresh Na<sub>2</sub>SO<sub>3</sub> solution.

When reading figure stable, press 'Zero STD' and then press 'ENTER', meter will automatically record zero oxygen value, zero oxygen calibration finishes.

#### 4.6.4.1.2 Full Scale Calibration

Take out dissolved oxygen electrode from solution, clean with water and dry it with filter paper, put into vessel with distilled water (such as tri-angel beaker, high leg beaker), place in the air close to water or place in air and make sure there is no water on the surface of electrode.

Press 'STD' key, meter enters into dissolved oxygen calibration status.

When reading figure stable, press 'Full STD' then press 'ENTER', meter will automatically record full scale value, full scale calibration finishes.

#### Generally, it is not necessary to set salinity. The acquiescent value is 0.0g/L

There are two ways to set DO salinity. One is: when calibrate DO electrode, press 'SETUP' key beside salinity, meter shows an input window, you can input salinity directly. Another is: In view using parameter of DO, you can also set salinity.

#### 4.6.4.1.3 Setup Air Pressure

The DO value measured by meter is in correspondence to atmosphere. The acquiescent air pressure is 101.3kPa. Before measuring, please select suitable air pressure.

Generally, it is not necessary to set air pressure. The acquiescent value is 103.1kPa.

There are also two ways to set DO air pressure. One is: when calibrate DO electrode, press 'SETUP' key beside air pressure, meter shows an input window, you can input air pressure directly. Another is: In view using parameter of DO, you can also set air pressure.

# 5. Maintenance

- 5.1. Meter must have good grounding and should avoid corrosive gas.
- 5.2. Before turning on meter, you must check if power is connected properly.
- 5.3. Connect power apparatus and turn on the mete. If the screen is not bright, please check if power apparatus is working or if you utilize the original apparatus. Required power apparatus is: 9VDC, above 800mA, inside is positive and outside is negative.
- 5.4. The socket should be kept clean, dry and could not be touched with acid, alkali, salty solution.
- 5.5. Meter can be used stably for long time. After measuring the sample, all electrodes should be immersed in distilled water.
- 5.6. When meter is not be used, please connect it with short circuit to avoid any damage of the meter.
- 5.7. For the measurement of ion select electrode, please refer to relevant material and obey strictly.
- 5.8. Improper use of conductive electrode will cause abnormal working condition of meter. Please keep conductive electrode immersed in solution totally. Please take care that the place to install conductive electrode should be: avoid to be installed in 'dead' corner, it should be installed in the place with good water cycling.
- 5.9. For the measuring of high pure water, it should be done in close and flow status. The direction of water flow is water can flow into but the flow shouldn't be too speedy.
- 5.10. If conductivity exceeds 3000μS/cm, bright electrode could not measure correctly. At this time, please exchange to platinum black electrode to measure.
- 5.11. When DO electrode is not be used, please keep the electrode in boiled and cooled distilled water. Please avoid immerse electrode into Na<sub>2</sub>SO<sub>3</sub> solution. In case the electrode is soaked with Na<sub>2</sub>SO<sub>3</sub> solution, meter will be damaged.
- 5.12. When the output of DO electrode which is newly installed with electrolyte and membrane is

low and not able to calibrate, if connect DO electrode to meter and turn on meter for 5 minutes, meter's displayed reading figure still not reach your needed value, there are two possibility: One is: membrane is not stick tightly to golden negative pole. Another is: the surface of golden negative pole is not wet. You can knock the oxygen electrode on table or chair. If reading figure increase, it means the function of electrode recovered. Otherwise, you have to exchange membrane.

- 5.13. After 5 minutes connecting with electric and be polarized, if the zero oxygen index of DO electrode is higher than technical condition, it is probably caused by the damage of negative pole. Check if there is hollow and holes on the surface of golden negative pole and also check if the surrounding of golden negative pole is disconnected with basement.
- 5.14. If all above conditions are not exist but meter still do not work, please contact with our company.

# 6. Complete Set of Meter

6.1. MP-2006 Multi-parameter Analyzer	1
6.2. DJS-1C Platinum Black Electrode	1
6.3. E-201-C-9 Complex Electrode	1
6.4. DO-958-S Dissolved Oxygen Electrode	1
6.5. T-818-B-6 Temperature Sensor	1
6.6. REX-2 Multi-function electrode supporter	1
6.7. Short circuit plug for electrode input (inserted in pH electrode socket wh	nen ex factory.)
6.8. pH standard buffer solution	5packages for each
6.9. Power apparatus (9VDC,800mA, inside – positive, outside - negative)	1
6.10.Quality Certification	1
6.11 One set of accessories please refer to packing list	

# 7. Appendix

#### 7.1 Compound of Fluorine ion solution

# • Standard Solution:

Precisely weigh 4.20g dissolve it in distilled water, add water to 1000mL and store in plastic bottle. This solution is  $1\times10^{-1}$  mol/L  $\,$  F $^ _{\circ}$ 

 $1\times10^{-2}$  mol/L F<sup>-</sup>: take  $1\times10^{-1}$  mol/L F<sup>-</sup> solution 100ml, dilute it to 1000ml;  $1\times10^{-3}$  mol/L F<sup>-</sup>: take  $1\times10^{-2}$  mol/L F<sup>-</sup> solution 100ml, dilute it to 1000ml;

#### Total Ion Strength Adjustment Agent (TISAB)

Weigh  $58.8g \text{ Na}_3\text{C}_6\text{H}_5\text{O}_7 \cdot 2\text{H}_2\text{O}$  and  $85g \text{ Na}_3\text{NO}_3$ , dissolve these in distilled water, adjust the pH value of solution to  $5 \sim 6$ , then dilute the solution to 1000mL

# Solution for Testing

pF2( $1\times10^{-2}$  mol/L F<sup>-</sup>): take  $1\times10^{-1}$  mol/L F<sup>-</sup> solution 10ml, add total ion strength adjustment agent (TISAB) 20ml, dilute the solution with distilled water to 100mL pF3( $1\times10^{-3}$  mol/L F<sup>-</sup>): take  $1\times10^{-2}$  mol/L F<sup>-</sup> solution 10ml, add total ion strength adjustment agent (TISAB) 20ml, dilute the solution with distilled water to 100mL pF4( $1\times10^{-4}$  mol/L F<sup>-</sup>): take  $1\times10^{-3}$  mol/L F<sup>-</sup> solution 10ml, add total ion strength adjustment agent (TISAB) 20ml, dilute the solution with distilled water to 100mL

# 7.2 Relationship Chart between DO Concentration and Temperature

$Temp(\Box)$	C <sub>S</sub> (mg/L)	$\Delta C_{\rm S}({\rm mg/L})$	$Temp(\square)$	C <sub>S</sub> (mg/L)	$\Delta C_{\rm S}({\rm mg/L})$
0.0	14.46	0.0925	21.0	8.90	0.0467
1.0	14.22	0.0890	22.0	8.73	0.0453
2.0	13.82	0.0857	23.0	8.57	0.0440
3.0	13.44	0.0827	24.0	8.41	0.0427
4.0	13.09	0.0798	25.0	8.25	0.0415
5.0	12.74	0.0771	26.0	8.11	0.0404
6.0	12.42	0.0745	27.0	7.96	0.0393
7.0	12.11	0.0720	28.0	7.82	0.0382
8.0	11.81	0.0697	29.0	7.69	0.0372
9.0	11.53	0.0675	30.0	7.56	0.0302
10.0	11.26	0.0653	31.0	7.43	
11.0	11.01	0.0633	32.0	7.30	
12.0	10.77	0.0614	33.0	7.18	
13.0	10.53	0.0595	34.0	7.07	
14.0	10.30	0.0577	35.0	6.95	
15.0	10.08	0.0559	36.0	6.84	
16.0	9.86	0.0543	37.0	6.73	
17.0	9.66	0.0527	38.0	6.63	
18.0	9.46	0.0511	39.0	6.53	
19.0	9.27	0.0496	40.0	6.41	
20.0	9.08	0.0481			

This chart is digested from ISO5813-1983 water analyze – measurement of dissolved oxygen in water – iodine titration

#### 7.3 Notice for Ordering

- a. When meter ex factory, the attached electrode is DSJ-1C platinum black electrode (electrode constant is 1.00)
- b. According to measuring range of conductivity, when user place order, please select suitable conductive electrode in reference to conductivity measuring chapter in this instructions.

# Packing List

# MP-2006

# Multi-Parameter Analyzer

No	Name and Specification	Amount
1	MP-2006 Multi-Parameter Analyzer	1
2	DJS-1D Platinum Black Electrode	1
3	65-1C pH Combination Electrode	1
4	DO-958-S Dissolved Oxygen Electrode	1
5	T-818-B-6 Temperature Sensor	1
6	Multi-function Electrode Support Stand	1
7	pH Standard buffer solution pH4,7,10	5bags for each
8	Universal power supply	
	(9V DC, 800mA, inside- positive and outside- negative )	
9	TP-16 Print Connecting Line	1
10	RS-232 Communication Line	1
11	Q9 Short Circuit Plug	1
12	Instruction Manual	1