FRIGEDOR REG

3001214

**UNIDAD REFRIGERADORA PARA BAÑOS
IMMERSION COOLERS FOR BATHS**

INDICE

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GENERAL INFORMATION

- 1) Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.
- 2) Do not install or use the equipment without reading this handbook before.
- 3) This handbook must always be attached to the equipment and it must be available for all users.
- 4) If you have any doubts or enquiries, please contact with your supplier or J.P. Selecta's technical service.
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY EQUIPMENT TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 89/655/CEE, the manufacturer is not responsible for the damage that can occur.
- 7) Do not use the equipment with liquids which can give off vapours capable of making explosive mixtures.

PACKING LIST

The standard equipment consists of the following components:

Description	Code	Quantity
FRIGEDOR REG	3001214	1
Pt100 probe	43068	1
Instruction Manual	80040	1

ACCESSORIES

Insulated tanks: Thermally insulated to withstand low temperatures. Insulated with injected foam of expanded-polyurethane to avoid condensations between -40°C and +90°C.

Code	Capacity litres	Tank dimensions cm			Overall dimensions cm			Kg.
		Depth	Length	Width	Height	Length	Width	
3001214	9	15	29	22	24	40	34	6

TECHNICAL FEATURES

Voltage supply 115-230V 50/60 Hz according to the characteristics plate indications.

Code		3001214
Regulation range °C		-20 / +20
Overall dimensions (cm)	Height	41
	Length	21
	Width	34
Cooling power (W)		at 20°C 276
Power (HP)		1/5
Consumption (W)		285
Weight (Kg)		14

EQUIPMENT DESCRIPTION

The immersion cooler for baths FRIGEDOR REG are designed for temperatures up to -20°C.

Desktop design with a sealed compressor, a ventilated condenser, a spiral cooling coil, with temperature regulator and a Pt100 probe.

COIL DIMENSIONS:

- * Line length: 900 mm.
- * Coil length: 150 mm.
- * Ø coil: 45 mm.

INSTALLATION

Place the equipment near a suitable mains connection.

CAUTION!!! IMPORTANT FOR YOUR SAFETY

Be sure that the voltage supply is the same as the one indicated on the characteristic plate of the equipment.

Do not use the unit to cool highly inflammable liquids.

Do not use the equipment if it is NOT earthed.

If you have to change the plug, bear in mind the following:

Blue cable:	Neutral.
Brown cable:	Phase.
Yellow/green cable:	Earth.

MAINTENANCE

Before removing the side cover, disconnect the equipment from the mains.



The guarantee will not have any effect if the equipment is manipulated by somebody who has not been expressly authorized by J.P.SELECTA, s.a.u.

CLEANING:

Remove the metallic casing and clean the condenser wings with the help of a brusher or compressed air. It is recommended to do it periodically.

OPERATION

CONTROL PANEL:

1. Main switch
2. Temperature regulator
3. Temperature probe plug
4. Pt100 temperature probe
5. Cooling coil
6. Key to scroll parameters
7. Keys to modify parameters
8. Key to exit from the menus
9. Running compressor led
10. Set temperature display
11. Current temperature display

CONTROL PANEL:

- *Place the refrigerator coils inside the liquid container that must be cool down.*
- *Connect the Pt100 temperature probe (4) in the connector (3).*
- *Place the probe inside the test liquid. It is preferable that it is not in contact with the bucket or the coil. You can use a bracket and a clamp to do this.*
- *Activate the control panel switch with operating light indicator (green).*
- *Select the operating temperature with the keys (7). The temperature selection is shown in the display (10).*

TEMPERATURE CORRECTION:

You can correct the temperature displayed on the screen if contrasted with a more accurate thermometer.

- *Press the key (6) scroll parameters. The display (11) will show "SHIFT" and the display (10) will show the correction value in grades whose default value is "0".*
- *Press keys (7) to modify the value.*
- *Press key (6) or (8) to validate and exit from the menu.*

Example 1:

If the standard thermometer indicates 10°C while the Frigedor unit temperature control indicates 9°C. The value to be entered is 1.

Example 2:

If the standard thermometer indicates 10°C while the Frigedor unit temperature control indicates 11°C. The value to be entered is -1.

Note: The temperature controller is configured to display tenths of a degree, as the compressor timing system ON / OFF does not guarantee a greater precision.

GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use of causes beyond the control of J.P. SELECTA, s.a.u.

Any manipulation of the equipment by unauthorized personnel by J.P. SELECTA, s.a.u. cancels the guarantee automatically.

SPARE PARTS

To guarantee the safety of the equipment, the spare parts must be purchased from J.P.SELECTA, s.a.u.

Description	Code
Schuko cable	7002
Handle	15490
Switch	20017
Cooling unit set	24040
Jebe clamps	27515
Cables Top	28004
Pt100 probe	43071
Temperature regulator	16167