ESTUFA REFRIGERADA HOT-COLD S

**REFRIGERATED CABINET HOT-COLD S**

2101618

# Indice

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# Index

<b>1 Packing list</b> .....	<b>12</b>
1.1 Standard package content .....	12
<b>2 Overview</b> .....	<b>12</b>
<b>3 Technical specification</b> .....	<b>12</b>
<b>4 Safety</b> .....	<b>12</b>
4.1 Safety marks .....	12
4.2 Important information mark .....	13
4.3 Personnel training .....	13
4.4 Proper uses of the equipment .....	13
4.5 Modifications .....	13
<b>5 Installation</b> .....	<b>13</b>
5.1 Suitable location .....	13
5.2 Connection to power supply .....	13
<b>6 Operation</b> .....	<b>14</b>
6.1 Status indication on the control panel .....	14
6.2 Set temperature .....	14
6.3 Set working time .....	15
6.4 Setting the temperature control .....	15
6.5 Setting the «defrost» operation .....	17
6.6 Setting the overtemperature alarm .....	18
6.7 Setting time .....	18
6.8 Power-Off recovery function .....	18
6.9 Adjustment of the "Offset" temperature reading .....	19
<b>7 Warranty</b> .....	<b>19</b>

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# 1 Packing list

## 1.1 Standard package content

- Refrigerated cabinet
- 2x try
- Power cord
- 4x Fuse
- Instruction manual

## 2 Overview

El HOT-COLD S is a refrigerated cabinet with forced air circulation. Suitable for use in pharmaceuticals, agriculture research, ... Its working range goes from 5°C to 65°C and it fits requirements for plants growing, culture, ferment, DBO test, etc...

## 3 Technical specification

Mains voltage:	230V 50/60Hz (See plate) Other voltage available (Ask)
Electrical power consumption:	700W (See plate)
Temperature range:	5°C ... 65°C
Display resolution:	0.1°C
Temperature precision:	+0.2°C (@37°C)
Temperature fluctuation:	±0.5 (Time variation)
Temperature uniformity:	From +0.3 to 2 (space variation)
Room ambient temperature:	5°C ... 35°C
Timed operation range:	1 minute to 100 hours
Power cord:	Phase, Neutral, Earth (section 0.75mm <sup>2</sup> )
Weight (Net):	70 Kg
Dimension int.: Wide x High x Deep	500 x 650 x 500mm
Dimension ext.: Wide x High x Deep	630 x 1280 x 630mm
Capacity:	160 Liters

## 4 Safety

For a safety use of the equipment, follow these instructions:

### 4.1 Safety marks

Safety marks identify hazard situations

#### Electrical hazard

- Grounding outlet to avoid accidental electric shock or fire disaster.
- Do not attempt to repair the product yourself. Do not open the power box to avoid electric shock.
- Do not pull out the plug when it is operating.
- Do not drag the wire to plug the outlet.
- When replacing the fuse. Replacing the fuse without cutting off the power supply will probably result in an electric shock.

#### 4.2 Important information mark

- Information important to obtain good results or for an optimal operation of the equipment.
- Information important to extend the life of the equipment or to avoid the degradation of some of its components.

#### 4.3 Personnel training

This equipment must only be operated by personnel trained and suitable qualified to work in a chemical laboratory and its common hazards.

This equipment must only be operated by personnel who have read and understood this manual or have been qualified on its manipulation and operation.

#### 4.4 Proper uses of the equipment



This equipment is for indoor laboratory use.

The improper use of this unit (not following manual instructions) can endanger the protection prepared by the equipment.

This equipment is not designed and not manufactured to work in an explosive ambient area.

#### 4.5 Modifications

The operation modification or the manipulation of the safety systems unauthorized by the manufacturer, could expose the operator to hazards not indicated in this manual.

## 5 Installation

Once the unit is unpacked, visually check that all the components are not damaged. Keep the original package for some days.

### 5.1 Suitable location

Place the unit in a flat, levelled and stable surface and suitable for the unit weight.

No longer that 1.5m from a water tap, waste output and power supply plug.

To work comfortably, let a free space of 50cm on each side and not least than 20cm on the rear part.

Do not expose the unit directly to sun light

Do not place the unit directly to air flow

### 5.2 Connection to power supply

Locate a power supply socket (rated to 16A) near to the unit.

Use the power cord supplied or another of similar specification.

For safety reasons, the socket must have ground terminal, since the unit must work grounded.

Before to plug it on the socket, check that its voltage match the voltage marked on the unit plate.

Do not manipulate the power cord with wet hands.

## 6 Operation

The control panel indicates the unit status and allows configure its on working mode.

Mains switch is on the right hand side wall. Press to "ON" position to power on the unit

Beside the mains switch there the lights ON/OFF switch. (Tube)

The unit can operate in a Continuous mode if no working time is entered, or in a TIMED mode if a working is programmed (See 5.3).

In order to optimize the unit performance, and get a precise temperature control, 3 working modes can be configured: (See 5.4)

- Refrigeration continuously runs.
- Refrigeration intermittently runs.
- Refrigeration stopped.

### 6.1 Status indication on the control panel

ALARM

Overtemperature alarm. Actual temperature value is above temperature set.

HEATING

Heater working.

COOLING

Refrigeration working.

°C

Display shows temperature °C


H.m.

Display shows time hh.mm

### 6.2 Set temperature






- Press  Lamp °C lights flashing.  
Display shows the preset temperature.

- Press  or  To change the temperature set.

**Note:** Keep press for 5s to accelerate changing rate.

- Press  to store the changes.  
Display shows the actual temperature on the oven.

### 6.3 Set working time

- Press  Lamp hh.mm lights flashing.  
Displays shows time (format 24h).
  - Press  Displays shows the preset working time.
  - Press  or  To change working time.
- Note:** Keep press for 5s to accelerate changing rate.
- Press  Display shows the time when the oven will stops  
(Format 24H hh.mm)
  - Press  To store the changed value.  
Display shows working time


**Note:**

- Max working time is 99 hours and 59 minutes.
- For a **continuous operation** enter 00.00
- When the unit works timed (not continuous) the decimal point flashes on the display.

### 6.4 Setting the temperature control

In order to get the optimal unit performance, temperature control must be properly set.

#### 6.4.1 Refrigeration continuously working

- Suitable to work from 4°C to 10°C
  - Temperature precision is  $\pm 0.1^\circ\text{C}$
  - Compressor is continuously working.
- Press  Keep 4s until «0» appears on the display.
  - Press  To set value «5».
  - Press  To validate the input.  
On the display appears «SCL» and alternately the refrigeration coefficient.
  - Press  or  to set value «0.5»
  - Press  to store the changes.




**Important note:**

On this working mode, automatic **defrost must be activated.**

### 6.4.2 Refrigeration working intermittently

- Suitable to work from 10°C to ambient +5°C
- Temperature precision is  $\pm 0.5^{\circ}\text{C}$
- Compressor works in a 2 minutes interval.
- Press  Keep 4s until «0» appears on the display
- Press  to set value «5».
- Press  To validate the input.  
On the display appears «SCL» and alternately the refrigeration coefficient.
- Press  | 0 |  to set value «0.0».
- Press  to store the changes.

### 6.4.3 Working with refrigeration stopped

- Suitable to work from ambient +5°C to 65°C
- Temperature precision is  $\pm 0.1^{\circ}\text{C}$
- Compressor remains stopped.
- Press  Keep 4s until «0» appears on the display
- Press  to set value «5».
- Press  To validate the input.  
On the display appears «SCL» and, alternately the refrigeration coefficient.
- Press  | 0 |  to set value «3»
- Press  to store the changes.

**Note:**

The working mode affects the precision of temperature





## 6.5 Setting the «defrost» operation

### 6.5.1 Terminology used

- **Defrosting period:**  
Time between two consecutive defrosting phases.
- **Defrosting time:**  
Duration of a defrosting action.

### 6.5.2 Configuring the defrost operation


- Press  Keep 4s until «0» appears on the display
- Press  To set value «5».
- Press  To validate the input.  
On the display appears «SHS» and alternately the preset defrosting period.
- Press  |  Set the defrost period. (minutes)
- Press  To validate the input.  
On the display appears «SHT» and alternately the preset defrosting time.
- Press  |  Set the defrost time (minutes)
- Press  to store the changes.

#### Notes about the defrosting settings:

- Because of the variance of the ambient temperature, the setting temperature and the time of continuous refrigeration, the defrosting timing and defrosting time should be set according to the accumulation of frost inside the evaporation chamber.
- If defrosting operation is end and the accumulation of frost inside the evaporation chamber is still severe, shorten the defrosting timing or extend the defrosting time appropriately.
- When the equipment works continuously under low temperature for a long period of time, the speed of frost accumulation will accelerate gradually because of the humidity increase inside the evaporation chamber. So, shorten the defrosting timing or extend the defrosting time appropriately to solve this problem.





#### Important suggestion:

For those equipment working under low temperature for a long period of time, it is good for the equipment to be dried thoroughly every 10 days.

To do it, set the refrigeration coefficient to be 3 (compressor stopped), set the temperature to be the highest value (65°C), and leave it working with the door open for 4 hours.

## 6.6 Setting the overtemperature alarm

- Press  Keep 4s until «0» appears on the display
- Press  To set value «3».
- Press  To validate the input  
On the display appears «AHL» and alternately the **upper** alarm temperature.
- Press  |  To set the **upper** alarm temperature.
- Press  To validate the input  
On the display appears «ALL» and alternately the **lower** alarm temperature.
- Press  |  To set the **lower** alarm temperature.
- Press  to store the changes.

### Overtemperature alarm effect:

When the deviation between the real temperature and the preset one is beyond the upper or lower limit value:




- The audio and visual alarm will be activated.
- The red alarm indicator turns on and the buzzer keeps alarming.
- The heating or refrigerating output will continue.

### Attention:

- The upper limit alarm value must be higher than the temperature set.
- The lower limit alarm value must be lower than the temperature set.

Or else, the temperature exceeding alarm will be always activated.

## 6.7 Setting time

- Press  Keep 4s until «0» appears on the display.
- Press  to set value «4».
- Press  To validate the input.  
On the display appears «Hr» and alternately the preset hour.
- Press  |  to set the hour.
- Press  To validate the input.  
On the display appears «Mu» and alternately the preset minute.
- Press  |  to set the minutes.
- Press  to store the changes.


## 6.8 Power-Off recovery function

The HOT-COLD S is provided with power-off recovery function.

When the external power is recovered after a sudden cut off, if the temperature drop inside is not bigger than 3, the equipment can automatically recover to run according to the originally designed program.

No need to set this function.

### 6.9 Adjustment of the "Offset" temperature reading

- Press  Keep until «0» appears on the display.
- Press  to set value «3».
- Press  To validate the input.  
On the display appears «CPO» and alternately the offset value.
- Press  or  to set correction value.
- Press  To validate the input.

## 7 Warranty

The warranty covers the defective components or faulty assembling of the unit for 1 year since the purchasing date.

Keep purchasing bill as the only valid document for the purchasing date.

The warranty does not cover damages caused for a unit malfunction. The unit must be suitably qualified before being used for valuable samples processing.

Return good are not accepted without previous J.P. SELECTA, s.a.u. approval.

J.P: SELECTA, s.a.

---

---