



WEIGHING DISHES "ACCU-WEIGH"

Anti static. Reproducible tare weighing.

Contamination free. Large selection of sizes.

Common features:

Made from high density polystyrene with flat base, round edge, contamination free surface and hydrophilic. Flexible in form for easy manipulation of the sample either liquid or solid.

Manufactured in two forms:

- Diamond form, colours white and black.
- Square form, colour white.

Part No.	Form	Capacity	Height / Width / Depth			Weight	Colour
5726004	Diamond	5	3.8	31	54	0.25	White
5726005	Diamond	5	3.8	31	54	0.25	Black
5726006	Diamond	30	13.5	55	85	1.00	White
5726007	Diamond	30	13.5	55	85	1.00	Black
5726008	Diamond	100	20	92	132	2.6	White
5726009	Diamond	100	20	92	132	2.6	Black
5726010	Square	7	8.5	44	44	0.65	White
5726011	Square	100	24	80	80	2.2	White
5726012	Square	250	24	134	134	6.1	White

Comes in packets of 100 units.