

## FEATURES

Fixed volume pipettes from 1 to 1000 µl, model dependent.  
Simple and easy to use.  
Discard tip ejection for all models.  
Comes complete with calibration tool.

MODELS	Part No.	Capacity µl	Accuracy %	Precision %
Dispen-Fix	5904009	5	±2.0	< 0.8
Dispen-Fix	5904010	10	±1.6	< 0.5
Dispen-Fix	5904011	20	±1.0	< 0.3
Dispen-Fix	5904012	25	±1.0	< 0.3
Dispen-Fix	5904013	50	±0.8	< 0.3
Dispen-Fix	5904014	100	±0.8	< 0.25
Dispen-Fix	5904015	200	±0.8	< 0.15
Dispen-Fix	5904017	500	±0.8	< 0.2
Dispen-Fix	5904018	1000	±0.8	< 0.13

# Variable volume dispensing pipettes "Dispen-Var"

WITH VALIDATION TEST ACCORDING TO SERIAL NUMBER.

AUTOCLAVABLE UP TO 121 °C.

## FEATURES

Variable volume from 1 to 5000 µl, model dependent.  
Simple and easy to use adjustable dispense volume.  
Digital set volume selection.  
Discard tip ejection for all models up to 1000 µl.  
Comes complete with calibration tool.


MODELS	Part No.	Capacity µl	Divisions µl	Precision %
Dispen-Var	5900010	1 - 10	0.1	< 1.5 - < 0.4
Dispen-Var	5900020	2 - 20	0.1	< 1.5 - < 0.8
Dispen-Var	5900100	10 - 100	1	< 1.0 - < 0.3
Dispen-Var	5900200	20 - 200	1	< 1.0 - < 0.5
Dispen-Var	5901000	100 - 1000	10	< 0.8 - < 0.28
Dispen-Var	5905006	1000 - 5000	10	< 0.4 - < 0.2



NEW DESIGN

## ACCESSORIES. PIPETTE TIPS AND RACKS. Accurate and precise, made of polypropylene.

Models N° image	Part No.	Range µl	Colour	Qty. tips x part No.	Suitable for pipette part No.	Volumes pipette µl
1	Tips 5490040	0.1- 10	Transparent	1000	5904009	5
	Rack 5490035	0.1- 10	Transparent	96	5904010	10
					5900010	1 - 10
2	Tips 5490041	2 - 20	Transparente	1000	5900020	2 - 20
	Rack 5490036	2 - 20	Transparente	96		
3					5904011	20
					5904012	25
	Tips 5490042	5 - 200	Transparent	1000	5904013	50
	Rack 5490037	5 - 200	Transparent	96	5904014	100
					5904015	200
					5900100	10 - 100
4	5490044	101 - 1000	Blue	1000	5900200	20 - 200
	5490034	101 - 1000	Blue	96	5904017	500
					5904018	1000
5					5901000	100 - 1000
	5490047	1000 - 5000	Blue	500	5905006	1000 - 5000
	5490039	1000 - 5000	Blue	500		



**Pipette rack.**  
Made of AISI 304 stainless steel. No. of places: 6 dosing pipettes up to 28 cm in length.  
Part No. 1001215



Rack for pipette tips.