



Microtiter Strip Reader "Reader M-2000"

STORAGE MEMORY FOR UP TO 59 METHODES INCLUDING: FACTORS, STANDARDS OR CALIBRATION CURVES WITH UP TO 10 STANDARDS. ADDITIONAL CALCULATION FEATURES.

APPLICATIONS

Clinical diagnostics, E.I.A., Veterinary, quality control.

FEATURES

Spectral range: 400 to 700 nm.

interference filter of 10 nm band pass.

Supplied complete with a filter of 450 nm (other wave-lengths available consult for options.)

Measuring range: -0.2 to 3.5 O.D.

Photometric accuracy: >1%.

Photometric precision: ±1%.

Light source: Long life LED.

Detector: solid state.

59 method memories.

Curve point editing feature.

Software: ABS. Concentration and polynomial curve fitting with up to 10 standards.

RS-232 output to connect a printer or to a PC.

Interface RS-232.



MODEL

Part No.	Height / Width / Depth cm			Power W	Weight Kg
4120013	11	18	28	7	2

ACCESSORY

Ink printer (not thermal). Paper 2 1/4" (56 mm) wide roll. Interface RS232. Includes power transformer and interface cables. (5V-3A)

Model

Part No.	Height / Width / Depth cm			Weight Kg
4120113	4	16	10	1



SPARE

Lamp of 4 mm Ø 5 V.

Part No. **4313030**