



Circulation pump “Ciragit”

FEATURES

The unit can be adjusted to fit most tanks or baths, can be submerged and held in position by an adjustable height support clamp.

Circulation pump allowing external circulation of liquids.
ON/ OFF illuminated power switch with light when on.

SPECIFICATIONS

OUTPUT water:
20 °C , 12 l/ min.- 150 mbar

OUTPUT silicon oil:
50 °C 4 l/min.- 130 mbar

Dimensions:
9.5 x 12.5 x 23 cm

Power: 30 W

Weight: 2.5 Kg.

Part No. **3000402**

