



Peristaltic pump "Percom N-M II"

ELECTRONIC SPEED CONTROL WITH DIGITAL DISPLAY,
WITH MULTI HEAD ADAPTER FOR TUBES FROM 4 TO 12 mm Ø (ext.). MAXIMUM WALL THICKNESS 1,6 mm.
ADJUSTABLE FLOW RATES FROM 0 TO 1100 ml/min.

FEATURES

Small footprint bench top model
Quiet with minimal vibrations.
Speed control 0,5 to 350 r.p.m (reversible).

CONTROL PANEL

Push button Stop/Start/switch.
Digital speed display.
Push button increase speed control.
Push button decrease speed control.


BACK

External control connection using a foot switch for Start / Stop, supplied as an accessory.

MODEL

Part No.	Height / Width / Depth (exterior) cm	Pump rate of H ₂ O ml/min	Power W	Weight Kg
3001301	19 13 23	0 to 1100	48	4.4

Graph of pump rate/r.p.m. dependent on the internal tube Ø.



Graph showing with liquids: Distilled water. These characteristics vary according to the density of the liquid.

NEW DESIGN

