



Peristaltic Pump "Percom S-Z BIG"

SPEED VARIABLE.

WITH MULTI-ADAPTER HEAD FOR 16 TO 20 MM Ø (EXT.) TUBES.

MAXIMUM WALL THICKNESS: 3,3 MM.

ADJUSTABLE FLOW FROM 220 TO 10000 ML/H.

FEATURES

High flow desktop model.

Easy-handling aluminum head.

Conveyor handle

Silent and vibration-free.

Selectable speed between 20 and 460 r.p.m.

CONTROL PANEL

Start-stop button.

Speed variator.

MODEL

Part No.	High/ Width /Depth (exterior) cm	Flow in H2O ml/hour	Power W	Weight Kg
5910155	23 21 40	220 to 10000	60	15,5

ACCESSORY Standard silicone tube. Ø Inside: 13.1 mm. Ø outside: 19.5 mm.

Part No. **5910156**

NEW

