



Membrane vacuum and pressure pump "Vacum-Pres"

FEATURES

Small footprint, bench top.

Quiet and vibration free.

Maximum working temperature: 40 °C.

CONTROL PANEL

Mains power switch.

Vacuum gauge, scale 0 to -1 bar.

Pressure gauge, scale 0 to 2.5 bar.

Control for pressure/vacuum.

Vacuum in.

Pressure out.

MODEL

Part No.	VACUUM		Pressure		Height / Width / Depth			Power W	Weight Kg
	Pump Rate /min.	Vacuum limit	Pressure Rate /min.	Max. Pressure	(exterior) cm				
3001002	1.08	-0.6 bar=440 mm Hg	1.3	1.2 bar	19	14	19	35	3.2

