

Free of lubrication • Free of maintenance • Quiet

APPLICATIONS

Vacuum filtration.
 Vacuum distillation.
 Rotating evaporators.
 Gas extraction and transport.
 Vacuum drying.

COMUN FEATURES

Constructed in aluminium.
 Membrane of polifluoropreno, resistant to the corrosion and long-life.
 Refrigeration system that allows a continuous operation of 24 hours.

MODEL “MV 50”



MODEL “MV 100”



CONTROL PANEL

Vacuum and pressure gauges.
 Adjustable pressure and vacuum.
 Tensile handle for easy transport.
 Button On-Off.

CONTROL PANEL

Vacuum indicator.
 Vacuum regulator.
 Button On-Off.

| MODELS | Part No. | Vacuum Rate m ³ /h | Vacuum Limit mbar | Pressure bar | Heigh / Width / Depth (exterior) cm | Power W | Power supply | Weight Kg |
|--------|----------|----------------------------------|----------------------|-----------------|--|------------|-----------------|--------------|
| MV 50 | 5910100 | 1,8 | 200 | 2 | 23,5 16 21 | 160 | 220V/50Hz | 8 |
| MV 100 | 5910101 | 3,6 | 200 | 2,75 | 23,5 16 30 | 160 | 220V/50Hz | 9 |

They are supplied with a 1 meter of silicone tube Ø 8 x 12mm.