



## Magnetic stirrer "Agimicro" without heating

### FEATURES

Maximum stir volume: 2 Litres.  
Adjustable stir speed from 40 to 1400 r.p.m.  
Durable moulded ABS plastic body.  
Excellent chemical resistance.

### MODEL

Part No.	Height cm	Ø Plate cm	Power W	Weight Kg
<b>7001638</b>	<b>5</b>	<b>12</b>	<b>3</b>	<b>0.5</b>

Comes complete with a 8 Ø x 32 mm P.T.F.E. coated stir bar.

### CONTROL PANEL

Mains switch with luminous "on" indicator.  
Analogue control of speed in r.p.m.

### ACCESSORIES

Support bar and electrode clamp (1 of 13 and 2 of 16 mm Ø). Used with pH and conductivity meters.  
Part No. **7001639**

