



Electronic temperature controller "Electemp-B"

ADJUSTABLE TEMPERATURE FROM 0 °C TO 310 °C.

FEATURES

Regulates the power to the heating element by pulsing electronic circuit, free of noise and interferences

APPLICATIONS

Regulates the power of temperature in ovens baths, radiators, hotplates, etc.

CONTROL PANEL

Temperature display indicator.
Set temperature display indicator.
Set value increase/decrease button.

BACK OF CONTROLLER

Output of power (Schuko plug system).
Protection fuse.
Connection for a PT100 probe (included)

MODEL

Part No.	Height / Width / Depth cm (exterior)	Max. power W (at 230 Vav)	Max. power W (at 110 Vav)	Weight Kg
3000896	15 15.5 19	2500	1700	1.2

NEW
DESIGN



"Electemp-B" application to a hotplate, showing the unit connected with an electronic contact thermometer.