



## Electronic contact thermometer "Sensoterm II"

DIGITAL DISPLAY, PT 1000 PROBE, AISI 304 STAINLESS STEEL.  
STABILITY  $\pm 1$  °C,  $\pm 1$  DIGIT. RESOLUTION  $\pm 1$  °C.



ONLY SUITABLE FOR INSTRUMENTS  
THAT CAN USE AN ELECTRICAL CONTACT  
THERMOMETER.

### FEATURES

Advantages include a separate temperature probe that maintain the main body away from steam and other vapours. Safer than a glass contact thermometer, eliminating the risk of broken it and a mercury spillage. Main body made of ABS plastic.

### TECHNICAL DATA

Temperature range: from  $-50$  °C to  $300$  °C.  
Temperature probe:  $3.5$  mm  $\varnothing$  x  $250$  mm immersion.  
Temperature Display:  $8$  mm high LCD display.  
Ambient working range:  $0$  to  $60$  °C and a relative humidity of  $80$  %.

### MODEL

Part No.	Height / Width / Depth cm (exterior)	Weight Kg
<b>1001009</b>	<b>13 7 2.4</b>	<b>0.250</b>

