



Rectangular hotplate "Placfin"

ADJUSTABLE HOT PLATE TEMPERATURE
FROM AMBIENT +5 °C UP TO 130 °C. STABILITY ± 3 °C.
SUITABLE FOR THERMAL TESTS (BUT NOT BOILING).

FEATURES

Temperature control by a hydraulic thermostat.

Protected thermometer from 0 °C to 130 °C.

Hotplate surface protected by chemically treated aluminium, with heating elements distributed throughout the hotplate's surface.

Enamelled external case, resistant against corrosive chemicals.

CONTROL PANEL

Hydraulic thermostat control, synchronised with an indicator lamp when in operation.

MODEL

Part No.	Width / Depth (plate) cm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
7000497	23 22	11 22 30	1500	3.7



ACCESSORY

Slide support adapter made of AISI 304 stainless steel. Capacity: up to 18 slide of 76 x 255 mm (in a horizontal position).

Part No. **7000498**

