




Hotplate-Rectangular “Combiplac”

PLATE TEMPERATURE ADJUSTABLE UP TO 400 °C.

FEATURES

Electronic energy temperature control.

The non warping steel hotplate has heating elements embedded throughout the hotplate surface maintaining the heat distribution. 

The main case is made of AISI 304 stainless steel.

CONTROL PANEL

Main power switch with “ON” indicator lamp.

Electronic energy control from 0 to 100 %.

An indicator lamp shows when the unit is hot.

Connector at the back for the external digital thermometer, “Sensoterm” or “Electemp-TFT”.

See Regulators and Controllers (pages 293 and 294).

MODELS

Part No.	Width / Depth (plate) cm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
3000156	20 40	12 22 55	2300	9
3000718	30 30	13 33 46	2800	13



ACCESSORIES

Sand bath adapters. Made entirely of AISI 310 stainless steel, removable for easy cleaning.

Part No.	Hotplate Part No.	Height/Width/Depth (Sand tray dims) cm
6000719	3000156	5.5 20 40
6000720	3000718	5.5 31 31

