



## Cooling plate "Plac-Center"

DIGITAL DISPLAY WITH ANALOGUE TEMPERATURE CONTROLLER.  
ADJUSTABLE TEMPERATURES FROM -10 TO 5 °C.

### APPLICATIONS

Designed to keep material in a solid state.

### FEATURES

External case made of AISI 304 stainless steel. Cooling plate made of black anodised aluminium.  
Anti vibration, hermetically sealed CFC free compressor.

### CONTROL PANEL

Illuminated Mains ON/OFF switch.

Analogue electronic temperature control with digital display (2 digits).

Resolution:  $\pm 1$  digit.

Temperature probe: PTC.

### MODEL

Part No.	Width/Depth (plate) cm	Height/Width/Depth (exterior) cm	Power HP	Power W	Weight Kg
<b>4000501</b>	<b>30 15</b>	<b>27 33 50</b>	<b>1/3</b>	<b>160</b>	<b>15</b>

