



Infrared tip sterilizer "Sterilbio"

APPLICATIONS

Sterilization of culture tips, needles and small instruments. Incinerates organic material at 900 °C for 5/8 seconds.

FEATURES

Made from ceramic elements with a quality quartz tube interior.
Does not need purged oxygen and can be used in anaerobic environments.
Eliminates problems of aerosols and bacterial contamination.
Support at the side for 4 Kollé handles, needles, lancet or dissecting forceps etc.

MODEL

Part No.	Height/Width/Depth (exterior) cm			Ø/Length (tube interior) cm		Power W	Weight Kg
3000788	20	13	18	1.4	15	110	1.5

