



# Incubator for Petri capsules

NATURAL CONVECTION.

MICROPROCESSOR REGULATION AND TEMPERATURE DIGITAL CONTROL.

FOR ADJUSTABLE TEMPERATURES FROM AMBIENT +5°C TO 60°C.

STABILITY: ±0,1°C TO 37°C. HOMOGENEITY: ±0,1°C TO 37 °C. SETPOINT ERROR: ±0,1°C. RESOLUTION: 0,1°C.

*Small size. Culture visual control. Transportable.*

## APPLICATIONS

Specially designed for bacteria and fungi cultures in Petri capsules at the same temperature of human body.

## FEATURES

Culture surface 320 x 220 mm (Inner height: 20mm)

Culture visual monitoring.

Transparent cover.

Easy access to samples.

**Approximate capacity:** (single level) (mm)

15 Petri capsules of Ø55.

10 Petri capsules of Ø80.

7 Petri capsules of Ø90.

6 Petri capsules of Ø100.

3 Petri capsules of 120x120.

2 Petri capsules of Ø140.

## MODEL

Part No.	Margin T° °C	Height / Width / Depth (Exterior) cm	Height / Width / Depth (interior) cm	Power W.	Weight Kg
<b>4002629</b>	Amb +5 to 60	9 45 24	2 32 22	200	3.5

