



Refrigerated ovens with fixed humidity Hotcold Humidity F & J

WITH FORCED AIR CIRCULATION
DIGITAL TEMPERATURE ELECTRONIC REGULATION
FOR ADJUSTABLE TEMPERATURES FROM 6°C TO 20°C.



SAFETY:

DIN 12880 STANDARDS. ADJUSTABLE SAFETY THERMOSTAT INCORPORATED WHICH DISCONNECTS HEATING IN CASE OF THE HEATER'S REGULATOR FAILURE, WITH MANUAL RESET.

ACCURACY - TEMPERATURE TABLE:

STABILITY:..... ±0,5 °C
HOMOGENEITY: ±1 °C
SETPOINT ERROR: ±1 °C
RESOLUTION 0,1 °C

APPLICATIONS

Refrigerated oven for conservation and products storage which require temperatures below ambient and a controlled fixed humidity.

FEATURES

- Colour TFT touch screen.
- Exterior case, door and interior enclosure covered in stainless steel AISI 304.
- Door with lock and easy to replace gasket, with automatic return mechanism and manual lockout. Resistant to chemical agents.
- Tropicalized hermetic compressor group.
- Automatic defrosting with water evaporation due to tempered electric defrost.
- Supports for shelves adjustable in height.
- Front foot supports adjustable in level.
- Integrated interior light.
- Entry port for external probes introduction.
- Supplied with wheels.



HOTCOLD F

TEMPERATURE RANGES- HUMIDITY

TEMPERATURE	HUMIDITY	PRECISION
6 °C:	33% HUMIDITY	±3 %.
8 °C:	31% HUMIDITY	±3 %.
12 °C:	30% HUMIDITY	±3 %.
16 °C:	29% HUMIDITY	±3 %.
20 °C:	26% HUMIDITY	±3 %.

HOTCOLD J

TEMPERATURE RANGES - HUMIDITY

TEMPERATURE	HUMIDITY	ACCURACY
6 °C:	60% HUMIDITY	±3 %.
8 °C:	57% HUMIDITY	±3 %.
12 °C:	53% HUMIDITY	±3 %.
16 °C:	50% HUMIDITY	±3 %.
20 °C:	47% HUMIDITY	±3 %.

CONTROL PANEL

- Model with 4.3 inches TFT touch screen.
- Main switch.
- Visual alarm.
- Probe error detection.
- Self Diagnostics.
- Over humidity alarm.
- Safety thermostat (TS) by software and mechanical.



MODEL

HOTCOLD	Part No.	Range °C	Humidity Range	Capacity litres	Height / Width / Depth (interior) cm	Height / Width / Depth (exterior) cm	Trays number	Power W	Weight Kg
F	2101517	6 +20	26% a 33%	600	138 62 68	210 73 84	3	1350	125
J	2101518	6 +20	47% a 60%	600	138 62 68	210 73 84	3	1350	125

SPARES

Tray + 2 guides.
Part No. 1001858

Supplied with 3 steel trays, PVC laminated.

ACCESSORIES

- Data logger with digital screen "View 2"
- Working temperature from -25°C to +85 °C.
- Humidity 0-100% RH
- IP65 protection.
- Capacity up to 30000 logs.
- Two programmable alarms.
- Log intervals of 1 second up to 10 days.
- Data download option on standby or operating.
- Display mode in °C or °F.
- Part No. 2101508



OPTIONAL COMPLEMENTS

- ETHERNET thermometer module**
Temperature Record, with control of programmable warnings and graphic representation. It must be installed at the factory.
ETHERNET connection is needed at the place of use.
Part No. 2101627
- Data logger with digital display**
PT100 type probe with two meters of cable. Measuring range of -100 °C. at +200 °C. With GSM module for data transmission.
Part No. 2101520