

# COMECTA Precise analytical balance “FA-2204B”

COMECTA

## FEATURES

- Digital 7 segment screen.
- Time to stability: 8 seconds.
- Adjustable time and sensitivity control.
- ON/OFF automatic zero adjustment.
- Manual external 220g calibration.
- Tare range by subtraction up to the maximum capacity.
- Piece counting function.
- Density measurement in the base of the unit.
- Units of measurement in g, ct, and oz.
- Level adjustment (bubble).
- Data interface: RS-232 for connection to a printer.
- Operable ambient temperature 5 to 40 °C
- Power supply: 110/220 V-50/60Hz.
- Supplied complete with cover.



MODELS	Part No.	Maximum weight capacity g	Precision $\pm$ g	Linearity $\pm$ g	Height / Width / Depth (exterior) cm	Height / Width / Depth (Draft shield) cm	$\emptyset$ plate cm	Weight Kg
FA-2204B	5830039	220	0.0001	0.0002	35 21 34	24 19 17.5	8	8.5