




Hotplate double circular “Duplac”

FOR PLATE TEMPERATURES UP TO 400 °C.

FEATURES

Hot plate with non warping cast steel plate with embedded heating elements throughout the hot plate surface. 

Independent plate temperature control is maintained by pulsed energy controllers. The exterior is epoxy covered steel.

CONTROL PANEL

Two temperature energy controllers.

Two indicator lamps function when the heaters are on.

MODEL

Part No.	Plate Ø cm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
4000634	19	15 47 24	1750	8.5

