



50 Litre distilled water reservoir

FEATURES

Recommended for the storage of distilled water. Made of high density polyethylene. 50 litre capacity, with handles and dispenser tap. 15cm Ø top lid with stopper. Dimensions: 68 cm high x 35 cm Ø. Part No. **0106006**

ACCESSORY

Water reservoir trolley.

Robust steel trolley with epoxy coating. Lockable wheels, Inner section for holding jars and carboys etc. Fits containers 38 cm Ø. Height 48 cm. Part No. **5903049**

