



Water distiller, specially for autoclaves "Dest-4"

DISTILLATION CAPACITY: 1.2 LITRES/HOUR. INTERNAL CAPACITY 4 LITRES.

FEATURES

Case made from thermo resistant plastic. Internal stainless steel. Cooling by forced air through a condenser. Without glass heater elements, Easy to use. Safety cut out if over temperature due to insufficient water.

Supplied complete with a 4 litres plastic bottle and three dispenser tap with filter. Note: for daily use it is recommended to change the filter every 3 months.

MODEL

Part No.	Distillation capacity l/h	Conductivity $\mu\text{s}/\text{cm}$	\emptyset	Height cm	Power W	Weight Kg
4001729	1.2	5	29	39	750	3.5

No external additional water required.



SPARE. Dispensing filter tip. Part No. **4001730**