



## Distillation still "L-3"

**DISTILLATION CAPACITY: 3 LITRE PER HOUR.**

### FEATURES

Made entirely of borosilicate glass.

Functions: automatic with continuous water production. Produces distilled water of a high purity, free of pyrogens and metal ions.

Equipped with quartz heater elements and protection device using a safety thermostat that is activated in the case of low water pressure to the condenser or boiler unit. The unit automatically reactivates when the water pressure is restored.

Suitable for wall mounting.

### MODEL

Part No.	Distillation capacity l/h	Water consumption l/min	Conductivity at 20 °C	Height / Width / Depth (exterior) cm	Power W	Weight Kg
<b>4903000</b>	<b>3</b>	<b>1</b>	<b>3 µs/cm</b>	<b>34.5 56 19</b>	<b>2200</b>	<b>4.5</b>

