



Water distillers “R-4 Reser” and “R-8 Reser”

**DISTILLATION CAPACITY: 4 AND 8 LITRES/HOUR.
STORAGE TANK: 8 LITRES.**

FEATURES

- External case epoxy covered.
- Cooling and heating elements made of stainless steel.
- Easy to dismantle for cleaning.
- Thermostat that disconnects the heater when there's lack of water supply.
- Monitoring function that disconnects the distiller when the storage tank is complete full, and resumes as it empties.
- Distilled water dispenser located in the front part.
- Distilled water and water supply records of 3/4 inch.
- Desktop model with device to hang it on the wall.

CONTROL PANEL

- Illuminated power switch.
- Lack of water flow indicator lamp.



MODEL

Part No.	Distillation capacity l/h	Water storage tank l	Water consumption l/min	Conductivity at 20° C	Height / Width / Depth (exterior) cm	Power W	Weight Kg
4903005	4	8	1	1,5 µS/cm	50 63 32	1x3000	24
4903017	8	16	1,4	1,5 µS/cm	60 71 34	2x3000	34