



## Water distiller “D-4 Large” and “D-10 Large”

DISTILLATION CAPACITY: 4 AND 10 LITRES/HOUR.

### FEATURES

Equipment entirely made of borosilicate glass.

Automatic operation.

Continuous production of high purity distilled water, free of dissolved gases, metals, salts and pyrogens.

Spiral heating system.

Double safety system:

1. By a presostat that turns the heater off in case of lack of water supply. 2. Water level device that avoids the computer to function until the heating element is covered with water.

### CONTROL PANEL

Power switch.

Lack of water flow indicator lamp.

### MODEL

| Part No.       | Distillation capacity l/h | Water consumption l/min | Conductivity at 20° C | Height / Width / Depth (exterior) cm | Power W | Weight Kg |
|----------------|---------------------------|-------------------------|-----------------------|--------------------------------------|---------|-----------|
| <b>4903001</b> | 4                         | 1                       | 2 µs/cm               | 69 28 30                             | 3000    | 6,3       |
| <b>4903015</b> | 10                        | 2                       | 1,5 µS/cm             | 96 36 40                             | 6600    | 23        |

NEW

