



## Water distillers, “AC-L4” and “AC-L8”

**DISTILLATION CAPACITY: 3.5 AND 8 LITRES/HOUR.**

### FEATURES

Tank and heating element made from stainless steel.

External case made from epoxy covered steel.

Easily dismantled for cleaning and service. All seals made from silicon.

In the case of over temperature or low water supply, a safety thermostat cuts power to the heater.

Water supply tube 3/4” mm.

Distilled water tube connection Ø 12 mm.

Voltage: 230 V / 50/60 Hz

### CONTROL PANEL

Illuminated power switch.

Safety switch thermostat with rearm system.

MODELS	Part No.	Distillation capacity l/h	Water consumption l/min	Conductivity at 20 °C	Height / Width / Depth cm	Power W	Weight Kg
AC-L4	4903007	4	1	2 µs/cm	44 25 23	3000	8
AC-L8	4903008	8	1,4	2 µs/cm	61 26 26	6000	14

