

COMECTA Pocket Balance "Pocket Scale"

FEATURES

Capacity up to 200 g, measuring in g, oz and qty.

Digital LCD display.

Stabilization time ± 3 seconds.

Tare by subtraction to maximum weight.

Piece counting function and automatic reset to zero.

Automatic calibration, external.

Batteries 3 x 1.5 V.

Lid to cover the weighing pan and operation menu.

MODEL

Part No.	Max. weight capacity g	Precision \pm g	Height / Width / Depth (exterior) mm	Width / Depth (plate) mm	Weight Kg
5830024	200	0.1	24 74 120	74 71	0.13

Comes with protective cover.

