



Heating mantle for flasks and funnels “Fibroman-0”

WITH TRIPOD SUPPORT STAND.

APPLICATIONS

Suitable for round bottom flasks, decanting funnels and reaction flasks with a bottom drain.

FEATURES

- Woven mineral fibre mantle.
- Epoxy coated exterior.
- On/Off switch.
- Safety earth.

MODELS

Part No.	Flask capacity ml	Ø /Height (Exterior dims.) cm		Ø cm Bottom hole	Tripod Height cm	Power W	Weight Kg
3049110	1000	24	10.5	8	34	410	1.5
3049120	2000	28.5	12	8	36	530	2.5
3049150	5000	37.5	17.5	9	38	840	3.5
3491100	10000	43.5	20	9	43	1400	5
3491200	20000	48	25	12	36	2300	13

Due to its capacity the part no. 3491200, 20000 ml capacity is manufactured in a cylindrical case. See accessories (page 230) and controller (pages 293 and 294).

