



Wax dispenser "Dispenser"

ADJUSTABLE TEMPERATURE FROM AMBIENT +5 TO 100 °C.
STABILITY ± 1.5 °C.

FEATURES

External metal case coated with epoxy. Internal tank made of aluminium coated with black indelible PTFE. Independently heated tap.

CONTROL PANEL

Illuminated Mains switch.

Indicator lamp shows when the heater is on.

Adjustable temperature thermostat controller.

Adjustable tap temperature controller.

MODEL

Part No.	Capacity litres	Ø / Height (usable) cm	Ø / Height (exterior) cm	Power W	Weight Kg
4000490	4	20 15	28 30	400	5

