



# Magnetic Stirrers “Agimatic” and “Agiman” analogue control

FOR SPEEDS FROM 60 TO 1600 rpm..

MAXIMUM STIR VOLUME: 10 LITRES.

PLATE MADE FROM AISI 304 STAINLESS STEEL WITH REMOVABLE SAFETY RING. PLATE Ø 14.5 cm.



## “Agimatic-S” without heating

### CONTROL PANEL

1. Mains switch illuminated.
2. Electronic stir speed control.



### MODEL

	Part No.	Maximum stir volume litres	Height/width/Depth (exterior) cm	Power W	Weight Kg
<b>AGIMATIC-S</b>	<b>7000270</b>	up to 10	14,5 19 26	90	3,5

Supplied with a PTFE coated stir bar of Ø 8 x 40 mm long.

