




## “Agimatic-E” with heating

ELECTRONIC TEMPERATURE CONTROL OF SPEED AND POWER TO THE HEATER.  
FROM 50 °C TO 350 °C PRECISION  $\pm 2$  °C.

### CONTROL PANEL

1. Switch on button with indicator lamp for heating and for stirring speed.
2. Electronic stir speed control.
3. Electronic temperature control of the heater plate from 50 °C to 350 °C synchronised with an indicator lamp that functions when the plate is hot.
4. Illuminated heater indicator when the heater is “on”. Precision  $\pm 2$  °C.



### MODEL

Part No.	Maximum stir volume litres	Height/width/Depth (exterior) cm	Power W	Weight Kg
<b>AGIMATIC-E 7000273</b>	<b>up to 10</b>	<b>14,5 19 26</b>	<b>550</b>	<b>3,5</b>

Supplied with a PTFE coated stir bar of  $\varnothing 8 \times 40$  mm long.

