



# Magnetic Stirrers “Agimatic-ED” and “Agimatic-ED-C” digital control



FOR SPEEDS FROM 60 TO 1600 RPM.

FOR ADJUSTABLE TEMPERATURES FROM 50 °C TO 350 °C (200 °C WITH EXTERNAL PROBE). PRECISION  $\pm 2$  °C.

MAXIMUM STIR VOLUME: 10 LITRES.

PLATE Ø IS 14,5 cm.

WITH RENOVABLE SAFETY RING AGAINST SPILLAGE.



## “Agimatic-ED” with heating

ELECTRONIC TEMPERATURE CONTROL OF SPEED AND POWER TO THE HEATER.

SILUMIN UPPER PLATE COVERED BY AISI 304 STAINLESS STEEL.

### CONTROL PANEL

1. Digital display indicating stir speed or temperature.
2. Push button, parameter selector, Stir speed or Temperature.
3. Push button, increase parameter value.
4. Push button, decrease parameter value.
5. Indicator lamp, safety thermostat “on”
6. Push button Stop-Start

Precision  $\pm 2$  °C.



### MODEL

Part No.	Maximum stir volume litres	Height/width/Depth (exterior) cm	Power W	Weight Kg
<b>AGIMATIC-ED 7000280</b>	<b>up to 10</b>	<b>14,5 19 26</b>	<b>550</b>	<b>3,5</b>

Supplied with a PTFE coated stir bar of Ø 8 x 40 mm long.

