



## “Agimatic-ED-C” with heating

ELECTRONIC TEMPERATURE CONTROL OF SPEED AND POWER TO THE HEATER.

SILUMIN UPPER PLATE WITH CERAMIC SURFACE.

### CONTROL PANEL

1. Digital display indicating stir speed or temperature.
  2. Push button, parameter selector, Stir speed or Temperature.
  3. Push button, increase parameter value.
  4. Push button, decrease parameter value.
  5. Indicator lamp, safety thermostat “on”
  6. Push button Stop-Start
- Precision  $\pm 2$  °C.



### MODEL

Part No.	Maximum stir volume litres	Height/width/Depth (exterior) cm	Power W	Weight Kg
<b>AGIMATIC-ED-C 7000281</b>	<b>up to 10</b>	<b>14,5 19 26</b>	<b>550</b>	<b>3,5</b>

Supplied with a PTFE coated stir bar of  $\varnothing 8 \times 40$  mm long.

### Ceramic surface:

Excelent conductivity and heat distribution.

Unalterable surface.

Higher acceleration to reach temperature.

Extreme resistance against chemical agents and scratches.



Complements for Agimatic-ED and Agimatic ED-C.  
See accessories on pages 31 and 32.