



Multiple position magnetic stirrers "Multimatic"



COMMON FEATURES

External plate made from AISI 304 stainless steel with engraved stirrer positions.
 Independent stirring speed control for each position.
 Front panel protected against splashes.

Models 5-S and 5-N: Suitable for 5 flasks, 600 ml max. (max. Ø 90 mm).
 9-S and 9-N: Suitable for 9 flasks, 300 ml max. (max. Ø 75 mm).
 Epoxy coated external case.

APPLICATIONS

Specially recommended for bio-pharmaceutical diagnostic assays.



Stirrers "Multimatic" 5-S and 9-S without heating FOR 5 AND 9 POSITIONS.



CONTROL PANEL

Illuminated mains switch.
 Individual position stir speed control.

MODELS

	Part No.	Height / Width / Depth (exterior) cm	Width / Depth (plate) cm	Speed control r.p.m.	Maximum stir volume ml	Power W	Weight Kg
MULTIMATIC 5-S	7000951	11 28 35	28 33	Up to 1600	5 x 600	16	7
MULTIMATIC 9-S	7000952	11 28 35	28 33	Up to 1600	9 x 300	27	9

Comes complete with a 5 or 9 mm stir bar, model dependent covered in PTFE of 8 Ø x 32 mm long.

ACCESSORIES

Transparent methacrylate bath.
 Fits onto the non heating stir plate of the "Multimatic" 5-S y 9-S
 The bath can withstand temperatures up to 60 °C.
 See section immersion Thermostats (page 94).

Part No.	Capacity litres	Height / Width / Depth (exterior) cm	Weight Kg
1001007	14	16 30 55	4

