



Powerful magnetic stirrer “Agimatic-H” without heating

FOR LOW VISCOSITY VOLUMES FROM 1 TO 60 LITRES.

USING REDUCED VOLUMES THE UNIT CAN STIR MEDIUM TO HIGH VISCOSITY LIQUIDS.

Maximum stirring power

FEATURES

Electronic speed control from 50 to 980 r.p.m, will maintain a constant stir speed irrespective of the liquid viscosity.

Digital speed display.

AISI 304 stainless steel top plate with an epoxy coated metal case.

CONTROL PANEL

Main switch with a luminous “on” indicator.

Electronic speed control.

Digital speed display in r.p.m.

MODEL

	Part No.	Speed control r.p.m.	Maximum stir volume in litres	Height/Width/Depth (exterior) cm	Power W	Weight Kg
AGIMATIC-HS	7001005	50 to 980	up to 25	12 28 34	35	8
AGIMATIC-HL	7001006	50 to 980	up to 60	14 44 38	40	12

Model HS comes complete with a 16 Ø x 127 mm stir bar and model HL with a 19 Ø x 150 mm stir bar.

