




Complete KIT stirrer Agimatic-ED with support rod and temperature sensor PT100. With heating

SUPPLIED COMPLETE WITH SUPPORT ROD, TEMPERATURE SENSOR PT100, CLAMP AND STIR BAR.
ELECTRONIC TEMPERATURE CONTROL OF SPEED AND POWER TO THE HEATER.

SILUMIN UPPER PLATE COVERED BY AISI 304 STAINLESS STEEL.

CONTROL PANEL

1. Digital display indicating stir speed or temperature.
 2. Push button, parameter selector, Stir speed or Temperature.
 3. Push button, increase parameter value.
 4. Push button, decrease parameter value.
 5. Indicator lamp, safety thermostat "on"
 6. Push button Stop-Start
- Precision ± 2 °C.



MODEL

	Part No.	Maximum stir volume litres	Height/width/Depth (exterior) cm	Power W	Weight Kg
KIT	7001512	up to 10	14,5 19 26	550	3,5

Supplied with a PTFE coated stir bar of $\varnothing 8 \times 40$ mm long.

INCLUDED

Support rod with thread of 12 mm \varnothing x 500 mm high, stirrer dependent, attaches to the main body of the stirrer.

Temperature sensor Pt100. Controls the sample temperature in the digital hotplate stirrers.

Immersion Probe: 4mm \varnothing x 130 mm long.

Cable connector 150 cm long.

Clamp with bosshead, "AUTOMATIC" Chrome coated alloy with screw closure.

For \varnothing from 0 to 25mm.

