



## Desiccator for materials

WITH HYDROMETER CONTROL.

### APPLICATIONS

Cabinet with protection against humidity and dust for anhydrous, biological and chemical preservation of samples.

### FEATURE

Made of robust transparent 12mm thick methacrylate. The door has a silicon seal and magnetic catch.

Volume: 55 Litres.

Dimensions 50 cm high x 38 cm wide x 29 cm deep.

Supplied complete with three perforated shelves and a stainless steel AISI 304 tray to hold desiccating material.

Part No. **1001403**

