



Glass drying oven "Dryglass"

FAN ASSISTED AIR CIRCULATION.

ADJUSTABLE TEMPERATURE FROM 40 °C TO 170 °C.

SAFETY:

EN.61012 STANDARD OVER TEMPERATURE SAFETY CUT OUT FITTED.
DIN 12880. STANDARD (CLASS 2 AND 3.1) ADJUSTABLE SAFETY THERMOSTAT FITTED.

FEATURES

- Hydraulic thermostat for temperature control.
- Air circulation by turbo fan.
- Inner chamber made of AISI 304 stainless steel with shelf runners.
- Removable tempered glass sliding doors.
- Ventilation port for steam.
- Epoxy coated external case.

STANDARD EQUIPMENT

- 2 shelves and 4 shelf guides.

CONTROL PANEL

- Dual heating power selector switch.
- Mains indicator lamp.
- Hydraulic thermostat for temperature control.
- Locking system of thermostat knob.
- Heater "ON" operation indicator lamp.
- Analogue thermometer.
- Adjustable over temperature safety thermostat, that cuts off the power if the control thermostat fails, manual reset with "on" indicator lamp.



MODEL

Part No.	Capacity litres	Height / Width / Depth (interior) cm			Height / Width / Depth (exterior) cm			Shelf Positions	Power W	Weight Kg
2000381	126	45	70	40	66	94	54	8	3000	65

ACCESSORIES

Accessories must be factory installed.



Part No.

2000002 Timer switch 0-120 minutes.

2000009 24 hour programmer with continuous on/off cycling up to every 15 minutes.

SPARES

Part No.

2000081 Shelf guides x 2.

2000091 Shelf.

Each shelf requires 2 guides.